

This Page Is Inserted by IFW Operations  
and is not a part of the Official Record

## **BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images may include (but are not limited to):

- BLACK BORDERS
- TEXT CUT OFF AT TOP, BOTTOM OR SIDES
- FADED TEXT
- ILLEGIBLE TEXT
- SKEWED/SLANTED IMAGES
- COLORED PHOTOS
- BLACK OR VERY BLACK AND WHITE DARK PHOTOS
- GRAY SCALE DOCUMENTS

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning documents *will not* correct images,  
please do not report the images to the  
Image Problem Mailbox.**

**THIS PAGE BLANK (USPTO)**

**PCT**WORLD INTELLECTUAL PROPERTY ORGANIZATION  
International Bureau

## INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification <sup>6</sup> : <b>C07K 14/245, C12N 9/24, 9/52, G01N 33/53</b>	<b>A1</b>	(11) International Publication Number: <b>WO 96/16082</b>
		(43) International Publication Date: <b>30 May 1996 (30.05.96)</b>

(21) International Application Number: **PCT/SE95/00761**(22) International Filing Date: **21 June 1995 (21.06.95)**(30) Priority Data:  
**9404072-2**      **24 November 1994 (24.11.94)**      **SE**(71) Applicant (for all designated States except US): **ASTRA AKTIEBOLAG [SE/SE]; S-151 85 Södertälje (SE).**

(72) Inventors; and

(75) Inventors/Applicants (for US only): **BALGANESH, Tanjore, Soundararajan [IN/IN]; 5th Main, Postal Colony II Stage, Sanjaynagar, Bangalore 560094 (IN). TOWN, Christine, Mary [GB/SE]; Högloftsvägen 31, S-151 68 Södertälje (SE).**(74) Agent: **ASTRA AKTIEBOLAG; Patent Dept., S-151 85 Södertälje (SE).**(81) Designated States: **AM, AT, AU, BB, BG, BR, BY, CA, CH, CN, CZ, DE, DK, EE, ES, FI, GB, GE, HU, IS, JP, KE, KG, KP, KR, KZ, LK, LR, LT, LU, LV, MD, MG, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, TJ, TM, TT, UA, UG, US, UZ, VN, European patent (AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, ML, MR, NE, SN, TD, TG), ARIPO patent (KE, MW, SD, SZ, UG).****Published***With international search report.*(54) Title: **NOVEL POLYPEPTIDES**

(57) Abstract

The present invention relates to variants of Penicillin Binding Proteins (PBP), which proteins are involved in bacterial peptidoglycan biosynthesis. Disclosed are also DNA molecules coding for the said PBP variants, as well as vectors and cells harbouring such DNA molecules. The invention is also related to processes for assaying and designing therapeutically useful compounds which have high affinity to PBP, which processes utilize the said PBP variants.

**FOR THE PURPOSES OF INFORMATION ONLY**

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AT	Austria	GB	United Kingdom	MR	Mauritania
AU	Australia	GE	Georgia	MW	Malawi
BB	Barbados	GN	Guinea	NE	Niger
BE	Belgium	GR	Greece	NL	Netherlands
BF	Burkina Faso	HU	Hungary	NO	Norway
BG	Bulgaria	IE	Ireland	NZ	New Zealand
BJ	Benin	IT	Italy	PL	Poland
BR	Brazil	JP	Japan	PT	Portugal
BY	Belarus	KE	Kenya	RO	Romania
CA	Canada	KG	Kyrgyzstan	RU	Russian Federation
CF	Central African Republic	KP	Democratic People's Republic of Korea	SD	Sudan
CG	Congo	KR	Republic of Korea	SE	Sweden
CH	Switzerland	KZ	Kazakhstan	SI	Slovenia
CI	Côte d'Ivoire	LJ	Liechtenstein	SK	Slovakia
CM	Cameroon	LK	Sri Lanka	SN	Senegal
CN	China	LU	Luxembourg	TD	Chad
CS	Czechoslovakia	LV	Latvia	TG	Togo
CZ	Czech Republic	MC	Monaco	TJ	Tajikistan
DE	Germany	MD	Republic of Moldova	TT	Trinidad and Tobago
DK	Denmark	MG	Madagascar	UA	Ukraine
ES	Spain	ML	Mali	US	United States of America
FI	Finland	MN	Mongolia	UZ	Uzbekistan
FR	France			VN	Viet Nam
GA	Gabon				

## NOVEL POLYPEPTIDES

### TECHNICAL FIELD

5 The present invention relates to variants of Penicillin Binding Proteins (PBP), which proteins are involved in bacterial peptidoglycan biosynthesis. Disclosed are also DNA molecules coding for the said PBP variants, as well as vectors and cells harbouring such DNA molecules. The invention is also related to processes for assaying and designing therapeutically useful  
10 compounds which have high affinity to PBP, which processes utilize the said PBP variants.

### BACKGROUND ART

15 Bacteria and most other unicellular organisms possess a cell wall, which comprises a cross-linked polysaccharide-peptide complex called peptidoglycan. Peptidoglycan biosynthesis consists of three stages: (1) synthesis of precursors (sugar nucleotides) in the cytosol, (2) precursor transfer across the membrane and formation of the polysaccharide chain,  
20 and (3) cross-linking of individual peptidoglycan strands in the cell wall.

In the latter stage of peptidoglycan biosynthesis, new bonds must be made between nascent glycan strands and existing peptidoglycan. The newly synthesized chains are about 10 disaccharides long and are extended by  
25 transglycosylase enzymes to a final glycan strand of between 100 and 150 disaccharide units. The peptidoglycan is crosslinked by the action of transpeptidases which link the terminal D-ala of one glycan strand to a free  $\epsilon$ -amino group on a diaminopimelic acid residue on an adjacent region.

30

A number of antibiotics inhibit bacterial growth by interfering with the formation of the peptidoglycan layer. The cross-linking reaction is the

target for action of two important classes of such antibiotics, the penicillins and the cephalosporins. Penicillin is thought to react irreversibly with the transpeptidase that catalyses cross-linking.

- 5 The penicillin interactive proteins fall into three groups: the  $\beta$ -lactamases, the Low Molecular Weight-Penicillin Binding Proteins (PBPs), which mainly include the carboxypeptidases, and the High Molecular Weight-Penicillin Binding Proteins. Penicillin Binding Proteins are those enzymes which have been shown to bind radiolabelled penicillin G. In *Escherichia*  
10 *coli* such proteins are called e.g. PBP 1A and PBP 1B, both belonging to the class High Molecular Weight-PBPs. PBP 1A and 1B, which are known to be membrane bound proteins, maintain cell integrity and control peptidoglycan side wall extension during growth. Inactivation of either PBP 1A or PBP 1B can be tolerated by the bacteria while the deletion of  
15 both the genes, designated *ponA* and *ponB*, is lethal (Yousif et al., 1985).

PBP 1B is known to be a bifunctional enzyme possessing both transpeptidase and transglycosylase activity (Ishino et al., 1980). PBP 1A is believed to be bifunctional since it can substitute for PBP 1B. The  $\beta$ -lactam  
20 antibiotics, such as penicillin, inhibit only the transpeptidase activity of these proteins.

The transglycosylase reaction is inhibited by e.g. moenomycin, which is a phosphoglycolipid used as a growth promoter in animal nutrition and  
25 which has been shown to possess broad spectrum bactericidal activity. The enzyme transglycosylase has been shown to be present in *Escherichia coli*, *Staphylococcus aureus*, *Bacillus megaterium* and *Bacillus subtilis*. This suggests that interference of peptidoglycan biosynthesis by inhibition of transglycosylase could be a lethal event in all clinically important  
30 pathogens.

The putative transglycosylase domain of PBP 1B has been assigned to the N-terminal 478 amino acids (Nakagawa et al., 1987). This regions includes three conserved stretches of amino acids between the N-terminal half of both PBP 1A and 1B and could represent residues involved in the transglycosylase activity.

Preparation of Penicillin Binding Protein 2A from *Staphylococcus aureus* is disclosed in EP-A-0505151.

## 10 DISCLOSURE OF THE INVENTION

There is a growing number of reports of bacteria which are resistant to antibiotics. There is consequently a need for new compounds which inhibit bacterial growth by means of binding Penicillin Binding Proteins. The present invention provides PBP variants which facilitate processes for assaying and designing therapeutically useful compounds which have high affinity to PBPs.

Accordingly, it is an object of the invention to provide polypeptides which are water-soluble active derivatives of bacterial bifunctional Penicillin Binding Proteins, said Penicillin Binding Proteins being bound to the cell membrane when expressed in a bacterial cell and being capable of exhibiting both transglycosylase and transpeptidase activities and said derivatives lacking a membrane anchoring sequence but retaining the capability to exhibit one or both of said enzymic activities. The "bacterial cell" mentioned above is preferably an *Escherichia coli* cell or a *Streptococcus pneumoniae* cell.

The soluble PBP variants according to the invention retains transglycosylase activity, indicating that soluble variants of PBP, devoid of membrane anchoring sequences, can recognize lipid linked substrate and polymerise the disaccharide into repeating units. It can thus be assumed

that other analogues of PBP lacking residues involved in membrane attachment would be enzymatically functional.

5 Molecules interacting with the penicillin interactive region of soluble PBP variants could be assumed to be capable of interacting identically with wild-type PBPs. Consequently, the soluble PBP variants according to the invention can be used for identifying compounds which are interacting with wild-type Penicillin Binding Proteins.

10 It is furthermore well known that membrane-bound proteins are very difficult to crystallize. The soluble enzymatically active PBP variants can be used for crystallisation and will thereby facilitate a rational design, based on X-ray crystallography, of therapeutic compounds inhibiting High Molecular Weight-PBPs.

15 A further object of the invention is to provide polypeptides which are truncated water-soluble derivatives of bacterial bifunctional Penicillin Binding Proteins, said Penicillin Binding Proteins being bound to the cell membrane when expressed in a bacterial cell and being capable of  
20 exhibiting both transglycosylase and transpeptidase activities and said derivatives lacking the membrane anchoring sequence but retaining the capability to exhibit the transglycosylase activity. The "bacterial cell" mentioned above is preferably an *Escherichia coli* cell.

25 Alignment of amino acid sequences of High Molecular Weight-Penicillin Binding Proteins, and the compilation of the motifs involved in the penicillin binding of  $\beta$ -lactamases and carboxypeptidase, have suggested the C-terminal half of PBP 1A and 1B to be the functional domain of the transpeptidase activity and includes the penicillin binding domain. In  
30 addition, Nakagawa et al. (1987) showed that a truncated *ponB* gene encoding the N-terminal 478 amino acids of PBP 1B is capable of the transglycosylase reaction.



On the basis of these findings, it has been suggested that the high molecular weight PBP 1A and 1B proteins are two domain-proteins, with the N-terminal half forming the transglycosylase domain and the C-terminal half the transpeptidase domain. The two domains have been predicted by computer analysis to be joined by a linker or hinge region which does not structurally or enzymatically contribute to the function of the protein. The linker region of *E.coli* PBP 1B has been predicted to be from position 545-559 while that for *E.coli* PBP 1A around position 501.

The monofunctional truncated variants of PBP according to the invention will, when used in x-ray crystallography, facilitate obtaining structural information of the transglycosylase domain of penicillin binding proteins. In addition, the reduced size of the monofunctional variant will facilitate crystallization.

In a preferred form, a water-soluble polypeptide according to the invention has an amino acid sequence which is identical to SEQ ID NO: 2, 4, 6, 12 or 13 in the Sequence Listing.

The observation that deletion of the *ponA* and *ponB* genes is lethal (Yousif et al., 1985) does not address the question of essentiality of the transglycosylase activity of the encoded PBP 1A proteins, since the deletion results in the loss of both transglycosylation and transpeptidation activities. In addition, this experiment does not address the possibility that the transglycosylase enzyme activity can be contributed by a Penicillin Binding Protein other than PBP 1A or PBP 1B. It is also possible that hitherto undescribed Penicillin Binding Proteins and/or other proteins that contribute to the transglycosylase activity exist.

Alignment of the amino acids forming the putative transglycosylase domain of PBP 1A and 1B reveals three stretches of 9 out of 12 (Region 1), 9/10 (Region 2) and 8/10 (Region 3) amino acids identical within the

N-terminal half of these two proteins (Broome-Smith et al., 1985) (Fig. 14). The same 3 regions are identically conserved among two other recently described protein sequences; *Streptococcus pneumoniae* PBP 1A (Martin et al., 1992) and a 94 kDa protein from *Haemophilus influenzae* (Tomb et al., 1991). The conservation of these residues in such diverse species suggests their critical requirement in either maintaining structural aspects of the protein, or in the transglycosylation reaction itself.

The overlapping functional transglycosylase and transpeptidase activities of PBP 1A and 1B also suggests conservation of the catalytic centres and that molecules designed to interact with the catalytic centre of PBP 1A would be reactive also with PBP 1B.

The functional transglycosylase activity of the expressed protein can be studied either in a direct *in vitro* assay using appropriate substrates, or in an assay measuring the ability of the protein to complement the deletion of the corresponding genes in the chromosome. It has been shown that a plasmid with a gene encoding the wild type product (PBP 1A or PBP 1B) is capable of maintaining the viability of the *E.coli* cell (Yousif et al., 1985). This trans-complementation technique can be utilized to assess the functional nature of the mutant gene(s) encoding PBPs with mutations inactivating one of the enzymic (transglycosylation or transpeptidation) functions. The ability of such mutant products to complement in trans the deletion of the chromosomal *ponA* and *ponB* genes would define the essential requirement of the individual enzymic functions.

There is consequently a need for research tools which will make it possible to study the effects of specific inactivation of the transglycosylase activity of Penicillin Binding Proteins.

Consequently, a further aspect of the invention is a polypeptide which is a transglycosylase deficient derivative of a bacterial bifunctional penicillin

binding protein, said penicillin binding protein being bound to the cell membrane when expressed in a bacterial cell and being capable of exhibiting both transglycosylase and transpeptidase activities and said derivative lacking the capability to exhibit transglycosylase activity but retaining the capability to exhibit transpeptidase activity. The "bacterial cell" mentioned above is preferably an *Escherichia coli* cell.

The transglycosylase deficient PBP variants can advantageously be used in X-ray crystallography for the purpose of obtaining structural information of the activity sites of PBPs. Structural analysis of crystal form of soluble transglycosylase deficient PBP variants could allow delineation of the catalytic region and facilitate the design of molecules capable of specifically inhibiting the transglycosylase activity.

In a preferred form, the transglycosylase deficient polypeptide according to invention is a polypeptide which is lacking transglycosylase activity because of a mutation or deletion in the second conserved region of the gene coding for said polypeptide.

In a further preferred form, the transglycosylase deficient polypeptide according to the invention has an amino acid sequence which is identical to SEQ ID NO: 7, 8, 9, or 10 in the Sequence Listing.

The conventional purification procedure employed for the enrichment of penicillin binding proteins has been the use of a "penicillin" affinity. The binding of the protein to penicillin is covalent and requires harsh conditions to elute the bound protein. This may lead to a certain degree of inactivation of the enzymic activity of the protein. There is consequently a need for alternate affinity matrices for the efficient purification of the proteins.

Included in the invention is consequently a polypeptide comprising (a) a first polypeptide which is a PBP variant according to the invention; and (b) an additional polypeptide which allows binding to an affinity matrix; there being a cleavage site between said polypeptides.

5

The "additional polypeptide" mentioned above can preferably be glutathione-S-transferase or a polypeptide substantially similar to glutathione-S-transferase. Such an additional polypeptide will enable rapid purification of the protein using Glutathione Sepharose<sup>®</sup> affinity matrix. In another preferred form, the additional polypeptide is a polypeptide rich in histidine residues, which residues will confer on the protein the ability to bind to an Ni affinity column. The additional polypeptide can be fused either to the N-terminus or the C-terminus of the soluble/membrane bound PBP.

15

The ability of the fusion proteins to bind to an affinity matrix allows immobilisation of the protein. Such immobilised proteins can be used for analysis of competitive binding of different ligands to the bound active protein, and thus for screening of compounds binding to the enzymic domain of interest.

20

The polypeptides according to the invention are not to be limited strictly to any one of the sequences shown in the Sequence Listing. Rather the invention encompasses polypeptides carrying modifications like substitutions, small deletions, insertions or inversions, which polypeptides nevertheless have substantially the biochemical activities of the PBP variants which amino acid sequence is disclosed in the Sequence Listing. Included in the invention are consequently also polypeptides, the amino acid sequence of which is at least 90% homologous, preferably at least 95% homologous, with the amino acid sequence of any of the PBP variants according to the invention.

25

30

A further object of the invention is to provide isolated and purified DNA molecules which have nucleotide sequences coding for any one of the PBP variants according to the invention.

5 In a preferred form of the invention, the said DNA molecules have nucleotide sequences identical to SEQ ID NO: 1, 3 or 5 in the Sequence Listing. However, the DNA molecules according to the invention are not to be limited strictly to any of the sequences shown in the Sequence Listing. Rather the invention encompasses DNA molecules carrying modifications  
10 like substitutions, small deletions, insertions or inversions, which nevertheless encode proteins having substantially the biochemical activities of the PBP variants according to the invention.

Included in the invention is also a DNA molecule which nucleotide  
15 sequence is degenerate, because of the genetic code, to the said nucleotide sequence coding for a PBP variant according to the invention. The natural degeneracy of the genetic code is well known in the art. It will thus be appreciated that the DNA sequences shown in the Sequence Listing are only examples within a large but definite group of DNA sequences which  
20 will encode the PBP variants which amino acid sequences are shown in the Sequence Listing.

A further aspect of the invention is a replicable expression vector which carries and is capable of mediating the expression of a DNA molecule  
25 according to the invention. In the present context the term "replicable" means that the vector is able to replicate in a given type of host cell into which it has been introduced. Examples of vectors are viruses such as bacteriophages, cosmids, plasmids and other recombination vectors. Nucleic acid molecules are inserted into vector genomes by methods well  
30 known in the art. A vector according to the invention can preferably be one of the plasmids listed in Table 1 below.

Included in the invention is also a host cell harbouring a vector according to the invention. Such a host cell can be a prokaryotic cell, a unicellular eukaryotic cell or a cell derived from a multicellular organism. The host cell can thus e.g. be a bacterial, yeast or mammalian cell. The methods employed to effect introduction of the vector into the host cell are well-known to a person familiar with recombinant DNA methods.

A further aspect of the invention is a process for production of a polypeptide which is a derivative of penicillin binding protein, comprising growing a host cell according to the invention in or on a culture medium for expression of the polypeptide and optionally recovering the polypeptide. An appropriate host cell may be any of the cell types mentioned above, and the medium used to grow the cells may be any conventional medium suitable for the purpose.

The High Molecular Weight-Penicillin Binding Proteins have been shown to be anchored to the membrane, but the majority of the protein is within the periplasmic space of the cell (Edelman et al. 1987). Thus PBP derivatives, devoid of the membrane signal / anchoring sequences, are forced to fold into their native state in a heterologous environment, namely the cytosol. This often leads to misfolding, and the majority of the expressed protein aggregates into an inactive form referred to as inclusion bodies.

It has now surprisingly been found that high yields of an active water-soluble PBP variant can be obtained by regulated transcription of the gene encoding the said PBP variant. Such regulated transcription involves (i) using a suboptimal concentration of the inducer isopropyl thiogalactoside (IPTG); and (ii) culturing the cells expressing the PBP variant at reduced temperature. A cumulative effect of these factors contributes to the overall recovery of the active soluble protein. Consequently, lower rates of expression, achieved through the mentioned combination of (i) sub-optimal

de-repression of promoter systems and (ii) increased generation time by lowering of the temperature of cultivation, will enhance the solubility of proteins lacking the membrane anchoring segment.

5 A further important aspect of the invention is a process for the production of a water soluble polypeptide according to the invention which comprises culturing *Escherichia coli* cells harbouring an expression vector wherein a DNA coding sequence for said polypeptide is under the control of an isopropyl thiogalactoside (IPTG) - inducible promoter, said culturing being  
10 carried out in the presence of a sub-optimal concentration of IPTG for induction of the said promoter and at a temperature in the range of 20 to 24°C, preferably 22°C. The concentration of IPTG can preferably be approximately 0.01 mM.

15 In the case of expression of *ponAdel23*, a gene encoding a PBP variant according to the invention, such regulated transcription by (i) controlled de-repression of the T7 promoter by using sub-optimal concentration of the inducer IPTG and (ii) reducing the growth rate by culturing at 22°C, resulted in yields of the active protein which reached nearly 50% of the  
20 total induced protein of interest. The growth and induction conditions were critical for the efficient recovery of the soluble protein, as growth at higher temperatures or induction with higher concentrations of IPTG resulted in the majority of the protein becoming inactive and forming inclusion bodies.

25 I will be appreciated that this method for controlled expression is applicable to other inducible promoter systems, e.g. the tac system, where the inducer is IPTG and the host is a lac Y negative host.

30 A route to obtain relevant structural information on the active site configuration of an enzyme is the production and characterisation of monoclonal antibodies capable of inhibiting the enzymic reaction. The

antibodies inhibiting the activity represent molecules which block or compete with the substrate for entry into the active site pocket, or can represent molecules which can prevent structural transitions required for catalytic activity. In either case, these antibodies can be used as a tool to quantitate interaction of the target enzyme with binding of radiolabelled inhibitory compounds to judge the affinity of interaction provided the affinity of the inhibiting antibody is known. A further use of mapping the epitopes recognised by the inhibitory antibodies is the ability to delineate residues forming the active site.

Consequently, a further aspect of the invention is a method of identifying an antibody capable of binding a bacterial bifunctional penicillin binding protein which includes the step of employing a polypeptide according to the invention in an antibody binding assay and selecting antibodies that bind to the polypeptide.

Also included in the invention are monoclonal antibodies directed to a PBP variant according to the invention. Such a monoclonal antibody is prepared using known hybridoma technology by fusing antibody-producing B-cells from immunized animals with myeloma cells and selecting the resulting hybridoma cell line producing the desired antibody.

Another aspect of the invention is a method of assaying for compounds which bind to a penicillin binding protein, said method comprising (a) contacting a polypeptide which is a PBP variant according to the invention with a compound to be investigated; and (b) detecting whether said compound binds to the said PBP variant.

For example, a method of assaying for compounds which bind to a penicillin binding protein can comprise (a) culturing host cells according to the invention; (b) lysing the said cells and isolating the crude cell extract; (c) exposing the said cell extract to potential inhibitors of a penicillin



binding protein; (d) introducing an agent, known to bind a penicillin binding protein, to the said cell extract; (e) removing the unbound fraction of said agent; and (f) assaying the presence of said agent remaining in the cell extract.

5

Another method of assaying for compounds which bind to a penicillin binding protein could comprise (a) exposing a polypeptide which is a PBP variant according to the invention, immobilised on a solid support, to a potential inhibitor of a penicillin binding protein; (b) exposing an agent,  
10 known to bind a penicillin binding protein, to the immobilised polypeptide; (c) removing the unbound fraction of said agent; and (d) assaying the presence of said agent bound to the immobilised polypeptide.

15

In a preferred form, the said method is a method of assaying for compounds which bind to the transglycosylase domain of a penicillin binding protein, said method comprising (a) exposing the transglycosylase domain of a polypeptide according to the invention, with the proviso that the polypeptide is not a transglycosylase deficient PBP variant, said polypeptide being immobilised on a solid support, to a potential inhibitor  
20 of the transglycosylase activity of a penicillin binding protein; (b) exposing an agent, known to bind the transglycosylase domain of a penicillin binding protein, to the immobilised polypeptide; (c) removing the unbound fraction of said agent; and (d) assaying the presence of said agent bound to the immobilised polypeptide.

25

Antibodies specific for transpeptidase can be immobilised on a BIAcore sensor chip surface. The BIAcore system, wherein "BIA" stands for "Biospecific Interaction Analysis", is available from Pharmacia Biosensor, Sweden. Protein binding to the immobilised antibody is detected by the  
30 output RU-signal. Screening for TP inhibitors will be possible by a competitive assay wherein soluble protein is preincubated with test compounds. Binding of a test compound to the protein will result in a

decrease in protein binding to TP specific antibody. In the same way, monoclonal antibodies specific for transglycosylase can be used in screening for TG inhibitors.

- 5 In a similar way, ampicillin or modified moenomycin can be coupled to the surface and used in an indirect competitive assay whereby protein is preincubated with test ligand prior to introduction in the BIAcore.

10 Consequently, yet another method of assaying for compounds which bind to a penicillin binding protein could comprise (a) exposing a polypeptide which is a PBP variant according to the invention to a potential inhibitor of a penicillin binding protein; (b) exposing the polypeptide to an agent, known to bind a penicillin binding protein, which agent is immobilised on a solid support; and (c) assaying the presence of polypeptide bound to the  
15 immobilised agent.

In a preferred form, the said method is a method of assaying for compounds which bind to the transglycosylase domain of a penicillin binding protein, said method comprising (a) exposing the transglycosylase  
20 domain of a polypeptide according to the invention, with the proviso that the polypeptide is not a transglycosylase deficient PBP variant, to a potential inhibitor of a penicillin binding protein; (b) exposing the said polypeptide to an agent, known to bind to the transglycosylase domain of a penicillin binding protein, which agent is immobilised on a solid support; and (c) assaying the presence of polypeptide bound to the immobilised  
25 agent.

The "agent known to bind a penicillin binding protein" referred to above can e.g. be a monoclonal antibody or a labelled antibiotic compound such  
30 as [<sup>3</sup>H]ampicillin.

A further aspect of the invention is a method of determining the protein structure of a penicillin binding protein, characterized in that a polypeptide which is a PBP variant according to the invention is utilized in X-ray crystallography.

5

Some of the features of the preferred PBP variants according to the invention are summarised in Table 1 below. The plasmids listed in the Table have been deposited under the Budapest Treaty at the National Collection of Industrial and Marine Bacteria Limited (NCIMB), Aberdeen,

10

Scotland, UK. The date of deposit is 28 June 1994.

TABLE 1

Example no.	Features	Plasmid (pARC)	Deposit no. (NCIMB)	Fig.	SEQ ID NO:
<b>Soluble variants</b>					
1.1	<i>E.coli</i> PBP 1A with aa 1-23 deleted	0558	40666	3	1, 2
2.1	<i>E.coli</i> PBP 1B with aa 65-87 deleted	0559	40667	9	3, 4
3.1	<i>S.pneumoniae</i> PBP 1A with aa 1-38 deleted	0512	40665	12	5, 6
<b>Transglycosylase deficient variants</b>					
4.1	<i>E.coli</i> PBP 1B with glutamines 270-271 substituted to alanines	0438	40661		7
	<i>E.coli</i> PBP 1B with glutamines 270-271 substituted to leucines	0468	40662		8
	<i>E.coli</i> PBP 1B with aa 264-271 deleted	0469	40663		9
4.2	<i>E.coli</i> PBP 1A with glutamines 123-124 substituted to alanines	0571	40668	19	10
<b>Truncated variants</b>					
5.1	aa 1-553 of <i>E.coli</i> PBP 1B	0592	40669	21	11
	aa 1-553 of <i>E.coli</i> PBP 1B, with aa 65-87 deleted	0593	40670	22	12
5.2	aa 210-368 of <i>E.coli</i> PBP 1B	0392	40659	23	13
<b>Fusion proteins</b>					
6.1	<i>E.coli</i> PBP 1A with 23 aa deletion, ligated to glutathione-S-transferase	0499	40664	24	
6.2	<i>E.coli</i> PBP 1A with 23 aa deletion, ligated to histidine stretch	0400	40660	25	

## EXAMPLES OF THE INVENTION

5 In the following examples, the terms "standard protocols" and "standard procedures" are to be understood as protocols and procedures found in an ordinary laboratory manual such as the one by Sambrook, Fritsch and Maniatis (1989).

### EXAMPLE 1

10

#### 1.1. Construction of gene encoding soluble form of *E.coli* PBP 1A

15 The possible amino acid residues involved in the membrane anchoring region of PBP 1A was deduced following the computer program described by Kyte & Dolittle (1982). The predicted hydrophobicity of the N-terminal 60 amino acid is shown in Fig. 1. Based on this hydrophobicity profile, it was predicted that the N-terminal 23 amino acids were strongly implicated to contribute to the membrane anchoring domain of the protein, but may not entirely encompass the membrane anchoring domain. This region was  
20 then putatively designated as the region involved in "membrane anchoring".

25 The plasmid pBS98, harbouring the native *ponA* gene (encoding wild type PBP 1A), was obtained from Prof. B.S. Spratt, Microbial Genetics Group, School of Biological Sciences, University of Sussex, Brighton, UK. The construction of pBS98 is described in Broome-Smith et al. (1985). Plasmid DNA from cells harbouring pBS98 was made following standard protocols.

30 Oligonucleotide primers for use in the polymerase chain reaction (PCR) were synthesized in Applied Biosystems Model 380 A. The 5'-oligonucleotide primer used was TG-82:

NcoI  
↓

5' -TCG ACC ATG GGC CTA TAC CGC TAC ATC G-3'

M G L Y R Y I

5                                    23   24   25   26   27   28   29 (Amino acid No.)

TG-82 incorporates the following characteristics: (1) it allows construction of mutant *ponA* gene whose encoded product would have the 24th amino acid (glycine) of the wild type PBP 1A as the second amino acid of the expressed mutant protein; and (2) it introduces DNA sequences recognized by the restriction enzyme *NcoI*. This introduces the codon ATG which corresponds to the first amino acid of the mutant PBP 1A when expressed in suitable systems.

The 3'-oligonucleotide primer used was TG-64:

5' -CGC GGA TCC GAA TCA CAA CAA TTC CTG TGC-3'

↑  
*BamHI*

TG-64 has the following characteristics: (1) it introduces a termination codon following the 850th amino acid of the structural protein of PBP 1A; (2) it introduces a site for the restriction enzyme *BamHI* to facilitate cloning into suitable expression vectors.

Using these primers, PCR was carried out using pBS98 DNA as template following standard protocols. A DNA fragment of approximately 2.5 kb was amplified. The fragment was digested with the restriction enzyme *NcoI* followed by digestion with *BamHI*. This 2.5 kb *NcoI* - *BamHI* DNA fragment was then ligated to the vector pBR329 (Covarrubias et al., 1982) previously cut with *NcoI* and *BamHI*. Ligation of the two DNA fragments were carried out using standard protocols and the ligation mixture transformed into *E.coli* DH 5α. The transformed cells were plated on LB agar plates with 50 µg/ml ampicillin. Following overnight incubation at 37°C, individual ampicillin resistant colonies were tested for their

tetracycline sensitivity as insertion into the *NcoI* - *BamHI* region renders the plasmid chloramphenicol and tetracycline sensitive. A recombinant plasmid bearing the 2.5 kb insert was designated pARC0488.

5 The *NcoI* - *BamHI* 2.5 kb DNA fragment was released from pARC0488 and ligated to *NcoI* - *BamHI* cleaved and purified pARC038 (Fig. 2). The plasmid pARC038 is a derivative of pET11d (Studier et al., 1990) in which the *EcoRI* and *PstI* sites were made blunt ended with T4 exonuclease and  
10 the *EcoRI* - *PstI* 0.75 kb DNA fragment replaced with a blunt ended kanamycin resistance cartridge (Pharmacia Biochemicals). The ligation mixture was transformed into competent cells of *E.coli* BL 26 (DE3). The transformation mix was plated on LB agar with 50 µg/ml kanamycin. Mini-prep plasmid DNA was made from several kanamycin resistant colonies and screened by restriction endonuclease mapping using standard  
15 procedures.

One of the colonies harbouring plasmid with expected structure (Fig. 3) was labelled pARC0558 (NCIMB 40666). The DNA sequence of the mutant *ponA* gene labelled as *ponAdel23* is shown as SEQ ID NO: 1. The amino  
20 acid sequence of the soluble PBP 1Adel23 is shown as SEQ ID NO: 2.

### 1.2. Expression of *ponAdel23*

25 *E.coli* BL 26 (DE3) cells (obtained from Dr. J.J. Dunn, Biology Dept., Brookhaven National Lab., Long Island, NY, USA) harbouring pARC0558 were grown in LB with 50 µg/ml kanamycin till an O.D. at 600 nm of 0.6 and induced with 0.01 mM isopropyl thiogalactoside (IPTG) for 6 hours.

30 Following 6 hours of induction, cells were harvested and broken by passing through a French press. After centrifugation at low speed to remove unbroken cells and debris, the cytosolic (soluble) fraction was obtained by either of the following two methods: (1) following a procedure

described Page et al. (1982) in which the pellet, membrane and soluble proteins are separated by sucrose gradient centrifugation; or (2) by spinning the obtained supernatant at 200,000 x g for 90 minutes, whereafter the supernatant obtained is taken as the cytosolic / soluble protein  
5 fraction.

### 1.3. Penicillin binding of expressed PBP 1Adel23

The obtained cytosolic fraction was tested for the presence of mutant PBP  
10 1A by following the method of Rojo et al. (1984). This procedure involves using [<sup>125</sup>I]cephradine as the labelled penicillin as it is specific for PBP 1A. Mutant PBP 1Adel23 capable of binding the labelled cephradine could be demonstrated in the cytosolic fraction. Approximately 50% of the expressed mutant protein fractionated as a soluble protein, while the remaining 50%  
15 fractionated into the inclusion body and/or into the membrane associated fractions. Consequently, enhanced levels of active mutant PBP 1Adel23 were obtained since the cells were induced with sub-optimal concentration of IPTG and the since cultures were grown at 22°C. The penicillin binding profile of the soluble PBP 1Adel23 is shown in Fig. 4.

20

### 1.4. Purification of soluble PBP 1Adel23

The cell pellet of *E.coli* BL26 (DE3) / pARC0558 obtained following 6 hours of induction at 22°C was washed twice with buffer A (30 mM Tris-Cl, pH  
25 8.0; 10 mM EDTA; 10 µg/ml leupeptin; 10 µg/ml aprotinin; 5 mM DTT) and resuspended in the same buffer. The cell suspension was passed through a French press at 1200 psi. The lysate was spun at 10,000 rpm for 10 minutes and the obtained supernatant centrifuged at 200,000 x g for 45 minutes. The obtained supernatant was then adjusted to 30% saturation  
30 with ammonium sulphate. The mixture was centrifuged at 12,000 rpm for 10 min and the pellet resuspended in buffer A containing 1 M NaCl. The dissolved pellet was then treated with Cephradine-Affigel 10 matrix.



Cephadrine was conjugated to Affigel 10 following the instructions of the manufacturers (Biorad Laboratories, USA). The soluble PBP 1Adel23 containing fraction, dissolved in buffer A containing 1 M NaCl, was incubated 16 hrs at 4°C. with cephradine-affigel 10 beads. The beads were then washed with Buffer A containing 1 M NaCl until the absorbance at 280 nm was nearly zero. Elution of PBP 1Adel23 was monitored by assaying for penicillin binding activity in the wash. This activity was measured using [<sup>125</sup>I]cephradine prepared as described in Rojo et al. (1984). Bound PBP 1Adel23 was eluted from the beads using 1 M hydroxylamine (pH 8.5) at 25°C for 120 minutes. This fraction was concentrated by ultrafiltration using YM 30 filters (Amicon, USA) in Buffer A with 0.25 M NaCl. The ultrafiltration also resulted in the removal of hydroxylamine. The purified fraction containing >85% of the protein species corresponding to PBP 1Adel23 showed both penicillin binding and transglycosylase enzyme activities. The protein profile as seen by Coomassie Brilliant Blue staining and the [<sup>125</sup>I]cephradine / penicillin binding profile of the different fractions, obtained during the various stages of purification, are shown in Fig. 5. The N-terminal amino acid sequence of the soluble PBP 1Adel23 was confirmed by sequencing the purified protein.

#### 1.5. Transglycosylase activity of soluble PBP 1Adel23

The transglycosylase activity of the soluble PBP 1Adel23 protein was measured using essentially the method described by Ishino et al. (1980). The substrate for the detection of the enzymic activity were essentially prepared and purified following the protocols described by Heijenoort et al. (1992). The concentration dependent transglycosylase activity of PBP 1Adel23 measured as the amount of peptidoglycan formed was compared to the amounts of peptidoglycan formed by different concentrations of the membrane bound form of native PBP 1A. As seen in Fig. 6, the peptidoglycan polymerizing efficiency of the mutant soluble PBP 1Adel23

was nearly identical to the enzymic activity of the membrane bound form of the protein.

5 It has consequently been found that the elimination of the 23 amino acid residue stretch does not interfere with the ability of the protein to assume its native structure capable of both the enzymatic activities, i.e. the transglycosylase and the transpeptidase activities.

## EXAMPLE 2

10

### 2.1. Construction of gene encoding soluble form of *E.coli* PBP 1B

15 The *ponB* gene encoding PBP 1B was obtained on a plasmid pBS96 from Prof. B.S. Spratt, Microbial Genetics Group, School of Biological Sciences, University of Sussex, Brighton, UK. The construction of pBS96, as well as the nucleotide sequence of the wild-type *ponB* gene and the derived amino acid sequence, are described in Broome-Smith et al. (1985).

20 The hydropathy plot of the N-terminal approximately 150 amino acids as derived using the method of Kyte and Doolittle (1982) is shown in Fig. 7. Analysis of the hydropathicity plot indicated that the amino acids at positions 65 to 87 of the PBP 1B sequence contributed largely to the hydrophobicity of the N-terminus and can be putatively assigned to be the membrane anchoring domain of the protein. In addition,  $\beta$ -lactamase studies of Edelman et al. (1987) had indicated that amino acids C-terminal to amino acid position 87 were present in the periplasmic space of the *E.coli* cell and that amino acids N-terminal to position 65 of PBP 1B were within the cytoplasm of the cell.

30 The strategy employed to construct a mutant *ponB* gene encoding a soluble form of PBP 1B is shown in Fig. 8. Initially a DNA fragment of approximately 200 bp of the 5'-end of the *ponB* gene was amplified by

PCR, from the *ponB* gene on the plasmid pBS96 (Broome-Smith et al., 1985). The oligonucleotide primers used were 5'-primer TG-77 (5'-GAA AAA CCA TGG CCG GGA ATG ACC-3') which includes a *NcoI* restriction enzyme site which also coincides with the start ATG codon of the sequence, and 3'-primer TG-84 (5'-AAG TCG CGA GCC GCG TTT GCC AC-3') which includes a site for the restriction enzyme *NruI* and encodes for amino acids corresponding to position 64 of the PBP 1B sequence.

Step 1: The PCR amplified fragment following restriction with the enzymes *NcoI* and *NruI* was cloned into the *NcoI* - *NruI* sites of the cloning vector pBR 329 (Covarrubias et al., 1982). Ligation, transformation and screening were carried out using standard protocols and the recombinant plasmid with the expected structure labelled pARC0547 (Fig. 8) was obtained.

Another DNA fragment of approximately 1.2 kb was amplified by PCR using primer sequences corresponding to amino acid 87 to 480. This DNA fragment encodes the C-terminal half of the TG domain of PBP 1B. The primers used were 5'-primer TG-79 (5'-CGG ATA TCG ATC AAA AAA TTC GTA GCC G-3') which included the nucleotide sequence for the cleavage site for the restriction enzyme *EcoRV*, and 3'-primer TG-80 (5'-GCG GAT CCT TAG TCG ACG ACC ACA ATC GCA G-3') which included the sequence for *BamHI* cleavage.

Step 2: The PCR amplification of this fragment was done using the *ponB* gene on pBS96 (Broome-Smith et al., 1985) DNA as template. The amplified fragment was cloned into the *EcoRV* - *BamHI* sites of pBR 329 (Covarrubias et al., 1982) using standard protocols. The recombinant plasmid obtained was labelled pARC0534 (Fig. 8).

Step 3: The 200 bp *NcoI* - *NruI* fragment cloned in pARC0547 was excised as a *NcoI* - *NruI* fragment and cloned into *NcoI* - *EcoRV* cleaved pARC0534 to obtain pARC0551 (Fig. 8).

5 The mutant *ponB* gene on pARC0551 has DNA sequences coding for the N-terminal 64 amino acids of PBP 1B fused to the nucleotide sequences encoding the amino acids 88 to 480. A 1.3 kb *Pst*I - *Bam*HI DNA fragment of pBS96 was then ligated to *Pst*I - *Bam*HI cleaved pARC0551 and the ligation mixture transformed into *E.coli* DH5 $\alpha$  using standard procedures. Individual transformants were then screened and colonies harbouring recombinant plasmid with the expected structure identified. The plasmid was labelled pARC0552. A *Nco*I - *Bam*HI fragment from pARC0552 encompassing the entire mutant *ponB* gene was then excised and ligated to the T7 expression vector pARC038 to obtain pARC0559 (NCIMB 40667; Fig. 9).

15 The 3'-end of the cloned fragment of Step 1 has the nucleotide sequence TCG (partial *Nru*I site sequence) while the 5'-end of the fragment cloned in Step 2 has the sequence ATC (partial *Eco*RV cleavage sequence). The junction nucleotide sequence which is the outcome of the fusion of TCG and ATC results in the introduction of the codons for serine and isoleucine. Thus the mutant *ponB* gene encodes a PBP 1B with the amino acid sequence 1 to 64 corresponding to the wild type PBP 1B fused to the sequence 87 to 844. The two stretches are joined by the amino acids serine and isoleucine.

20 The nucleotide sequence of the mutant *ponB* gene is shown as SEQ ID NO: 3 and the derived amino acid as SEQ ID NO: 4.

25

## 2.2. Expression of soluble PBP 1B

30 The plasmid DNA of pARC0559 was transformed into the T7 expression host *E.coli* BL 26 (DE3) and the restriction map profile of the transformed plasmid confirmed using standard procedures. *E.coli* BL26 (DE3)/pARC0559 were grown at 22°C and induced with 0.01 mM IPTG and the cells allowed to grow for 6 hours. Cells were then harvested and

broken by passage through a french press. The lysate was centrifuged at 10,000 rpm for 10 minutes and the supernatant obtained was centrifuged at 200,000 x g for 45 minutes in a Beckman ultracentrifuge.

## 5 2.3. Characterization of the expressed soluble PBP 1B

The obtained supernatant, i.e. the cytosolic / soluble fraction, was tested for the presence of the mutant PBP 1B using [<sup>125</sup>I]ampicillin as the radio-ligand. The [<sup>125</sup>I]ampicillin was prepared as described by Rojo et al.  
10 (1984) for the preparation of [<sup>125</sup>I]cephradine. The mutant PBP 1B was detected in the soluble fraction and bound radioactive ampicillin.

Soluble PBP 1B could also be purified using Ampicillin - Affigel beads by a procedure analogous to the one described in Section 1.4. The protein  
15 profile of the different fractions seen by Coomassie Blue staining and the binding of [<sup>125</sup>I]ampicillin of the enriched PBP 1B fraction is shown in Fig. 10.

The purified protein was enzymatically active in the peptidoglycan  
20 transglycosylase assay (Heijenoort et al., 1992) and bound penicillin with an affinity comparable to that of the membrane bound native PBP 1B.

## EXAMPLE 3

### 25 3.1. Construction of gene encoding soluble form of *Streptococcus pneumoniae* PBP 1A

The molecular architecture of the *S.pneumoniae* PBP 1A is predicted to be similar to that of *E.coli* PBP 1A and PBP 1B protein in the fact that the  
30 protein is anchored to the membrane via a N-terminal membrane anchoring sequence. The nucleotide sequence of the gene encoding native membrane bound *S.pneumoniae* PBP 1A and its derived amino acid

sequence are described in Martin et al., (1992). The hydropathicity profile of the N-terminal 100 amino acids as derived by the Kyte and Doolittle plot is shown in Fig. 11. A stretch of 38 amino acids contributed significantly to the hydrophobicity of this region and was assumed to be the membrane interacting domain. A mutant gene of *S.pneumoniae* PBP 1A was constructed by deleting the nucleotide sequence coding for the N-terminal 38 amino acids of *S.pneumoniae* PBP 1A.

Using standard PCR protocols, sequences encoding the wild type *S.pneumoniae* PBP 1A gene was amplified as a 2.5 kb DNA fragment from the chromosome of *S.pneumoniae* strain PM1 (obtained from S.A. Lacks, Biology Department, Brookhaven National Laboratory, Upton, New York, USA) (Lacks, 1968) using the primers designed based on the sequence reported by Martin et al. (1992) and the amplified fragment cloned into the pneumococcal vector pLS 101 (Balganesh and Lacks, 1984).

The mutant gene encoding a soluble form of *S.pneumoniae* PBP 1A was constructed by using of plasmid DNA harbouring the wild type gene as template and amplifying a 2.3 kb DNA fragment by using PCR following standard procedures. The sequence of the primers used were 5'-primer TG-24 (5'-TAC GTT ACC ATG GCT CCT AGC CTA TCC-3') and 3'-primer TG-25 (5'-GAC AGG ATC CTG AGA AGA TGT CTT CTC A-3').

The 5'-primer TG-24 includes the sequence for the restriction enzyme *NcoI* while the 3'-primer TG-25 includes the site for the restriction enzyme *BamHI*. The *NcoI* and *BamHI* digested PCR amplified DNA fragment was ligated to *NcoI* - *BamHI* cleaved pARC039. The plasmid pARC039 is a derivative of pET 8c (Studier et al., 1990) in which the gene coding for the  $\beta$ -lactamase has been replaced by a kanamycin resistance cartridge.

Following ligation and screening using standard protocols, the structure of the recombinant plasmid was confirmed by detailed restriction mapping

and transformed into the T7 expression host *E.coli* BL 21 (DE3) (Studier et al., 1990). The recombinant plasmid was labelled pARC0512 (NCIMB 40665) and is schematically represented in Fig. 12.

- 5 The nucleotide sequence of the mutant *S.pneumoniae* PBP 1A gene is shown as SEQ ID NO: 5 and the derived amino acid sequence is shown as SEQ ID NO: 6.

10 3.2. Expression and characterization of soluble form of *Streptococcus pneumoniae* PBP 1A

The gene coding for soluble *S.pneumoniae* PBP 1A was expressed by a procedure analogous to the one described in Section 1.2.. The cytosolic fraction of *E.coli* BL 21 (DE3)/pARC0512 was isolated and tested for the  
15 presence of the soluble form of the *S.pneumoniae* PBP 1A<sub>del38</sub>. The radioactive ligand used for the binding studies was [<sup>3</sup>H]benzyl penicillin (Amersham) which was prepared as described earlier. Approximately 50% of the expressed protein from the mutant gene was found to be in the soluble fraction and bound [<sup>125</sup>I]penicillin (Rojo et al., 1984) or  
20 [<sup>3</sup>H]penicillin (Amersham) when the culture was grown and induced at 22°C with 0.01 mM IPTG. The growth and induction conditions were critical for the efficient recovery of the soluble protein, as growth at higher temperatures or induction with higher concentrations of IPTG resulted in the majority of the protein becoming inactive and forming inclusion bodies.  
25 Optimum levels of soluble active protein was found following induction for 6-8 h. (Fig. 13).

The soluble *S.pneumoniae* PBP 1A<sub>del38</sub> protein could also be efficiently purified essentially following the protocol used for the purification of the  
30 soluble *E.coli* PBP 1B protein.

The efficiency of penicillin binding of the soluble PBP 1Adel38 was comparable to that of the native membrane bound *S.pneumoniae* PBP 1A.

#### EXAMPLE 4

5

##### 4.1. Transglycosylase deficient *E.coli* PBP 1B

The conserved amino acids within Region 2 (Fig. 14) were chosen for site-directed mutagenesis. Within this stretch of 10 amino acids three different mutations were constructed:

- (a) the glutamines at position 270 and 271 of the PBP 1B sequence were changed to alanines;
- (b) the glutamines at position 270 and 271 of the PBP 1B sequence were changed to leucines; and
- (c) a deletion of the nucleotide sequence encoding amino acids from position 264 to 271.

Mutants of the *ponB* gene were constructed essentially following the procedure of Kunkel et al. (1985). A 1.5 kb *EcoRI* - *Sall* fragment of the *ponB* gene of the plasmid pBS96 was excised and cloned into *EcoRI* - *Sall* cleaved M13mp19 following standard protocols.

- (a) The primer used for mutating the nucleotide sequence coding for glutamine residues 270 and 271 into a sequence coding for alanine residues was TG-21:

5'-ACG CTG ACG GCC GCT CTG GTG AAA-3'  
T L T A A L V K

- (b) The primer used for mutating the sequence coding for the glutamine residues 270 and 271 into leucine residues was TG-23:

5'-ACG CTG ACG CTA TTG CTG GTG AAA-3'



T L T L L L V K

(c) The primer used for creating a deletion of the nucleotides encoding amino acids at position 264 to 271, all of which are within the conserved Region 2, was TG-22:

5' -CGC ACG GTA CAG CTG GTG AAA AAC-3'

R T V Q L V K

260 261 262 263 272 273 274 (amino acid no.)

10

Following mutagenesis, the nucleotide sequence of the mutagenized *EcoRI* - *Sall* fragment was determined following the protocol of Sanger et al. (1977). The sequencing confirmed the nucleotide changes and also ruled out any extraneous changes. This mutated 1.5 kb DNA fragment was ligated back to *EcoRI* - *Sall* cleaved pBS96 and the ligated DNA transformed in to *E.coli* DH5 $\alpha$  cells following standard protocols. Kanamycin resistant transformants were analyzed for their plasmid profiles and the plasmid with the TG-21 mutation (a) was labelled pARC0438 (NCIMB 40661). The mutant protein is referred to as PBP 1B QQ-AA (SEQ ID NO: 7).

20

The plasmid with the mutation (b) introduced by TG-23 was labelled pARC0468 (NCIMB 40662). The mutant protein is referred to as PBP 1B QQ-LL (SEQ ID NO: 8)

25

The plasmid with the deletion (c) obtained using TG-22 was labelled pARC0469 (NCIMB 40663). The mutant protein is referred to as PBP 1Bdel8 (SEQ ID NO: 9).

The four plasmid DNAs of pBS96, pARC0438, pARC0468 and pARC0469 were individually transformed into *E.coli ponB:spc<sup>r</sup>* cells (Broome-Smith et al., 1985) in which a deleted *ponB* gene had been marked with spectinomycin resistance marker.

30

5 *E.coli ponB:spc<sup>r</sup>* cells having the individual plasmids pBS96, pARC0438 or pARC0469 were grown and membrane preparations made following the procedure described by Spratt (1977) and the profile of the penicillin binding proteins analyzed on a 8% SDS-PAGE following labelling with radioactive penicillin. The mutant proteins were initially analyzed for *in vivo* stability and localization into the membrane using anti-PBP 1B sera raised against purified membrane bound native PBP 1B (Fig. 15).

10 The mutant proteins were found to be localized to the membrane and no degraded protein fragments reacting with the antibody could be detected indicating no gross instability. In addition the mutant proteins bound penicillin with an affinity comparable to that of the wild type PBP 1B (Fig. 15).

15 After assaying for transglycosylase activity as described in Heijenoort et al. (1978), no activity could be detected in the membranes expressing the mutant proteins, while the membrane with the wild type PBP 1B showed transglycosylase activity. This defines the amino acids 263 to 271 as being critical for transglycosylase activity.

20 The ability of the mutant proteins to bind penicillin with an affinity comparable to that of the wild type suggests that the transpeptidase activity of the mutant proteins would also be comparable to that of the wild type. Knowing that the bifunctional protein PBP 1B expressed on a plasmid can in trans complement the deletions of both *ponA* and *ponB* (Yousif et al., 1985) the ability of the transglycosylase negative / transpeptidase positive proteins PBP 1B QQ-AA and PBP 1Bdel8 to complement the absence of chromosomally encoded PBP 1A and 1B was tested.

30 The wild type *ponB* and the mutant *ponB* genes were cloned into low copy vector pMAK 705 (Hamilton et al., 1989). The resulting plasmids were

designated pARC0462, (wild type *ponB*, Fig. 16), pARC0463, (*ponB*del8, Fig. 17) and pARC0470 (*ponB* QQ-AA, Fig. 18). The plasmids were individually transformed into *E.coli* del *ponA* (*E.coli* with a deletion of the *ponA* gene).

5 *E.coli* del *ponA* /pARC0462, *E.coli* del *ponA*/pARC0463 and *E.coli* del *ponA*/pARC0470 were used as recipients of the P1 phage for the transduction of the *ponB*:*spc*<sup>r</sup> marker. The transduction was performed as described by Miller (1972) The phage P1 lysate was made on *E.coli*  
10 *ponB*:*spc*<sup>r</sup> strain (Yousif et al., 1985). Following infection, the infected cells were plated on spectinomycin. Integration of the DNA fragment harbouring and *ponB*:*spc*<sup>r</sup> transduced into any of the recipients results in the inactivation of the chromosomal *ponB* gene rendering the chromosome *ponA*<sup>-</sup> and *ponB*<sup>-</sup>. This genotype being lethal for the cell, the *E.coli*  
15 spectinomycin resistant transductants can remain viable only if the plasmid encoded *ponB* or the *ponB* mutant can functionally complement in trans.

The following *E.coli* strains were subject to phage P1 transduction analysis of trans-complementation: (1) *E.coli* AMA 1004 which has chromosomally coded wild type *ponA* and *ponB*; (2) *E.coli* AMA 1004 which has a  
20 chromosomally inactivated *ponB* and is the host for the plasmid coded mutant *ponB* genes; (3) *E.coli* AMA 1004 host bearing the plasmid pARC0462 encoding the wild type *ponB* gene; (4) *E.coli* AMA 1004 host bearing the plasmid pARC0463 encoding PBP 1Bdel8; and (5) *E.coli* AMA 1004 host bearing the plasmid pARC0470 encoding PBP 1B QQ-AA.

25

**Results**(Number of Km<sup>r</sup> transductants /-ml)

	(1) <i>E.coli</i> AMA 1004	3.0 × 10 <sup>4</sup>
	(2) <i>E.coli</i> AMA 1004, <i>ponB:spc<sup>r</sup></i>	< 1
5	(3) <i>E.coli</i> AMA 1004, <i>ponB:spc<sup>r</sup></i> (PBP 1B wt)	1.1 × 10 <sup>4</sup>
	(4) <i>E.coli</i> AMA 1004, <i>ponB:spc<sup>r</sup></i> (PBP 1Bdel8)	< 1
	(5) <i>E.coli</i> AMA 1004, <i>ponB:spc<sup>r</sup></i> (PBP 1B QQ-AA)	< 1

10 A comparable number of transductants were obtained for an internal marker : *trp* transduction using the same P1 phage lysate.

15 The above results show that viable transductants could be obtained only with wild type PBP 1B, indicating that the TG<sup>-</sup> TP<sup>+</sup> product encoded by *ponB* QQ-AA or *ponBdel8* could not complement the loss of chromosomally encoded PBP 1A and 1B. However, as these mutant proteins bind penicillin and thus can be assumed to have transpeptidase activity, the inability to complement must be the absence of the transglycosylase enzymic activity. These results confirm the essential nature of the transglycosylase activity of PBP 1A or 1B for the viability of the *E.coli* cell.

20 The mutants described define the Region 2 to be involved in the transglycosylase activity of the protein. As this stretch of amino acids is conserved within the four high molecular weight penicillin binding proteins namely *E.coli* PBP 1A, 1B and *S.pneumoniae* 1A and the 94 kDa protein of *H.influenzae* (Fig. 14) it is reasonable to assume similar catalytic or structural involvement of this region in all the transglycosylase enzymes utilizing substrates similar to that used by PBP 1A and 1B of *E.coli*.

#### 4.2. Transglycosylase deficient *E.coli* PBP 1A

30

The conserved Region 2 was chosen for site-directed mutagenesis and the nucleotide sequence coding for glutamine at positions 123 and 124 of *E.coli*

PBP 1A was changed to a sequence coding for alanine by PCR mutagenesis as follows. The 5' half of the *ponA* gene was amplified as 2 fragments, the 5'-fragment corresponding to amino acid 1 to 123 (fragment A) and the 3'-fragment corresponding to amino acid 124 to 434 (fragment B).

5

The sequence of the 5'-primer used for the amplification of fragment A was TG-93 (5'- GCG CGG ACC ATG GTG AAG TTC GTA AAG TAT-3') while the 3'-primer used for the amplification of fragment A was TG-106 (5'-CAG TGC TGC AGT AAT GGT ACT TGC CCC TTG-3').

10

The 3'-primer for fragment A amplification included the sequence for the restriction enzyme *PstI* which allowed the conversion of the sequence encoding the glutamine residues in position 123 and 124 into a nucleotide sequence coding for alanine residues.

15

Fragment B was amplified with the 5'-primer TG-107 (5'-ATT ACT GCA GCA CTG GCG AGA AAC TTC TTC-3') and the 3'-primer TG-108 (5'-TCG CGA GAT ATC TGG CGG ATT GAT CGA CAC-3').

20

The 5'-primer for amplifying fragment B included the sequence for the restriction enzyme *PstI* overlapping the sequence with that of 3'-primer for amplifying fragment A. Ligation of the 3'-end of fragment A to the 5'-end of fragment B recreated the site for *PstI* and resulted in the change of the nucleotide sequence encoding glutamine 123 and 124 into alanine 123 and 124. The amplified fragments A and B were individually cloned into pBR 329, and corresponding clones pARC0565 and pARC0566 were obtained.

25

30

Fragment A and B obtained from pARC0565 and pARC0566 were ligated to obtain pARC0567. The *ponA* sequences were completed by introducing an *XhoI* - *BamHI* fragment of pARC0489 (which is identical to pARC0558 (Fig. 3) except for having additional *LacI* and *Lac* operator sequences) into pARC0567 to obtain pARC0568. The *MluI* - *BglII* fragment of pARC0568

which included the Q<sub>123</sub> - Q<sub>124</sub> to A<sub>123</sub> - A<sub>124</sub> mutated region was then used to replace the otherwise identical MluI - BglII fragment of pBS98 to obtain the plasmid pARC0571 (Fig. 19; NCIMB 40668). The mutant protein was labelled PBP 1A QQ-AA (SEQ ID NO: 10).

5

Expression studies on the mutant indicated that the mutant protein was localised to the membrane (as detected by anti PBP 1A antibodies) and bound penicillin with an affinity comparable to that of the native PBP 1A (Fig. 20).

10

An *in vivo* complementation assay, similar to that described in the previous section, was performed by checking the ability of mutant PBP 1A protein to complement in trans. The *in vivo* complementation was performed using phage P1 transduction and transducing *ponB*:*spc*<sup>r</sup> into the host *E.coli* (recipient) del *ponA* harbouring the plasmid encoding the mutant protein PBP 1A QQ-AA.

15

In order to carry out the complementation analysis the wild type *ponA* gene was cloned into the low copy vector pMAK 705 (Hamilton et al, 1989) to obtain pARC0583 and the mutant *ponA* gene encoding PBP 1A QQ-AA cloned into pMAK 705 to obtain pARC0582.

20

The following *E.coli* strains were subject to phage P1 transduction analysis of trans-complementation: (1) *E.coli* AMA 1004 which has chromosomally coded *ponA* and *ponB*; (2) *E.coli* AMA 1004 *ponA* which has a chromosomally inactivated *ponA* and is the host for the plasmid coded mutant *ponA* genes; (3) host bearing the plasmid pARC0583 encoding the wild type *ponA* gene; (4) host bearing the plasmid pARC0582 encoding PBP 1A QQ-AA.

25

30

**Results**(Number of  $\text{Sp}^{\text{r}}$  transductants / ml)

	(1) <i>E.coli</i> AMA 1004	$2.1 \times 10^3$
	(2) <i>E.coli</i> AMA 1004, <i>ponA</i>	< 1
5	(3) <i>E.coli</i> AMA 1004, <i>ponA</i> (PBP 1A wt)	$1.64 \times 10^3$
	(4) <i>E.coli</i> AMA 1004, <i>ponA</i> (PBP 1A QQ-AA)	< 1

A comparable number of transductants were obtained for internal marker : *trp* transduction using the same P1 phage lysate.

10

As shown above, no viable transductants could be obtained with *E.coli* del *ponA* / pARC0582 as recipient indicating that the mutant PBP 1A QQ-AA could not complement the absence of chromosomally encoded PBP 1A/1B. This indicates that the Q<sub>123</sub> and Q<sub>124</sub> of region 2 of PBP 1A also affects transglycosylase activity of the protein as the loss of the complementing function must be a reflection of the loss of transglycosylase activity. The transpeptidase activity of the protein is unaffected as tested by its affinity to bind penicillin.

15

20 These results argue in favour of the region 2 as a critical stretch of amino acids involved in the transglycosylase enzymic function and may be the explanation for the strong evolutionary conservation of this stretch of amino acids.

25 **EXAMPLE 5**5.1. Truncated *E.coli* PBP 1B

30

A mutant gene encoding the truncated PBP 1B consisting of the N-terminal 553 amino acids was constructed by PCR amplification using the 5'-primer TG-77 (5'-GAA AAA CCA TGG CCG GGA ATG ACC-3') and the 3'-primer TG-116 (5'- ATG GGA TCC TTA ATC ATT CTG CGG TGA-3').

The 5' end of the primer corresponded to the amino acid 553 in the wild type followed by the stop codon and a site for the restriction enzyme *Bam*HI. A fragment of 1.7 kb was amplified using pBS96 DNA as template. The PCR amplified fragment was cut with *Pst*I and *Bam*HI and cloned into  
5 *Pst*I-*Bam*HI restricted pARC0555 (pARC0555 has the full length *ponB* gene cloned as *Nco*I-*Bam*HI fragment into the expression vector pET11d. The *Nco*I site includes the initiation codon ATG) to obtain pARC0592 (NCIMB 40669; Fig. 21) The expressed protein (SEQ ID NO: 11) was shown to have transglycosylase activity, thus confirming the functional independence of  
10 this domain.

The soluble truncated PBP 1B, i.e. PBP 1B with N-terminal 553 amino acids but lacking the membrane anchoring hydrophobic domain from 65-87, was constructed by replacing the *Pst*I-*Bam*HI fragment of pARC0559 (Fig. 9)  
15 with the *Pst*I-*Bam*HI fragment of pARC0592 to obtain pARC0593 (NCIMB 40670; Fig. 22). The mutant *ponB* gene encodes the soluble form of PBP 1B and the expressed protein (SEQ ID NO: 12) was found to have transglycosylase activity.

## 20 5.2. Minimum substrate binding domain of truncated *E.coli* PBP 1B

Detailed computer analysis of the anatomy of the presumptive TG domain (aa 1-553) of PBP 1B indicated that aa 210-368 were probably sufficient for the binding of the lipid linked substrate and the transglycosylase reaction.  
25 This stretch of amino acids includes the 3 conserved domains Region I, II and III. The mutant gene encoding the truncated protein stretch 210-368 was constructed as follows.

A fragment of approx size 480 bp was amplified from pBS96 as substrate  
30 with the 5'-primer having the sequence TG-154 (5'-CAA TCC ATG GGT GAG CAG CGT CTG TTT G-3') were the initiation ATG codon is



immediately followed by the sequence encoding the 210th amino acid of PBP 1B.

5 The 3'-primer corresponded to the sequence TG-155 (5'-T CCA GAA TTC CAG TTT TGG GTT ACG-3') were the sequence encoded the amino acid 368 of PBP 1B followed by the nucleotide sequence that provides the restriction site for *EcoRI*, enabling fusion to sequences encoding an enterokinase site and a histidine stretch, which allows rapid purification of the protein on an Ni affinity column (cf. section 6.2 below).

10

A *NcoI-EcoRI* fragment was cloned into the plasmid pARC0400 that was restricted with *NcoI-EcoRI* to obtain the recombinant plasmid pARC0392 (NCIMB 40659; Fig. 23). The recombinant plasmid was transformed into *E.coli* BL26 (DE3) and a protein of approximately 17 kDa was detected largely in the soluble fraction after induction with IPTG.

15

Along similar lines the minimum substrate binding region of PBP 1A could be predicted to involve the stretch 62-220 in the wild type protein.

Production of this protein as a fusion with a histidine stretch allows high efficiency affinity purification of the expressed product using the Ni<sup>2+</sup> column. That the results will be similar to that obtained with truncated PBP 1B can be anticipated.

20

## EXAMPLE 6

25

### 6.1. N-terminal fusion of soluble *E.coli* PBP 1A to glutathione-S-transferase

Fusion of the *ponAdel23* gene at its 5'-end in frame to sequences coding for glutathione-S-transferase was made as described in the following section.

30

The vector chosen for the fusion gene construction was pGEX-3X obtained from Pharmacia Biochemicals. In order to fuse the 5'-initiation ATG of

*ponAdel23* in frame with the gene encoding glutathione-S-transferase a *Bam*HI site was introduced using a PCR primer whose sequence included the sequence for the restriction enzyme *Eco*RI. The 5'-primer used was TG-115:

5

```

5'-TCG AGG ATC CCC ATG GGC CTA TAC CGC TAC ATC G-3'
      =====
      EcoRI   BamHI

```

10 The 3'-primer used was TG-106, described in Section 4.2. The PCR amplified DNA Fragment A was digested with *Bam*HI and *Pst*I and cloned into the *Bam*HI - *Pst*I sites of the standard cloning vector pUC8 to obtain pARC0496. This Fragment A includes the N-terminal 102 amino acids of the PBP 1Adel23 protein. A *Bam*HI - *Mlu*I (site present within the fragment

15 A) 270 bp fragment obtained from Fragment A, a 2.2 kb *Mlu*I - *Eco*RI fragment which includes the rest of the portion of the *ponA* gene obtained from pARC0490 (pARC0490 has the wild type *ponA* gene cloned into the *Xba*I - *Bam*HI sites of the low copy vector pWKS29 (Fu Wang et al., 1991) facilitating the 3'-end of the *ponA* del 23 gene to be excised as an *Eco*RI

20 fragment) and a *Eco*RI - *Bam*HI cleaved pGEX-3X were ligated together and transformed into competent *E.coli* cells. Individual transformants were screened for recombinant plasmid and the plasmid with the expected structure was designated pARC0499 (NCIMB 40664; Fig. 24). The encoded fusion product on pARC0499 has the glutathione-S-transferase sequences at

25 its C-terminus linked to PBP 1Adel23 sequences via a Factor Xa cleavage recognition sequence.

Following induction with 1 mM IPTG, a fusion protein of expected size was found to be induced. The protein bound penicillin and was active in

30 the transglycosylase assay. Following cell lysis by passing the suspension through a French press, the cell free supernatant fraction was prepared as detailed in Section 1.4. for the purification of PBP 1Adel23. The supernatant fraction was passed through a Glutathione Sepharose® matrix

(Pharmacia Biochemicals) and the bound GST-PBP 1Adel123 was eluted with glutathione. The eluted protein was found to be 80% homogeneous. Free glutathione was removed by dialysis and the GST-PBP 1Adel 23 was cleaved with factor Xa.

5

PBP 1Adel23 thus purified was found to be active in both penicillin binding and the transglycosylase reactions.

## 6.2. C-terminal fusion of soluble *E.coli* PBP 1A to histidine stretch

10

Fusion of the *ponAdel23* gene at its 3'-end in frame to sequences encoding a stretch of 6 histidines was made as described below.

15

In the first step the *ponAdel23* gene was amplified using pBS98 DNA as template using the 5'-primer TG-115 (5'-TCG AGG ATC CCC ATG GGC CTA TAC CGC TAC ATC G-3') and the 3'-primer TTG-121 (5'-GTT AGA ATT CGA ACA ATT CCT GTG-3').

20

The 3'-primer introduced an *EcoRI* site at the 3' end of the *ponAdel23* gene while eliminating the translation stop codon. The PCR amplified modified *ponAdel23* gene fragment was digested with *PstI* and *EcoRI* to release a 930 bp 5'-end fragment and ligated to *PstI*-*EcoRI* digested pBR 329 to obtain the recombinant plasmid pARC0467.

25

In the next step, a double stranded synthetic oligonucleotide with sequences encoding the six histidines and the DNA sequence coding for amino acids recognised as the enterokinase cleavage site was synthesised and ligated to the newly created *EcoRI* site at the 3'-end of the *ponAdel23* gene on pARC0467. The synthetic oligonucleotides used were TG-122:

30

EcoRI

5' -AA TTC GAC GAC GAC GAC AAG CAC CAC CAC CAC CAC CAC TGA TAA G-3'

5

ENTEROKINASE

HISTIDINES

and TG 123 (5'-GAT CCT TAT CAG TGG TGG TGG TGG TGG TGC TTG  
TCG TCG TCG TCG-3').

10 The plasmid pARC0467 was linearised with *EcoRI* and the synthetic double  
stranded oligonucleotide ligated. Following ligation a *PstI* - *BamHI*  
(Fragment A) was released from the ligation mixture and cloned into the  
*PstI* - *BamHI* sites of pARC0558 (Fig. 3), to obtain pARC0400 (NCIMB  
40660; Fig. 25). The mutant *ponAdel23* fusion gene thus encoded a protein  
15 with the PBP 1Adel23 sequence fused to the amino acid sequence Asp-  
Asp-Asp-Asp-Lys fused to His-His-His-His-His-His at its C-terminus. The  
Asp-Asp-Asp-Asp-Lys sequence is recognised by the protease enterokinase  
and cleaves following the lysine residue. The six histidine residues confer  
on the protein the ability to bind to the metal nickel.

20

The recombinant plasmid pARC0400 was transformed in *E.coli* BL26(DE3)  
cells and induced under culture and temperature conditions identical to  
those used for the purification of PBP 1Adel23. The cells were lysed by  
passing through a French press. The lysate was centrifuged at 10,000 rpm  
25 for 10 min. The supernatant obtained after low speed centrifugation was  
then spun at 200,000 x g for 45 min and the supernatant obtained  
represented the cytosolic fraction. This fraction contained the protein  
encoded by the fusion gene and the recombinant fusion protein was  
labelled PBP 1Adel23EH. This protein PBP 1Adel23EH bound  
30 [<sup>125</sup>I]cephradine and was also active in transglycosylase assay. The soluble  
fraction was passed through a Ni affinity column and bound protein eluted  
in batches with increasing concentrations of imidazole essentially following  
the procedure described in "The Qia Expressionist" obtained from QIAGEN  
Inc. 9259 Eton Avenue, Chateworth, CA 91311 USA. The majority of PBP

1Adel23EH eluted with 250 mM imidazole and was approximately 85%  
homogenous. It was the only cephradine binding protein eluted from the  
column. Thus the ability of fusion protein to bind to the Ni column can be  
easily exploited both for efficient purification and immobilisation of the  
5 active protein.

## EXAMPLE 7

### 10 7.1. Use of cell extracts for enzyme assays and in screening

The crude cell extract made according to Example 6 can be analyzed for  
the ability to bind penicillin by reacting with [<sup>3</sup>H]ampicillin prepared  
according to Hackenbeck (1983). To adapt the procedure to large-scale  
screening, a 96 well microtitre plate is used to contain the reactions and the  
15 assay is performed using a Beckman Biomek robot. Crude cell extract is  
mixed with [<sup>3</sup>H]ampicillin for 15 min at 37°C. The proteins in the  
reaction are precipitated with TCA and collected on a glass filter,  
unbound ampicillin is washed off and filters counted in a scintillation  
counter. Alternatively, autoradiography can be used to assay the degree of  
20 binding of ampicillin.

Based on the above method, a competitive assay can be used to assess the  
ability of test compounds to bind to the transpeptidase site of a PBP  
variant. In this assay, the test compound is exposed to the crude cell  
25 extract for 15 min prior to the addition of ampicillin. A positive result is  
indicated by a reduction in the amount of radioactivity present on the glass  
filter.

### 30 7.2. Use of soluble immobilised protein in screening

Protein containing a histidine peptide which has been purified as described  
can be used for screening for compounds which inhibit transpeptidase

activity or transglycosylase activity. The purified full length or truncated protein is immobilised onto agarose gel to which Ni(II) has been coupled. Aliquots of the beads containing immobilised protein are then transferred to the wells of a microtitre plate, test compounds are added to the plate and incubated before unbound test substance is washed free. Compounds which bind to the transpeptidase site of the bifunctional protein can be detected by adding [<sup>3</sup>H]ampicillin to the reaction vessel and continuing essentially as described above. Alternatively monoclonal antibodies known to bind to the transpeptidase region can be used. Compounds which bind to the transglycosylase site can be assessed in a competitive assay by the use of monoclonal antibodies which bind to the transglycosylase region of the protein.

#### EXAMPLE 8

##### 8.1. Production of monoclonal antibodies to PBP 1A

The protocol for the production of monoclonal antibodies (mAbs) was essentially that described in "Antibodies - a laboratory manual" (ed. Harlow David Lane, Cold Spring Harbor, USA). Purified membrane bound PBP 1A was used as the immunogen. Balb-C mice, 6-8 weeks old were immunised with 50 µg of purified native PBP 1A in Freund's Complete Adjuvant. A booster injection of 20 µg PBP 1A in incomplete Freund's adjuvant was given intraperitoneally. Two weeks later the presence of serum antibodies was checked by ELISA using PBP 1A as the coated antigen. Mice with circulating antibodies were immunised intraperitoneally daily for 4 days with 20 µg of PBP 1A in saline and the mice sacrificed for isolating splenocytes for generating fusions.

The myeloma cell line used in fusion experiments was Sp 2/0-Ag 14 and these cells were fused with splenocytes from immunised mice at a ratio of 10 : 1. Fusion was carried out using standard protocols and antibody

production from the clones was monitored by ELISA against PBP 1A when the cells were > 90% confluent.

5 72 high producing clones were expanded to 24 well plates and the secreted antibody characterised using the following screens: (1) ELISA against membrane bound form of PBP 1A; (2) ELISA against soluble form of PBP 1A del 23; (3) Dot blot analysis against membrane bound PBP 1A to eliminate monoclonals reacting with the detergent solubilised purified PBP 1A protein only due to changes in the configuration during purification; 10 and (4) ELISA against membrane bound form of PBP 1B.

Based on these screens, a panel of 5 secreting clones were selected and subcloned twice to ensure monoclonality. Ascites with these hybridoma clones were raised following standard procedures and IgG was purified 15 from these ascites fluids, using Protein G-Sepharose<sup>®</sup> affinity chromatography as recommended by the manufacturers of Protein G-Sepharose<sup>®</sup> (Pharmacia Biochemicals).

20 These purified antibodies react specifically with PBP 1A in both the membrane bound and the soluble forms in ELISA, Dot blots and in Western blotting. Clones were obtained by a cloning procedure employing 3 cells / well. To ensure the monoclonality these clones were subcloned into 96 well microtitre plates by limiting dilution at 1 cell / well. The wells receiving one cell were carefully confirmed under the microscope and 25 allowed to grow with macrophage feeder layers so as to obtain progeny from a single hybrid cell. Following sub-cloning the secretion of mAbs to PBP 1A was again assayed in ELISA using full length PBP 1A. Finally two clones from each parent hybridoma were selected and one of them was expanded as ascites in pristine primed Balb/c mice. All the five clones 30 adapted to grow in peritoneal cavities and produced ascitic mAbs.

The ascitic mAbs were titrated against purified PBP 1A in ELISA. All the ascitic mAbs had a titre of  $> 5 \times 10^5$  in ELISA and recognised full length protein in western immunoblots. The ascitic mAbs were purified by protein-G affinity columns.

5

The immunoglobulin isotype of mAbs was determined by mouse Ig - isotype by ELISA using a kit obtained from Sigma chemicals USA. Four of the monoclonals belonged to IgG1 and one belonged to IgG2a immunoglobulin isotype.

10

Further characterization of mAbs was done by using full length membrane bound PBP 1A/1B in western blots. In addition the transglycosylase (TG) and transpeptidase (TP) domain specificity of mAbs was determined by using various truncated forms of the membrane-bound N-terminal of PBP 1A, N-terminal of PBP 1B and C-terminal of PBP 1B in Western immunoblots. Various full length and truncated membrane bound PBPs were expressed and the prepared membrane fractions were resolved on a SDS-PAGE. The proteins were transferred onto nitrocellulose membranes and subjected to western blot analysis using polyclonal *E.coli* PBP 1A antibodies and monoclonal antibodies.

15

20

Assessment of the penicillin binding inhibitory potential of the mAbs was determined essentially following the protocol described by den Blaauwen et al. (1990). The protein-G affinity purified mAbs was preincubated with PBP 1A followed by addition of [ $^3\text{H}$ ]benzyl penicillin or [ $^{125}\text{I}$ ]cephradine. Two of the mAbs competitively inhibited binding of the radiolabelled penicillin to PBP 1A.

25

30

Monoclonal antibodies specific for the TG domain of PBP 1A have been obtained by screening the secreted antibody of the original hybridoma clones to react with the protein representing the N-terminal 434 amino acids of PBP 1A in western blots. Antibody from clone TG-2 reacted with



the N-terminal truncated 434 amino acid analogue of PBP 1A but also inhibited (>80% inhibition) the transglycosylase activity of PBP 1A. This indicates that the antibody recognises sequences in the protein which are involved in (a) binding of the substrate; (b) catalysing the enzymic action; or (c) altering conformation of the protein allosterically. In either of the three possibilities, identification of compounds competing for the binding of TG-2 to PBP 1A would represent molecules interacting with identical sequences on PBP 1A. Thus the competitive binding assay could be used as a screening assay for the identification of the TG inhibitory compound.

#### BRIEF DESCRIPTION OF THE DRAWINGS

##### Figure 1

Hydropathicity profile of *E.coli* PBP 1A. The figure shows in expanded form the hydropathicity pattern of the N-terminal 55 amino acids of PBP 1A.

##### Figure 2

Schematic representation of the T7 translation fusion expression vector pARC038.

— vector sequences

== gene conferring kanamycin resistance  $Km^r$ , gene encoding the lactose repressor ( $lac I_q$ ), origin of replication (ori), T7 lac operator promoter, T7 phage terminator.

The direction of transcription of the different genes are shown by arrows. Relevant restriction enzyme sites are shown. Numbers next to the restriction site represent the nucleotide position taking the nucleotide at the upper twelve o'clock-position as zero.

**Figure 3**

Schematic representation of the vector pARC0558 encoding soluble PBP 1Adel 23 of *E.coli*.

- vector sequences
- 5    = mutant gene encoding PBP 1Adel23, kanamycin resistance  $Km^r$ ,  
lactose repressor (lac Iq) and the origin of replication ori.

**Figure 4**

- Expression of soluble PBP 1Adel23. Panel A represents the autoradiogram of the [ $^{125}I$ ]cephradine binding profile of the uninduced and induced cultures of *E.coli* BL 26 (DE3) harbouring pARC0558. Panel B represents the Coomassie Brilliant Blue staining protein profile of the same uninduced and induced cells. Lane (1): uninduced cytosol fraction; (2): uninduced membrane fraction; (3): induced cytosol fraction; (4): induced membrane fraction; (M): molecular weight markers.
- 10
- 15

**Figure 5**

- SDS-PAGE pattern of purified PBP 1Adel23. Panel A: Coomassie blue staining. Panel B: [ $^{125}I$ ]cephradine binding protein profile. Lanes (1): *E.coli* BL 26(DE3)/pARC0558 cytosolic fraction (200,000g supernatant); (2): 30% Ammonium sulphate supernatant fraction; (3): 30% Ammonium sulphate pellet fraction; (4): Cephradine affigel breakthrough fraction; (5): Molecular weight markers; (6-8): Cephradine affigel eluate.
- 20

**Figure 6**

- Transglycosylase activity profile of wild type PBP 1A and mutant PBP 1Adel23 using purified proteins.
- (▲—▲) represents activity of soluble PBP 1Adel23;
- (●—●) represents activity of membrane bound PBP 1A solubilised with octyl- $\beta$ -glucoside. X-axis represents the concentration of the proteins used in  $\mu$ g. Y-axis represents the quantities of peptidoglycan formed.
- 25
- 30

**Figure 7**

Hydropathicity profile of *E.coli* PBP 1B. The figure represents the expanded hydropathicity profile of the N-terminal 150 amino acids of *E.coli* PBP 1B.

**Figure 8**

Schematic representation of the cloning of the soluble transglycosylase domain of *E.coli* PBP 1B.

— vector sequences

== sequences encoding *ponB* gene fragments and  $\beta$ -lactamase

The *NcoI-NruI* fragment encoding the N-terminal 64 amino acids of PBP 1B was cloned into the *NcoI-EcoRV* sites of pARC0534 to obtain the plasmid pARC0551. This recombinant plasmid harbours the gene encoding amino acid 1 to 480 of PBP 1B with internal deletion of amino acid 65 to 87.

**Figure 9**

Schematic representation of pARC0559 encoding soluble PBP 1B.

— vector sequences

== sequences of the mutant *ponB* gene encoding the soluble form of PBP 1B (solPBP 1B), lactose repressor ( $\text{lac I}_Q$ ), kanamycin resistance ( $\text{Km}^r$ ) and the origin of replication (*ori*).

Arrows represent direction of transcription of the genes.

**Figure 10**

Purification of soluble PBP 1B. Panel A: SDS-PAGE, Coomassie blue staining of the different fractions. Panel B: [ $^{125}\text{I}$ ]ampicillin binding profile of the same fractions. Lanes (1) and (2): Cytosolic fraction of *E.coli* BL 26(DE3)/pARC0559 induced cells; (3): Breakthrough fraction of Ampicillin-Affigel column; (4): Molecular weight markers; (5) and (6): Eluted fraction from the Ampicillin-Affigel column.

**Figure 11**

Hydropathicity profile of *S.pneumoniae* PBP 1A. The figure shows the expanded profile of the hydropathicity profile of the N-terminal 100 amino acids of *S.pneumoniae* PBP 1A.

5 *Figure 12*

Schematic representation of the plasmid pARC0512 encoding soluble form of *S.pneumoniae* PBP 1A.

— represents vector sequences

10 == represents sequences of the gene encoding soluble PBP 1A of *S.pneumoniae* (sPBP 1A), kanamycin resistance Km<sup>r</sup> and the origin of replication (ori).

*Figure 13*

15 Penicillin binding profile of soluble *S.pneumoniae* PBP 1A. Host: *E.coli* BL 21(DE3)/pARC0512. Panel A: Coomassie Blue staining. Panel B: *In vivo* labelling with [<sup>3</sup>H]benzyl penicillin followed by SDS-PAGE. Lanes (1) and (2): Cytosolic fraction of cells induced at 22°C for 2 h and 20 h respectively; (3): Cytosolic fraction of cells induced at 30°C for 2 h; (4): Cytosolic fraction of cells induced at 37°C for 2h; (5): Molecular weight markers.

20

*Figure 14*

Amino acid alignment of conserved regions of the transglycosylase domain of high molecular weight penicillin binding proteins. The figure compares the conserved residues of the Regions 1, 2 and 3 among E.1A (*E.coli* PBP 1A), E.1B (*E.coli* PBP 1B), S.1A (*S.pneumoniae* PBP 1A), and H.inf (*Haemophilus influenzae* PBP 1A). (\*) indicates identical amino acid residues.

25

*Figure 15*

30 Analysis of membrane protein of *E.coli* cells harbouring plasmids with genes encoding mutant PBP 1B. Panel A: [<sup>3</sup>H]benzyl penicillin binding profile. Panel B: Western blotting with anti-PBP 1B sera. Lanes (1): Molecular weight markers; (2): Membrane fraction of *E.coli* JM 101/pBS96

- cells; (3): Membrane fraction of *E.coli* 900521 *ponB*:Spc<sup>r</sup> cells (This host lacks chromosomal encoded PBP 1B); (4): Membrane fraction of *E.coli* 900521 *ponB*:spc./pARC0438 cells; (5): Membrane fraction of *E.coli* 900521 *ponB*:spc/pARC0469; (6): Membrane fraction of *E.coli* 900521 *ponB*:spc/pARC0468.

Figure 16

Schematic representation of plasmid pARC0462 encoding wild type PBP 1B:

- 10 ——— vector sequences  
===== sequences of the *ponB* gene, replication origin (ori), chloramphenicol acetyl transferase (cm<sup>r</sup>) and portions of the lac Z multiple cloning site.

15 Figure 17

Schematic representation of plasmid pARC0463 encoding mutant *ponB* gene.

- vector sequences  
===== sequences of mutant the *ponB* gene encoding PBP 1Bdel8 amino acids, replication origin (ori), chloramphenicol acetyl transferase (cm<sup>r</sup>) and portions of the lac Z multiple cloning site.

Figure 18

Schematic representation of plasmid pARC0470 encoding mutant *ponB* gene.

- 25 ——— vector sequences  
===== sequences of mutant the *ponB* gene encoding PBP 1B Q<sub>271-272</sub> - A<sub>271-272</sub>, replication origin (ori), chloramphenicol acetyl transferase (cm<sup>r</sup>) and portions of the lac Z multiple cloning site.

30

Figure 19

Schematic representation of pARC0571 harbouring mutant *ponA* gene.

- vector sequences  
== sequences of mutant *ponA* gene (PBP 1A QQ-AA), kanamycin resistance  $Km^r$  origin of replication (ori).

5 *Figure 20*

[<sup>125</sup>I]Penicillin binding protein profile of wild type and mutant *E.coli* PBP 1A. Lane (1): *E.coli* AMA 1004 *ponB:spc<sup>r</sup>*/pBS 98 (w.t. *ponA*); (2): *E.coli* BL21 (DE3) *ponB:spc<sup>r</sup>*/pARC0570 (w.t. *ponA*); (3): *E.coli* AMA 1004 del *ponA*/pARC0571 (QQ-AA *ponA*); (4): *E.coli* AMA 1004 del *ponA*/pBS 98 (w.t. *ponA*); (5): Molecular weight markers.

*Figure 21*

Schematic representation of plasmid pARC0592.

- vector sequences  
15 == sequences of truncated *ponB* gene encoding for the N-terminal 553 amino acids of PBP 1B (hinge 1B), kanamycin resistance ( $Km^r$ ) and origin of replication (ori)

*Figure 22*

20 Schematic representation of plasmid pARC0593.

- vector sequences  
== sequences of mutant truncated *ponB* gene encoding a soluble form of the truncated N-terminal 553 amino acids of PBP 1B (soluble hinge 1B), kanamycin resistance  $Km^r$  and origin of replication (ori).

25

*Figure 23*

Schematic representation of plasmid pARC0392.

- vector sequences  
30 == sequences of mutant gene encoding truncated fragment of PBP 1B protein, representing amino acids 210-368 sequences fused in frame at its 3'-end to sequences encoding a enterokinase site followed by a

stretch of 6 histidines, kanamycin resistance  $Km^r$  and origin of replication (ori).

*Figure 24*

5 Schematic representation of plasmid pARC0499.

—— vector sequences

===== sequences of mutant *ponAdel23* gene fused at its 5'-end in frame to sequences encoding Glutathione-S-transferase encoding sequences,  $\beta$ -lactamase  $amp^r$  and origin of replication (ori).

10

*Figure 25*

Schematic representation of plasmid pARC0400.

—— vector sequences

===== sequences of mutant *ponAdel23* sequences fused in frame at its 3'-end to sequences encoding a enterokinase site followed by a stretch of 6 histidines, kanamycin resistance  $Km^r$  and origin of replication (ori).

15

20 REFERENCES

Balganesh, T.S. and Lacks, S. (1984): Gene 29, 221-230

den Blaauwen, T. et al. (1990): J. Bact. 172, 63-70

25

Broome-Smith, J.K. et al. (1985): Eur. J. Biochem. 147, 437-446

Covarrubias, L. & Bolivar, F. (1982): Gene 17, 79-89

30

Edelman, A. et al. (1987): Molecular Microbiology 1, 101-106

Fu Wang, R. and Kushner, S.R. (1991): Gene 100, 195-199

35

Hackenbeck et al. (eds.) The target of penicillin. W. de Gruyter publications, Berlin/New York 1983.

Hamilton, C.A. et al. (1989): J. Bacteriol. 171, 4617-4622

Heijenoort, Y. van, et al. (1978): FEBS Letters 89, 141-144

- Heijenoort, Y. van, et al. (1992): J. Bacteriol. 174, 3549-3557
- Ishino, F. et al. (1980): Biochem. Biophys. Res. Comm. 97, 287- 293
- 5 Kunkel, T.A. (1985): Proc. Natl. Acad. Sci. U.S.A. 82, 488-492
- Kyte and Doolittle (1982): J. Mol. Biol. 157, 105-132
- Lacks, S.A. (1968): Genetics 60, 685-706
- 10 Martin, C. et al. (1992): J. Bacteriol. 174, 4517-4523
- Miller, J.H. (ed.) (1972): Experiments in Molecular Genetics. Cold Spring Harbor Publications.
- 15 Nakagawa, J.S. et al. (1984): J. Biol. Chem. 259, 13937-13946
- Page, W.J. et al. (1982): J. Bacteriol. 151, 237-242
- 20 Rojo et.al. (1984): J. Antibiotics. 37, 389-393
- Sambrook, J., Fritsch, E.F. and Maniatis, T. (1989): Molecular Cloning: A laboratory manual, 2nd Ed., Cold Spring Harbor Laboratory, Cold Spring Harbor, NY
- 25 Sanger, et al. (1977): Proc. Natl. Acad. Sci. U.S.A. 74, 5463-5467
- Spratt, B.S. (1977): Eur. J. Biochem. 72, 341-352
- 30 Studier, F.W. et al. (1990): Methods in Enzymology 185, 61-89
- Tomb, F. et al. (1991): Gene 104, 1-10
- Yousif, S.Y. et al. (1985): J. Gen. Microbiol. 131, 2839-2845



## SEQUENCE LISTING

## (1) GENERAL INFORMATION:

## (i) APPLICANT:

(A) NAME: ASTRA AKTIEBOLAG  
 (B) STREET: Kvarnbergagatan 16  
 (C) CITY: Sodertalje  
 (E) COUNTRY: Sweden  
 (F) POSTAL CODE (ZIP): S-151 85  
 (G) TELEPHONE: +46-8 553 260 00  
 (H) TELEFAX: +46-8 553 288 20  
 (I) TELEX: 19237 astra s

(ii) TITLE OF INVENTION: Novel Polypeptides

(iii) NUMBER OF SEQUENCES: 13

## (iv) COMPUTER READABLE FORM:

(A) MEDIUM TYPE: Floppy disk  
 (B) COMPUTER: IBM PC compatible  
 (C) OPERATING SYSTEM: PC-DOS/MS-DOS  
 (D) SOFTWARE: PatentIn Release #1.0, Version #1.25 (EPO)

## (2) INFORMATION FOR SEQ ID NO: 1:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 2487 base pairs  
 (B) TYPE: nucleic acid  
 (C) STRANDEDNESS: double  
 (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: DNA (genomic)

## (vi) ORIGINAL SOURCE:

(A) ORGANISM: Escherichia coli  
 (B) STRAIN: DH5 alpha

## (vii) IMMEDIATE SOURCE:

(A) LIBRARY: PCR cloning  
 (B) CLONE: pARC 0558 Soluble PBP 1A del 23

## (ix) FEATURE:

(A) NAME/KEY: CDS  
 (B) LOCATION: 1..2487

## (ix) FEATURE:

(A) NAME/KEY: mat\_peptide  
 (B) LOCATION: 1..2484

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 1:

ATG GGC CTA TAC CGC TAC ATC GAG CCA CAA CTG CCG GAT GTG GCG ACA	48
Met Gly Leu Tyr Arg Tyr Ile Glu Pro Gln Leu Pro Asp Val Ala Thr	
1 5 10 15	
TTA AAA GAT GTT CGC CTG CAA ATT CCG ATG CAG ATT TAC AGC GCC GAT	96
Leu Lys Asp Val Arg Leu Gln Ile Pro Met Gln Ile Tyr Ser Ala Asp	
20 25 30	
GGC GAG CTG ATT GCT CAA TAC GGT GAG AAA CGT CGT ATT CCG GTT ACG	144
Gly Glu Leu Ile Ala Gln Tyr Gly Glu Lys Arg Arg Ile Pro Val Thr	
35 40 45	

TTG	GAT	CAA	ATC	CCA	CCG	GAG	ATG	GTG	AAA	GCC	TTT	ATC	GCG	ACA	GAA	192
L u	Asp	Gln	Ile	Pro	Pro	Glu	M t	Val	Lys	Ala	Phe	Ile	Ala	Thr	Glu	
	50					55					60					
GAC	AGC	CGC	TTC	TAC	GAG	CAT	CAC	GGC	GTT	GAC	CCG	GTG	GGG	ATC	TTC	240
Asp	Ser	Arg	Phe	Tyr	Glu	His	His	Gly	Val	Asp	Pro	Val	Gly	Ile	Phe	
65					70					75					80	
CGT	GCA	GCA	AGC	GTG	GCG	CTG	TTC	TCC	GGT	CAC	GCG	TCA	CAA	GGG	GCA	288
Arg	Ala	Ala	Ser	Val	Ala	Leu	Phe	Ser	Gly	His	Ala	Ser	Gln	Gly	Ala	
				85					90					95		
AGT	ACC	ATT	ACC	CAG	CAG	CTG	GCG	AGA	AAC	TTC	TTC	CTC	AGT	CCA	GAA	336
Ser	Thr	Ile	Thr	Gln	Gln	Leu	Ala	Arg	Asn	Phe	Phe	Leu	Ser	Pro	Glu	
			100					105					110			
CGC	ACG	CTG	ATG	CGT	AAG	ATT	AAG	GAA	GTC	TTC	CTC	GCG	ATT	CGC	ATT	384
Arg	Thr	Leu	Met	Arg	Lys	Ile	Lys	Glu	Val	Phe	Leu	Ala	Ile	Arg	Ile	
		115					120					125				
GAA	CAG	CTG	CTG	ACG	AAA	GAC	GAG	ATC	CTC	GAG	CTT	TAT	CTG	AAC	AAG	432
Glu	Gln	Leu	Leu	Thr	Lys	Asp	Glu	Ile	Leu	Glu	Leu	Tyr	Leu	Asn	Lys	
	130					135					140					
ATT	TAC	CTT	GGT	TAC	CGC	GCC	TAT	GGT	GTC	GGT	GCT	GCG	GCA	CAA	GTC	480
Ile	Tyr	Leu	Gly	Tyr	Arg	Ala	Tyr	Gly	Val	Gly	Ala	Ala	Ala	Gln	Val	
145					150					155					160	
TAT	TTC	GGA	AAA	ACG	GTC	GAC	CAA	CTG	ACG	CTG	AAC	GAA	ATG	GCG	GTG	528
Tyr	Phe	Gly	Lys	Thr	Val	Asp	Gln	Leu	Thr	Leu	Asn	Glu	Met	Ala	Val	
				165					170					175		
ATA	GCC	GGG	CTG	CCG	AAA	GCG	CCT	TCC	ACC	TTC	AAC	CCG	CTC	TAC	TCG	576
Ile	Ala	Gly	Leu	Pro	Lys	Ala	Pro	Ser	Thr	Phe	Asn	Pro	Leu	Tyr	Ser	
			180					185					190			
ATG	GAT	CGT	GCC	GTC	GCG	CGG	CGT	AAC	GTC	GTG	CTG	TCG	CGG	ATG	CTG	624
Met	Asp	Arg	Ala	Val	Ala	Arg	Arg	Asn	Val	Val	Leu	Ser	Arg	Met	Leu	
		195					200					205				
GAT	GAA	GGG	TAT	ATC	ACC	CAA	CAA	CAG	TTC	GAT	CAG	ACA	GCG	ACT	GAG	672
Asp	Glu	Gly	Tyr	Ile	Thr	Gln	Gln	Gln	Phe	Asp	Gln	Thr	Arg	Thr	Glu	
	210					215					220					
GCG	ATT	AAC	GCT	AAC	TAT	CAC	GCG	CCG	GAG	ATT	GCT	TTC	TCT	GCG	CCG	720
Ala	Ile	Asn	Ala	Asn	Tyr	His	Ala	Pro	Glu	Ile	Ala	Phe	Ser	Ala	Pro	
225					230					235					240	
TAC	CTG	AGC	GAA	ATG	GTG	CGC	CAG	GAG	ATG	TAT	AAC	CGT	TAT	GGC	GAA	768
Tyr	Leu	Ser	Glu	Met	Val	Arg	Gln	Glu	Met	Tyr	Asn	Arg	Tyr	Gly	Glu	
				245					250					255		
AGT	GCC	TAT	GAA	GAC	GGT	TAT	CGC	ATT	TAC	ACC	ACC	ATC	ACC	GCG	AAA	816
Ser	Ala	Tyr	Glu	Asp	Gly	Tyr	Arg	Ile	Tyr	Thr	Thr	Ile	Thr	Arg	Lys	
			260					265					270			
GTG	CAG	CAG	GCC	GCG	CAG	CAG	GCG	GTA	CGT	AAT	AAC	GTG	CTG	GAC	TAC	864
Val	Gln	Gln	Ala	Ala	Gln	Gln	Ala	Val	Arg	Asn	Asn	Val	Leu	Asp	Tyr	
		275					280					285				
GAC	ATG	CGC	CAC	GGC	TAT	CGC	GGC	CCG	GCA	AAT	GTG	CTG	TGG	AAA	GTG	912
Asp	Met	Arg	His	Gly	Tyr	Arg	Gly	Pro	Ala	Asn	Val	Leu	Trp	Lys	Val	
	290					295					300					
GGC	GAG	TCG	GCG	TGG	GAT	AAC	AAC	AAG	ATT	ACC	GAT	ACG	CTG	AAG	GCG	960
Gly	Glu	Ser	Ala	Trp	Asp	Asn	Asn	Lys	Ile	Thr	Asp	Thr	Leu	Lys	Ala	
305					310					315					320	

CTG CCA ACC TAT GGT CCG CTG CTG CCT GCC GCA GTC ACC AGC GCC AAT Leu Pro Thr Tyr Gly Pro Leu Leu Pro Ala Ala Val Thr S r Ala Asn 325 330 335	1008
CCT CAG CAA GCG ACG GCG ATG CTG GCG GAC GGG TCG ACC GTC GCA TTG Pro Gln Gln Ala Thr Ala Met Leu Ala Asp Gly Ser Thr Val Ala L u 340 345 350	1056
AGT ATG GAA GGC GTT CGC TGG GCG CGT CCT TAC CGT TCG GAT ACT CAG Ser Met Glu Gly Val Arg Trp Ala Arg Pro Tyr Arg Ser Asp Thr Gln 355 360 365	1104
CAA GGA CCG ACG CCG CGT AAA GTG ACC GAT GTT CTG CAA ACG GGT CAG Gln Gly Pro Thr Pro Arg Lys Val Thr Asp Val Leu Gln Thr Gly Gln 370 375 380	1152
CAA ATC TGG GTT CGT CAG GTT GGC GAT GCA TGG TGG CTG GCA CAA GTG Gln Ile Trp Val Arg Gln Val Gly Asp Ala Trp Trp Leu Ala Gln Val 385 390 400	1200
CCG GAA GTG AAC TCG GCG CTG GTG TCG ATC AAT CCG CAA AAC GGT GCC Pro Glu Val Asn Ser Ala Leu Val Ser Ile Asn Pro Gln Asn Gly Ala 405 410 415	1248
GTT ATG GCG CTG GTC GGT GGC TTT GAT TTC AAT CAG AGC AAG TTT AAC Val Met Ala Leu Val Gly Gly Phe Asp Phe Asn Gln Ser Lys Phe Asn 420 425 430	1296
CGC GCC ACC CAG GCA CTG CGT CAG GTG GGT TCC AAC ATC AAA CCG TTC Arg Ala Thr Gln Ala Leu Arg Gln Val Gly Ser Asn Ile Lys Pro Phe 435 440 445	1344
CTC TAC ACC GCG GCG ATG GAT AAA GGT CTG ACG CTG GCA AGT ATG TTG Leu Tyr Thr Ala Ala Met Asp Lys Gly Leu Thr Leu Ala Ser Met Leu 450 455 460	1392
AAC GAT GTG CCA ATT TCT CGC TGG GAT GCA AGT GCC GGT TCT GAC TGG Asn Asp Val Pro Ile Ser Arg Trp Asp Ala Ser Ala Gly Ser Asp Trp 465 470 475 480	1440
CAG CCG AAG AAC TCA CCA CCG CAG TAT GCT GGT CCA ATT CGC TTA CGT Gln Pro Lys Asn Ser Pro Pro Gln Tyr Ala Gly Pro Ile Arg Leu Arg 485 490 495	1488
CAG GGG CTG GGT CAG TCG AAA AAC GTG GTG ATG GTA CGC GCA ATG CGG Gln Gly Leu Gly Gln Ser Lys Asn Val Val Met Val Arg Ala Met Arg 500 505 510	1536
GCG ATG GGC GTC GAC TAC GCT GCA GAA TAT CTG CAA CGC TTC GGC TTC Ala Met Gly Val Asp Tyr Ala Ala Glu Tyr Leu Gln Arg Phe Gly Phe 515 520 525	1584
CCG GCA CAA AAC ATT GTC CAC ACC GAA TCG CTG GCG CTG GGT TCA GCG Pro Ala Gln Asn Ile Val His Thr Glu Ser Leu Ala Leu Gly Ser Ala 530 535 540	1632
TCC TTC ACC CCA ATG CAG GTG GCG CGC GGC TAC GCG GTC ATG GCG AAC Ser Phe Thr Pro Met Gln Val Ala Arg Gly Tyr Ala Val Met Ala Asn 545 550 555 560	1680
GGC GGC TTC CTG GTG GAC CCG TGG TTT ATC AGC AAA ATT GAA AAC GAT Gly Gly Phe Leu Val Asp Pro Trp Phe Ile Ser Lys Ile Glu Asn Asp 565 570 575	1728
CAG GGC GGC GTG ATT TTC GAA GCG AAA CCG AAA GTA GCC TGC CCG GAA Gln Gly Gly Val Il Phe Glu Ala Lys Pro Lys Val Ala Cys Pro Glu 580 585 590	1776

TGC	GAT	ATT	CCG	GTG	ATT	TAC	GGT	GAT	ACG	CAG	AAA	TCG	AAC	GTG	CTG	1824
Cys	Asp	Ile	Pro	Val	Ile	Tyr	Gly	Asp	Thr	Gln	Lys	Ser	Asn	Val	Leu	
	595						600					605				
GAA	AAT	AAC	GAT	GTT	GAA	GAT	GTC	GCT	ATC	TCC	CGC	GAG	CAG	CAG	AAT	1872
Glu	Asn	Asn	Asp	Val	Glu	Asp	Val	Ala	Ile	Ser	Arg	Glu	Gln	Gln	Asn	
	610					615					620					
GTT	TCT	GTA	CCA	ATG	CCG	CAG	CTG	GAG	CAG	GCA	AAT	CAG	GCG	TTA	GTG	1920
Val	Ser	Val	Pro	Met	Pro	Gln	Leu	Glu	Gln	Ala	Asn	Gln	Ala	Leu	Val	
625					630					635					640	
GCG	AAG	ACT	GGC	GCG	CAG	GAG	TAC	GCA	CCG	CAC	GTC	ATC	AAC	ACT	CCG	1968
Ala	Lys	Thr	Gly	Ala	Gln	Glu	Tyr	Ala	Pro	His	Val	Ile	Asn	Thr	Pro	
			645						650					655		
CTG	GCA	TTC	CTG	ATT	AAG	AGT	GCT	TTG	AAC	ACC	AAT	ATC	TTT	GGT	GAG	2016
Leu	Ala	Phe	Leu	Ile	Lys	Ser	Ala	Leu	Asn	Thr	Asn	Ile	Phe	Gly	Glu	
			660					665					670			
CCA	GGC	TGG	CAG	GGT	ACT	GGC	TGG	CGT	GCA	GGT	CGT	GAT	TTG	CAG	CGT	2064
Pro	Gly	Trp	Gln	Gly	Thr	Gly	Trp	Arg	Ala	Gly	Arg	Asp	Leu	Gln	Arg	
		675					680					685				
CGC	GAT	ATC	GGC	GGG	AAA	ACC	GGG	ACC	ACT	AAC	AGT	TCG	AAA	GAT	GCG	2112
Arg	Asp	Ile	Gly	Gly	Lys	Thr	Gly	Thr	Thr	Asn	Ser	Ser	Lys	Asp	Ala	
	690					695					700					
TGG	TTC	TCG	GGT	TAC	GGT	CCG	GGC	GTT	GTG	ACC	TCG	GTC	TGG	ATT	GGC	2160
Trp	Phe	Ser	Gly	Tyr	Gly	Pro	Gly	Val	Val	Thr	Ser	Val	Trp	Ile	Gly	
705					710					715					720	
TTT	GAT	GAT	CAC	CGT	CGT	AAT	CTC	GGT	CAT	ACA	ACG	GCT	TCC	GGA	GCG	2208
Phe	Asp	Asp	His	Arg	Arg	Asn	Leu	Gly	His	Thr	Thr	Ala	Ser	Gly	Ala	
				725					730					735		
ATT	AAA	GAT	CAG	ATC	TCA	GGT	TAC	GAA	GGC	GGT	GCC	AAG	AGT	GCC	CAG	2256
Ile	Lys	Asp	Gln	Ile	Ser	Gly	Tyr	Glu	Gly	Gly	Ala	Lys	Ser	Ala	Gln	
			740					745					750			
CCT	GCA	TGG	GAC	GCT	TAT	ATG	AAA	GCC	GTT	CTT	GAA	GGT	GTG	CCG	GAG	2304
Pro	Ala	Trp	Asp	Ala	Tyr	Met	Lys	Ala	Val	Leu	Glu	Gly	Val	Pro	Glu	
		755					760					765				
CAG	CCG	CTG	ACG	CCG	CCA	CCG	GGT	ATT	GTG	ACG	GTG	AAT	ATC	GAT	CGC	2352
Gln	Pro	Leu	Thr	Pro	Pro	Pro	Gly	Ile	Val	Thr	Val	Asn	Ile	Asp	Arg	
	770					775					780					
AGC	ACC	GGG	CAG	TTA	GCT	AAT	GGT	GGC	AAC	AGC	CGC	GAA	GAG	TAT	TTC	2400
Ser	Thr	Gly	Gln	Leu	Ala	Asn	Gly	Gly	Asn	Ser	Arg	Glu	Glu	Tyr	Phe	
	785				790					795					800	
ATC	GAA	GGT	ACG	CAG	CCG	ACA	CAA	CAG	GCA	GTG	CAC	GAG	GTG	GGA	ACG	2448
Ile	Glu	Gly	Thr		Gln	Pro	Thr	Gln	Gln	Ala	Val	His	Glu	Val	Gly	
				805					810					815		
ACC	ATT	ATC	GAT	AAT	GGC	GAG	GCA	CAG	GAA	TTG	TTG	TG				2487
Thr	Ile	Ile	Asp	Asn	Gly	Glu	Ala	Gln	Glu	Leu	Leu					
			820					825								

## (2) INFORMATION FOR SEQ ID NO: 2:

- (i) SEQUENCE CHARACTERISTICS:
  - (A) LENGTH: 828 amino acids
  - (B) TYPE: amino acid
  - (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: protein

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 2:

Met Gly Leu Tyr Arg Tyr Ile Glu Pro Gln Leu Pro Asp Val Ala Thr  
 1 5 10 15  
 Leu Lys Asp Val Arg Leu Gln Ile Pro Met Gln Ile Tyr Ser Ala Asp  
 20 25 30  
 Gly Glu Leu Ile Ala Gln Tyr Gly Glu Lys Arg Arg Ile Pro Val Thr  
 35 40 45  
 Leu Asp Gln Ile Pro Pro Glu Met Val Lys Ala Phe Ile Ala Thr Glu  
 50 55 60  
 Asp Ser Arg Phe Tyr Glu His His Gly Val Asp Pro Val Gly Ile Phe  
 65 70 75 80  
 Arg Ala Ala Ser Val Ala Leu Phe Ser Gly His Ala Ser Gln Gly Ala  
 85 90 95  
 Ser Thr Ile Thr Gln Gln Leu Ala Arg Asn Phe Phe Leu Ser Pro Glu  
 100 105 110  
 Arg Thr Leu Met Arg Lys Ile Lys Glu Val Phe Leu Ala Ile Arg Ile  
 115 120 125  
 Glu Gln Leu Leu Thr Lys Asp Glu Ile Leu Glu Leu Tyr Leu Asn Lys  
 130 135 140  
 Ile Tyr Leu Gly Tyr Arg Ala Tyr Gly Val Gly Ala Ala Ala Gln Val  
 145 150 155 160  
 Tyr Phe Gly Lys Thr Val Asp Gln Leu Thr Leu Asn Glu Met Ala Val  
 165 170 175  
 Ile Ala Gly Leu Pro Lys Ala Pro Ser Thr Phe Asn Pro Leu Tyr Ser  
 180 185 190  
 Met Asp Arg Ala Val Ala Arg Arg Asn Val Val Leu Ser Arg Met Leu  
 195 200 205  
 Asp Glu Gly Tyr Ile Thr Gln Gln Gln Phe Asp Gln Thr Arg Thr Glu  
 210 215 220  
 Ala Ile Asn Ala Asn Tyr His Ala Pro Glu Ile Ala Phe Ser Ala Pro  
 225 230 235 240  
 Tyr Leu Ser Glu Met Val Arg Gln Glu Met Tyr Asn Arg Tyr Gly Glu  
 245 250 255  
 Ser Ala Tyr Glu Asp Gly Tyr Arg Ile Tyr Thr Thr Ile Thr Arg Lys  
 260 265 270  
 Val Gln Gln Ala Ala Gln Gln Ala Val Arg Asn Asn Val Leu Asp Tyr  
 275 280 285  
 Asp Met Arg His Gly Tyr Arg Gly Pro Ala Asn Val Leu Trp Lys Val  
 290 295 300  
 Gly Glu Ser Ala Trp Asp Asn Asn Lys Ile Thr Asp Thr Leu Lys Ala  
 305 310 315 320  
 Leu Pro Thr Tyr Gly Pro Leu Leu Pro Ala Ala Val Thr Ser Ala Asn  
 325 330 335

Pro Gln Gln Ala Thr Ala Met Leu Ala Asp Gly S r Thr Val Ala Leu  
 340 345 350  
 Ser Met Glu Gly Val Arg Trp Ala Arg Pro Tyr Arg Ser Asp Thr Gln  
 355 360 365  
 Gln Gly Pro Thr Pro Arg Lys Val Thr Asp Val Leu Gln Thr Gly Gln  
 370 375 380  
 Gln Ile Trp Val Arg Gln Val Gly Asp Ala Trp Trp Leu Ala Gln Val  
 385 390 395 400  
 Pro Glu Val Asn Ser Ala Leu Val Ser Ile Asn Pro Gln Asn Gly Ala  
 405 410 415  
 Val Met Ala Leu Val Gly Gly Phe Asp Phe Asn Gln Ser Lys Phe Asn  
 420 425 430  
 Arg Ala Thr Gln Ala Leu Arg Gln Val Gly Ser Asn Ile Lys Pro Phe  
 435 440 445  
 Leu Tyr Thr Ala Ala Met Asp Lys Gly Leu Thr Leu Ala Ser Met Leu  
 450 455 460  
 Asn Asp Val Pro Ile Ser Arg Trp Asp Ala Ser Ala Gly Ser Asp Trp  
 465 470 475 480  
 Gln Pro Lys Asn Ser Pro Pro Gln Tyr Ala Gly Pro Ile Arg Leu Arg  
 485 490 495  
 Gln Gly Leu Gly Gln Ser Lys Asn Val Val Met Val Arg Ala Met Arg  
 500 505 510  
 Ala Met Gly Val Asp Tyr Ala Ala Glu Tyr Leu Gln Arg Phe Gly Phe  
 515 520 525  
 Pro Ala Gln Asn Ile Val His Thr Glu Ser Leu Ala Leu Gly Ser Ala  
 530 535 540  
 Ser Phe Thr Pro Met Gln Val Ala Arg Gly Tyr Ala Val Met Ala Asn  
 545 550 555 560  
 Gly Gly Phe Leu Val Asp Pro Trp Phe Ile Ser Lys Ile Glu Asn Asp  
 565 570 575  
 Gln Gly Gly Val Ile Phe Glu Ala Lys Pro Lys Val Ala Cys Pro Glu  
 580 585 590  
 Cys Asp Ile Pro Val Ile Tyr Gly Asp Thr Gln Lys Ser Asn Val Leu  
 595 600 605  
 Glu Asn Asn Asp Val Glu Asp Val Ala Ile Ser Arg Glu Gln Gln Asn  
 610 615 620  
 Val Ser Val Pro Met Pro Gln Leu Glu Gln Ala Asn Gln Ala Leu Val  
 625 630 635 640  
 Ala Lys Thr Gly Ala Gln Glu Tyr Ala Pro His Val Ile Asn Thr Pro  
 645 650 655  
 Leu Ala Phe Leu Ile Lys Ser Ala Leu Asn Thr Asn Ile Phe Gly Glu  
 660 665 670  
 Pro Gly Trp Gln Gly Thr Gly Trp Arg Ala Gly Arg Asp Leu Gln Arg  
 675 680 685  
 Arg Asp Ile Gly Gly Lys Thr Gly Thr Thr Asn Ser Ser Lys Asp Ala  
 690 695 700

Trp Phe Ser Gly Tyr Gly Pro Gly Val Val Thr Ser Val Trp Ile Gly  
 705 710 715 720  
 Phe Asp Asp His Arg Arg Asn Leu Gly His Thr Thr Ala Ser Gly Ala  
 725 730 735  
 Ile Lys Asp Gln Ile Ser Gly Tyr Glu Gly Gly Ala Lys Ser Ala Gln  
 740 745 750  
 Pro Ala Trp Asp Ala Tyr Met Lys Ala Val Leu Glu Gly Val Pro Glu  
 755 760 765  
 Gln Pro Leu Thr Pro Pro Pro Gly Ile Val Thr Val Asn Ile Asp Arg  
 770 775 780  
 Ser Thr Gly Gln Leu Ala Asn Gly Gly Asn Ser Arg Glu Glu Tyr Phe  
 785 790 795 800  
 Ile Glu Gly Thr Gln Pro Thr Gln Gln Ala Val His Glu Val Gly Thr  
 805 810 815  
 Thr Ile Ile Asp Asn Gly Glu Ala Gln Glu Leu Leu  
 820 825

## (2) INFORMATION FOR SEQ ID NO: 3:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 2472 base pairs
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: double
- (D) TOPOLOGY: linear

## (ii) MOLECULE TYPE: DNA (genomic)

## (vi) ORIGINAL SOURCE:

- (A) ORGANISM: Escherichia coli
- (B) STRAIN: DH5 alpha

## (vii) IMMEDIATE SOURCE:

- (A) LIBRARY: PCR cloning
- (B) CLONE: pARC 0556 Soluble PBP 1B

## (ix) FEATURE:

- (A) NAME/KEY: CDS
- (B) LOCATION: 1..2472

## (ix) FEATURE:

- (A) NAME/KEY: mat\_peptide
- (B) LOCATION: 1..2469

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 3:

ATG	GCC	GGG	AAT	GAC	CGC	GAG	CCA	ATT	GGA	CGC	AAA	GGG	AAA	CCG	ACG	48
Met	Ala	Gly	Asn	Asp	Arg	Glu	Pro	Ile	Gly	Arg	Lys	Gly	Lys	Pro	Thr	
1			5						10					15		
CGT	CCG	GTC	AAA	CAA	AAG	GTA	AGC	CGT	CGT	CGT	TAC	GAA	GAT	GAC	GAT	96
Arg	Pro	Val	Lys	Gln	Lys	Val	Ser	Arg	Arg	Arg	Tyr	Glu	Asp	Asp	Asp	
		20						25					30			
GAT	TAC	GAC	GAT	TAT	GAT	GAC	TAT	GAG	GAT	GAA	GAA	CCG	ATG	CCG	CGC	144
Asp	Tyr	Asp	Asp	Tyr	Asp	Asp	Tyr	Glu	Asp	Glu	Glu	Pro	Met	Pro	Arg	
		35					40					45				

AAA Lys	GGT Gly	AAG Lys	GGC Gly	AAA Lys	GGC Gly	AAA Lys	GGG Gly	CGT Arg	AAG Lys	CCT Pro	CGT Arg	GGC Gly	AAA Lys	CGC Arg	GGC Gly	192
50						55				60						
TCG Ser	ATC Ile	GAT Asp	CAA Gln	AAA Lys	ATT Ile	CGT Arg	AGC Ser	CGT Arg	ATT Ile	GAT Asp	GGC Gly	AAG Lys	GTC Val	TGG Trp	CAA Gln	240
65					70					75					80	
CTC Leu	GCT Ala	GCG Ala	GCA Ala	GTT Val	TAT Tyr	GGC Gly	CGA Arg	ATG Met	GTC Val	AAT Asn	CTT Leu	GAG Glu	CCA Pro	GAC Asp	ATG Met	288
				85					90					95		
ACC Thr	ATC Ile	AGC Ser	AAG Lys	AAC Asn	GAG Glu	ATG Met	GTG Val	AAG Lys	CTG Leu	CTG Leu	GAG Glu	GCG Ala	ACC Thr	CAG Gln	TAT Tyr	336
			100					105					110			
CGT Arg	CAG Gln	GTG Val	TCG Ser	AAA Lys	ATG Met	ACC Thr	CGT Arg	CCT Pro	GGC Gly	GAA Glu	TTT Phe	ACC Thr	GTG Val	CAG Gln	GCC Ala	384
		115					120					125				
AAC Asn	AGC Ser	ATT Ile	GAG Glu	ATG Met	ATT Ile	CGC Arg	CGT Arg	CCG Pro	TTT Phe	GAT Asp	TTC Phe	CCG Pro	GAC Asp	AGT Ser	AAA Lys	432
130						135					140					
GAA Glu	GGA Gly	CAG Gln	GTG Val	CGC Arg	GCG Ala	CGT Arg	CTG Leu	ACC Thr	TTT Phe	GAT Asp	GGC Gly	GAT Asp	CAT His	CTG Leu	GCG Ala	480
145					150					155				160		
ACG Thr	ATC Ile	GTC Val	AAT Asn	ATG Met	GAG Glu	AAC Asn	AAC Asn	CGT Arg	CAG Gln	TTC Phe	GGT Gly	TTC Phe	TTC Phe	CGT Arg	CTT Leu	528
				165					170					175		
GAT Asp	CCG Pro	CGT Arg	CTG Leu	ATC Ile	ACC Thr	ATG Met	ATC Ile	TCT Ser	TCG Ser	CCA Pro	AAC Asn	GGT Gly	GAG Glu	CAG Gln	CGT Arg	576
			180					185					190			
CTG Leu	TTT Phe	GTG Val	CCG Pro	CGC Arg	AGT Ser	GGT Gly	TTC Phe	CCG Pro	GAT Asp	TTG Leu	CTG Leu	GTG Val	GAT Asp	ACT Thr	TTG Leu	624
		195					200					205				
CTG Leu	GCG Ala	ACA Thr	GAA Glu	GAC Asp	CGT Arg	CAT His	TTT Phe	TAC Tyr	GAG Glu	CAT His	GAT Asp	GGA Gly	ATC Ile	AGT Ser	CTC Leu	672
	210					215					220					
TAC Tyr	TCA Ser	ATC Ile	GGA Gly	CGT Arg	GCG Ala	GTG Val	CTG Leu	GCA Ala	AAC Asn	CTG Leu	ACC Thr	GCC Ala	GGA Gly	CGC Arg	ACG Thr	720
225					230					235					240	
GTA Val	CAG Gln	GGT Gly	GCG Ala	AGT Ser	ACG Thr	CTG Leu	ACG Thr	CAA Gln	CAG Gln	CTG Leu	GTG Val	AAA Lys	AAC Asn	CTG Leu	TTC Phe	768
				245					250					255		
CTC Leu	TCC Ser	AGC Ser	GAG Glu	CGT Arg	TCT Ser	TAC Tyr	TGG Trp	CGT Arg	AAA Lys	GCG Ala	AAC Asn	GAA Glu	GCT Ala	TAC Tyr	ATG Met	816
			260					265					270			
GCG Ala	CTG Leu	ATC Ile	ATG Met	GAC Asp	GCG Ala	CGT Arg	TAC Tyr	AGC Ser	AAA Lys	GAC Asp	CGT Arg	ATT Ile	CTT Leu	GAG Glu	CTG Leu	864
		275					280					285				
TAT Tyr	ATG Met	AAC Asn	GAG Glu	GTG Val	TAT Tyr	CTC Leu	GGT Gly	CAG Gln	AGC Ser	GGC Gly	GAC Asp	AAC Asn	GAA Glu	ATC Ile	CGC Arg	912
	290					295					300					
GGC Gly	TTC Phe	CCG Pro	CTG Leu	GCA Ala	AGC Ser	TTG Leu	TAT Tyr	TAC Tyr	TTT Phe	GGT Gly	CGC Arg	CCG Pro	GTA Val	GAA Glu	GAG Glu	960
305					310					315					320	



CTA AGC CTC GAC CAG CAG GCG CTG TTA GTC GGT ATG GTG AAA GGG GCG L u Ser Leu Asp Gln Gln Ala Leu Leu Val Gly Met Val Lys Gly Ala 325 330 335	1008
TCC ATC TAC AAC CCG TGG CGT AAC CCA AAA CTG GCG CTG GAG CGA CGT S r Ile Tyr Asn Pro Trp Arg Asn Pro Lys Leu Ala Leu Glu Arg Arg 340 345 350	1056
AAT CTG GTG CTG CGT CTG CTG CAA CAG CAA CAG ATT ATT GAT CAA GAA Asn Leu Val Leu Arg Leu Leu Gln Gln Gln Ile Ile Asp Gln Glu 355 360 365	1104
CTC TAT GAC ATG TTG AGT GCC CGT CCG CTG GGG GTT CAG CCG CGC GGT Leu Tyr Asp Met Leu Ser Ala Arg Pro Leu Gly Val Gln Pro Arg Gly 370 375 380	1152
GGG GTG ATC TCT CCT CAG CCA GCC TTT ATG CAA CTG GTG CGT CAG GAG Gly Val Ile Ser Pro Gln Pro Ala Phe Met Gln Leu Val Arg Gln Glu 385 390 395 400	1200
CTG CAG GCA AAA CTG GGC GAT AAG GTA AAA GAT CTC TCC GGC GTG AAG Leu Gln Ala Lys Leu Gly Asp Lys Val Lys Asp Leu Ser Gly Val Lys 405 410 415	1248
ATC TTC ACT ACC TTT GAC TCG GTG GCC CAG GAC GCG GCA GAA AAA GCC Ile Phe Thr Thr Phe Asp Ser Val Ala Gln Asp Ala Ala Glu Lys Ala 420 425 430	1296
GCC GTG GAA GGC ATT CCG GCA CTG AAG AAA CAG CGT AAG TTG AGC GAT Ala Val Glu Gly Ile Pro Ala Leu Lys Lys Gln Arg Lys Leu Ser Asp 435 440 445	1344
CTT GAA ACT GCG ATT GTG GTC GTC GAC CGC TTT AGT GGT GAA GTT CGT Leu Glu Thr Ala Ile Val Val Val Asp Arg Phe Ser Gly Glu Val Arg 450 455 460	1392
GCG ATG GTC GGA GGT TCT GAG CCG CAG TTT GCG GGC TAC AAC CGT GCG Ala Met Val Gly Gly Ser Glu Pro Gln Phe Ala Gly Tyr Asn Arg Ala 465 470 475 480	1440
ATG CAG GCG CGT CGT TCG ATT GGT TCC CTT GCA AAA CCA GCG ACT TAT Met Gln Ala Arg Arg Ser Ile Gly Ser Leu Ala Lys Pro Ala Thr Tyr 485 490 495	1488
CTG ACG GCC TTA AGC CAG CCG AAA ATC TAT CGT CTG AAT ACG TGG ATT Leu Thr Ala Leu Ser Gln Pro Lys Ile Tyr Arg Leu Asn Thr Trp Ile 500 505 510	1536
GCG GAT GCG CCA ATT GCG CTG CGT CAG CCG AAT GGC CAG GTC TGG TCA Ala Asp Ala Pro Ile Ala Leu Arg Gln Pro Asn Gly Gln Val Trp Ser 515 520 525	1584
CCG CAG AAT GAT GAC CGT CGT TAT AGC GAA AGC GGC AGA GTG ATG CTG Pro Gln Asn Asp Asp Arg Arg Tyr Ser Glu Ser Gly Arg Val Met Leu 530 535 540	1632
GTG GAT GCG TTG ACC CGT TCG ATG AAC GTG CCG ACG GTA AAT CTG GGG Val Asp Ala Leu Thr Arg Ser Met Asn Val Pro Thr Val Asn Leu Gly 545 550 555 560	1680
ATG GCG CTG GGG CTG CCT GCG GTT ACG GAG ACC TGG ATT AAA CTG GGC Met Ala Leu Gly Leu Pro Ala Val Thr Glu Thr Trp Ile Lys Leu Gly 565 570 575	1728
GTA CCG AAA GAT CAG TTG CAT CCG GTT CCG GCA ATG CTG CTG GGG GCG Val Pro Lys Asp Gln Leu His Pro Val Pro Ala Met Leu Leu Gly Ala 580 585 590	1776

TTG AAC TTA ACG CCA ATC GAA GTG GCG CAG GCA TTC CAG ACC ATC GCC	1824
Leu Asn Leu Thr Pro Ile Glu Val Ala Gln Ala Phe Gln Thr Ile Ala	
595 600 605	
AGC GGT GGT AAC CGT GCA CCG CTT TCT GCG CTG CGT TCG GTA ATC GCG	1872
Ser Gly Gly Asn Arg Ala Pro Leu Ser Ala Leu Arg Ser Val Ile Ala	
610 615 620	
GAA GAT GGC AAA GTG CTG TAT CAG AGC TTC CCG CAG GCG GAA CGC GCT	1920
Glu Asp Gly Lys Val Leu Tyr Gln Ser Phe Pro Gln Ala Glu Arg Ala	
625 630 635 640	
GTT CCG GCG CAG GCG GCG TAT CTG ACA CTA TGG ACC ATG CAG CAG GTG	1968
Val Pro Ala Gln Ala Ala Tyr Leu Thr Leu Trp Thr Met Gln Gln Val	
645 650 655	
GTA CAA CGC GGT ACG GGT CGT CAG CTT GGG GCG AAA TAC CCG AAC CTG	2016
Val Gln Arg Gly Thr Gly Arg Gln Leu Gly Ala Lys Tyr Pro Asn Leu	
660 665 670	
CAT CTG GCA GGG AAA ACA GGG ACT ACC AAC AAT AAC GTA GAT ACC TGG	2064
His Leu Ala Gly Lys Thr Gly Thr Thr Asn Asn Asn Val Asp Thr Trp	
675 680 685	
TTT GCG GGC ATT GAC GGC AGC ACG GTG ACC ATC ACC TGG GTC GGC CGT	2112
Phe Ala Gly Ile Asp Gly Ser Thr Val Thr Ile Thr Trp Val Gly Arg	
690 695 700	
GAT AAC AAC CAG CCG ACC AAA CTG TAT GGT GCC AGC GGG GCA ATG TCG	2160
Asp Asn Asn Gln Pro Thr Lys Leu Tyr Gly Ala Ser Gly Ala Met Ser	
705 710 715 720	
ATT TAT CAG CGT TAT CTG GCT AAC CAG ACG CCA ACG CCG CTG AAT CTT	2208
Ile Tyr Gln Arg Tyr Leu Ala Asn Gln Thr Pro Thr Pro Leu Asn Leu	
725 730 735	
GTT CCG CCA GAA GAT ATT GCA GAT ATG GGC GTG GAC TAC GAC GGC AAC	2256
Val Pro Pro Glu Asp Ile Ala Asp Met Gly Val Asp Tyr Asp Gly Asn	
740 745 750	
TTT GTT TGC AGC GGT GGC ATG CGT ATC TTG CCG GTC TGG ACC AGC GAT	2304
Phe Val Cys Ser Gly Gly Met Arg Ile Leu Pro Val Trp Thr Ser Asp	
755 760 765	
CCG CAA TCG CTG TGC CAG CAG AGC GAG ATG CAG CAG CAG CCG TCA GGC	2352
Pro Gln Ser Leu Cys Gln Gln Ser Glu Met Gln Gln Gln Pro Ser Gly	
770 775 780	
AAT CCG TTT GAT CAG TCT TCT CAG CCG CAG CAA CAG CCG CAA CAG CAA	2400
Asn Pro Phe Asp Gln Ser Ser Gln Pro Gln Gln Gln Pro Gln Gln Gln	
785 790 795 800	
CCT GCT CAG CAA GAG CAG AAA GAC AGC GAC GGT GTA GCC GGT TGG ATC	2448
Pro Ala Gln Gln Glu Gln Lys Asp Ser Asp Gly Val Ala Gly Trp Ile	
805 810 815	
AAG GAT ATG TTT GGT AGT AAT TA	2472
Lys Asp Met Phe Gly Ser Asn	
820	

## (2) INFORMATION FOR SEQ ID NO: 4:

- (i) SEQUENCE CHARACTERISTICS:  
 (A) LENGTH: 823 amino acids  
 (B) TYPE: amino acid  
 (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: prot in

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 4:

M t Ala Gly Asn Asp Arg Glu Pro Ile Gly Arg Lys Gly Lys Pro Thr  
1 5 10 15  
Arg Pro Val Lys Gln Lys Val Ser Arg Arg Tyr Glu Asp Asp Asp  
20 25 30  
Asp Tyr Asp Asp Tyr Asp Asp Tyr Glu Asp Glu Glu Pro Met Pro Arg  
35 40 45  
Lys Gly Lys Gly Lys Gly Lys Gly Arg Lys Pro Arg Gly Lys Arg Gly  
50 55 60  
Ser Ile Asp Gln Lys Ile Arg Ser Arg Ile Asp Gly Lys Val Trp Gln  
65 70 75 80  
Leu Ala Ala Ala Val Tyr Gly Arg Met Val Asn Leu Glu Pro Asp Met  
85 90 95  
Thr Ile Ser Lys Asn Glu Met Val Lys Leu Leu Glu Ala Thr Gln Tyr  
100 105 110  
Arg Gln Val Ser Lys Met Thr Arg Pro Gly Glu Phe Thr Val Gln Ala  
115 120 125  
Asn Ser Ile Glu Met Ile Arg Arg Pro Phe Asp Phe Pro Asp Ser Lys  
130 135 140  
Glu Gly Gln Val Arg Ala Arg Leu Thr Phe Asp Gly Asp His Leu Ala  
145 150 155 160  
Thr Ile Val Asn Met Glu Asn Asn Arg Gln Phe Gly Phe Phe Arg Leu  
165 170 175  
Asp Pro Arg Leu Ile Thr Met Ile Ser Ser Pro Asn Gly Glu Gln Arg  
180 185 190  
Leu Phe Val Pro Arg Ser Gly Phe Pro Asp Leu Leu Val Asp Thr Leu  
195 200 205  
Leu Ala Thr Glu Asp Arg His Phe Tyr Glu His Asp Gly Ile Ser Leu  
210 215 220  
Tyr Ser Ile Gly Arg Ala Val Leu Ala Asn Leu Thr Ala Gly Arg Thr  
225 230 235 240  
Val Gln Gly Ala Ser Thr Leu Thr Gln Gln Leu Val Lys Asn Leu Phe  
245 250 255  
Leu Ser Ser Glu Arg Ser Tyr Trp Arg Lys Ala Asn Glu Ala Tyr Met  
260 265 270  
Ala Leu Ile Met Asp Ala Arg Tyr Ser Lys Asp Arg Ile Leu Glu Leu  
275 280 285  
Tyr Met Asn Glu Val Tyr Leu Gly Gln Ser Gly Asp Asn Glu Ile Arg  
290 295 300  
Gly Phe Pro Leu Ala Ser Leu Tyr Tyr Phe Gly Arg Pro Val Glu Glu  
305 310 315 320  
Leu Ser Leu Asp Gln Gln Ala Leu Leu Val Gly Met Val Lys Gly Ala  
325 330 335

Ser Ile Tyr Asn Pro Trp Arg Asn Pro Lys Leu Ala Leu Glu Arg Arg  
 340 345 350  
 Asn Leu Val Leu Arg L u Leu Gln Gln Gln Gln Ile Ile Asp Gln Glu  
 355 360 365  
 Leu Tyr Asp Met Leu Ser Ala Arg Pro Leu Gly Val Gln Pro Arg Gly  
 370 375 380  
 Gly Val Ile Ser Pro Gln Pro Ala Phe Met Gln Leu Val Arg Gln Glu  
 385 390 395 400  
 Leu Gln Ala Lys Leu Gly Asp Lys Val Lys Asp Leu Ser Gly Val Lys  
 405 410 415  
 Ile Phe Thr Thr Phe Asp Ser Val Ala Gln Asp Ala Ala Glu Lys Ala  
 420 425 430  
 Ala Val Glu Gly Ile Pro Ala Leu Lys Lys Gln Arg Lys Leu Ser Asp  
 435 440 445  
 Leu Glu Thr Ala Ile Val Val Val Asp Arg Phe Ser Gly Glu Val Arg  
 450 455 460  
 Ala Met Val Gly Gly Ser Glu Pro Gln Phe Ala Gly Tyr Asn Arg Ala  
 465 470 475 480  
 Met Gln Ala Arg Arg Ser Ile Gly Ser Leu Ala Lys Pro Ala Thr Tyr  
 485 490 495  
 Leu Thr Ala Leu Ser Gln Pro Lys Ile Tyr Arg Leu Asn Thr Trp Ile  
 500 505 510  
 Ala Asp Ala Pro Ile Ala Leu Arg Gln Pro Asn Gly Gln Val Trp Ser  
 515 520 525  
 Pro Gln Asn Asp Asp Arg Arg Tyr Ser Glu Ser Gly Arg Val Met Leu  
 530 535 540  
 Val Asp Ala Leu Thr Arg Ser Met Asn Val Pro Thr Val Asn Leu Gly  
 545 550 555 560  
 Met Ala Leu Gly Leu Pro Ala Val Thr Glu Thr Trp Ile Lys Leu Gly  
 565 570 575  
 Val Pro Lys Asp Gln Leu His Pro Val Pro Ala Met Leu Leu Gly Ala  
 580 585 590  
 Leu Asn Leu Thr Pro Ile Glu Val Ala Gln Ala Phe Gln Thr Ile Ala  
 595 600 605  
 Ser Gly Gly Asn Arg Ala Pro Leu Ser Ala Leu Arg Ser Val Ile Ala  
 610 615 620  
 Glu Asp Gly Lys Val Leu Tyr Gln Ser Phe Pro Gln Ala Glu Arg Ala  
 625 630 635 640  
 Val Pro Ala Gln Ala Ala Tyr Leu Thr Leu Trp Thr Met Gln Gln Val  
 645 650 655  
 Val Gln Arg Gly Thr Gly Arg Gln Leu Gly Ala Lys Tyr Pro Asn Leu  
 660 665 670  
 His Leu Ala Gly Lys Thr Gly Thr Thr Asn Asn Asn Val Asp Thr Trp  
 675 680 685  
 Phe Ala Gly Ile Asp Gly Ser Thr Val Thr Ile Thr Trp Val Gly Arg  
 690 695 700

[illegible]

(2) INFORMATION FOR SEQ ID NO: 5:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 2049 base pairs  
(B) TYPE: nucleic acid  
(C) STRANDEDNESS: double  
(D) TOPOLOGY: linear

## (ii) MOLECULE TYPE: DNA (genomic)

(vi) ORIGINAL SOURCE:

- (A) ORGANISM: Streptococcus pneumoniae  
(B) STRAIN: PM 1

(vii) IMMEDIATE SOURCE:

- (A) LIBRARY: PCR cloning  
(B) CLONE: pARC 0512 Soluble PBP 1A del 38

(ix) **FEATURE:**

- (A) NAME, KEY: CDS  
(B) LOCATION: 1..2049

(ix) FEATURE:

- ```
(A) NAME/KEY: mat_peptide
(B) LOCATION: 1..2046
```

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 5:

|     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| ATG | GCT | CCT | AGC | CTA | TCC | GAG | AGT | AAA | CTA | GTT | GCA | ACA | ACT | TCT | AGT | 48  |
| Met | Ala | Pro | Ser | Leu | Ser | Glu | Ser | Lys | Leu | Val | Ala | Thr | Thr | Ser | Ser |     |
| 1   |     |     |     | 5   |     |     |     |     | 10  |     |     |     |     | 15  |     |     |
|     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |
| AAA | ATC | TAC | GAC | AAT | AAA | AAT | CAA | CTC | ATT | GCT | GAC | TTG | GGT | TCT | GAA | 96  |
| Lys | Ile | Tyr | Asp | Asn | Lys | Asn | Gln | Leu | Ile | Ala | Asp | Leu | Gly | Ser | Glu |     |
|     |     |     | 20  |     |     |     |     | 25  |     |     |     |     | 30  |     |     |     |
|     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |
| CGC | CGC | GTC | AAT | GCC | CAA | GCT | AAT | GAT | ATT | CCC | ACA | GAT | TTG | GTT | AAG | 144 |
| Arg | Arg | Val | Asn | Ala | Gln | Ala | Asn | Asp | Ile | Pro | Thr | Asp | Leu | Val | Lys |     |
|     |     | 35  |     |     |     |     | 40  |     |     |     |     | 45  |     |     |     |     |

|     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| GCA | ATC | GTT | TCT | ATC | GAA | GAC | CAT | CGC | TTC | TTC | GAC | CAC | AGG | GGG | ATT | 192 |
| Ala | Ile | Val | Ser | Ile | Glu | Asp | His | Arg | Phe | Ph  | Asp | His | Arg | Gly | Il  |     |
|     | 50  |     |     |     |     | 55  |     |     |     |     | 60  |     |     |     |     |     |
| GAT | ACC | ATC | CGT | ATC | CTG | GGA | GCT | TTC | TTG | CGC | AAT | CTG | CAA | AGC | AAT | 240 |
| Asp | Thr | Ile | Arg | Ile | Leu | Gly | Ala | Phe | Leu | Arg | Asn | Leu | Gln | Ser | Asn |     |
| 65  |     |     |     |     | 70  |     |     |     |     | 75  |     |     |     |     | 80  |     |
| TCC | CTC | CAA | GGT | GGA | TCA | GCT | CTC | ACT | CAA | CAG | TTG | ATT | AAG | TTG | ACT | 288 |
| Ser | Leu | Gln | Gly | Gly | Ser | Ala | Leu | Thr | Gln | Gln | Leu | Ile | Lys | Leu | Thr |     |
|     |     |     |     | 85  |     |     |     |     | 90  |     |     |     |     | 95  |     |     |
| TAC | TTT | TCA | ACT | TCG | ACT | TCC | GAC | CAG | ACT | ATT | TCT | CGT | AAG | GCT | CAG | 336 |
| Tyr | Phe | Ser | Thr | Ser | Thr | Ser | Asp | Gln | Thr | Ile | Ser | Arg | Lys | Ala | Gln |     |
|     |     |     | 100 |     |     |     |     | 105 |     |     |     |     | 110 |     |     |     |
| GAA | GCT | TGG | TTA | GCG | ATT | CAG | TTA | GAA | CAA | AAA | GCA | ACC | AAG | CAA | GAA | 384 |
| Glu | Ala | Trp | Leu | Ala | Ile | Gln | Leu | Glu | Gln | Lys | Ala | Thr | Lys | Gln | Glu |     |
|     |     | 115 |     |     |     |     | 120 |     |     |     |     | 125 |     |     |     |     |
| ATC | TTG | ACC | TAC | TAT | ATA | AAT | AAG | GTC | TAC | ATG | TCT | AAT | GGG | AAC | TAT | 432 |
| Ile | Leu | Thr | Tyr | Tyr | Ile | Asn | Lys | Val | Tyr | Met | Ser | Asn | Gly | Asn | Tyr |     |
|     | 130 |     |     |     |     | 135 |     |     |     |     | 140 |     |     |     |     |     |
| GGA | ATG | CAG | ACA | GCA | GCT | CAA | AAC | TAC | TAT | GGT | AAA | GAC | CTC | AAT | AAT | 480 |
| Gly | Met | Gln | Thr | Ala | Ala | Gln | Asn | Tyr | Tyr | Gly | Lys | Asp | Leu | Asn | Asn |     |
| 145 |     |     |     |     | 150 |     |     |     |     | 155 |     |     |     |     | 160 |     |
| TTA | AGT | TTA | CCT | CAG | TTA | GCC | TTG | CTG | GCT | GGG | ATG | CCT | CAG | GCA | CCA | 528 |
| Leu | Ser | Leu | Pro | Gln | Leu | Ala | Leu | Leu | Ala | Gly | Met | Pro | Gln | Ala | Pro |     |
|     |     |     |     | 165 |     |     |     |     | 170 |     |     |     |     | 175 |     |     |
| AAC | CAA | TAT | GAC | CCC | TAT | TCA | CAT | CCA | GAA | GCA | GCC | CAA | GAC | CGC | CGA | 576 |
| Asn | Gln | Tyr | Asp | Pro | Tyr | Ser | His | Pro | Glu | Ala | Ala | Gln | Asp | Arg | Arg |     |
|     |     |     | 180 |     |     |     |     | 185 |     |     |     |     | 190 |     |     |     |
| AAC | TTG | GTC | TTA | TCT | GAA | ATG | AAA | AAT | CAA | GGC | TAC | ATC | TCT | GCT | GAA | 624 |
| Asn | Leu | Val | Leu | Ser | Glu | Met | Lys | Asn | Gln | Gly | Tyr | Ile | Ser | Ala | Glu |     |
|     |     | 195 |     |     |     |     | 200 |     |     |     |     | 205 |     |     |     |     |
| CAG | TAT | GAG | AAA | GCA | GTC | AAT | ACA | CCA | ATT | ACT | GAT | GGG | CTA | CAA | AGT | 672 |
| Gln | Tyr | Glu | Lys | Ala | Val | Asn | Thr | Pro | Ile | Thr | Asp | Gly | Leu | Gln | Ser |     |
|     | 210 |     |     |     |     | 215 |     |     |     |     | 220 |     |     |     |     |     |
| CTC | AAA | TCA | GCA | AGT | AAT | TAC | CCT | GCT | TAC | ATG | GAT | AAT | TAC | CTC | AAG | 720 |
| Leu | Lys | Ser | Ala | Ser | Asn | Tyr | Pro | Ala | Tyr | Met | Asp | Asn | Tyr | Leu | Lys |     |
| 225 |     |     |     |     | 230 |     |     |     |     | 235 |     |     |     |     | 240 |     |
| GAA | GTC | ATC | AAT | CAA | GTT | GAA | GAA | GAA | ACA | GGC | TAT | AAC | CTA | CTC | ACA | 768 |
| Glu | Val | Ile | Asn | Gln | Val | Glu | Glu | Glu | Thr | Gly | Tyr | Asn | Leu | Leu | Thr |     |
|     |     |     |     | 245 |     |     |     |     | 250 |     |     |     |     | 255 |     |     |
| ACT | GGG | ATG | GAT | GTC | TAC | ACA | AAT | GTA | GAC | CAA | GAA | GCT | CAA | AAA | CAT | 816 |
| Thr | Gly | Met | Asp | Val | Tyr | Thr | Asn | Val | Asp | Gln | Glu | Ala | Gln | Lys | His |     |
|     |     |     | 260 |     |     |     |     | 265 |     |     |     |     | 270 |     |     |     |
| CTG | TGG | GAT | ATT | TAC | AAT | ACA | GAC | GAA | TAC | GTT | GCC | TAT | CCA | GAC | GAT | 864 |
| Leu | Trp | Asp | Ile | Tyr | Asn | Thr | Asp | Glu | Tyr | Val | Ala | Tyr | Pro | Asp | Asp |     |
|     |     | 275 |     |     |     |     | 280 |     |     |     |     | 285 |     |     |     |     |
| GAA | TTG | CAA | GTC | GCT | TCT | ACC | ATT | GTT | GAT | GTT | TCT | AAC | GGT | AAA | GTC | 912 |
| Glu | Leu | Gln | Val | Ala | Ser | Thr | Ile | Val | Asp | Val | Ser | Asn | Gly | Lys | Val |     |
|     | 290 |     |     |     |     | 295 |     |     |     |     | 300 |     |     |     |     |     |
| ATT | GCC | CAG | CTA | GGA | GCA | CGC | CAT | CAG | TCA | AGT | AAT | GTT | TCC | TTC | GGA | 960 |
| Ile | Ala | Gln | Leu | Gly | Ala | Arg | His | Gln | Ser | Ser | Asn | Val | Ser | Phe | Gly |     |
| 305 |     |     |     |     | 310 |     |     |     |     | 315 |     |     |     |     | 320 |     |

|                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |      |
|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|------|
| ATT<br>11         | AAC<br>Asn        | CAA<br>Gln        | GCA<br>Ala        | GTA<br>Val<br>325 | GAA<br>Glu        | ACA<br>Thr        | AAC<br>Asn        | CGC<br>Arg        | GAC<br>Asp<br>330 | TGG<br>Trp        | GGA<br>Gly        | TCA<br>Ser        | ACT<br>Thr        | ATG<br>Met<br>335 | AAA<br>Lys        | 1008 |
| CCG<br>Pro        | ATC<br>11         | ACA<br>Thr        | GAC<br>Asp<br>340 | TAT<br>Tyr        | GCT<br>Ala        | CCT<br>Pro        | GCC<br>Ala        | TTG<br>Leu<br>345 | GAG<br>Glu        | TAC<br>Tyr        | GGT<br>Gly        | GTC<br>Val        | TAC<br>Tyr        | GAG<br>Glu<br>350 | TCA<br>Ser        | 1056 |
| ACT<br>Thr        | GCC<br>Ala        | ACT<br>Thr        | ATC<br>Ile        | GTT<br>Val<br>355 | CAC<br>His        | GAT<br>Asp        | GAG<br>Glu<br>360 | CCC<br>Pro        | TAT<br>Tyr        | AAC<br>Asn        | TAC<br>Tyr        | CCT<br>Pro<br>365 | GGG<br>Gly        | ACA<br>Thr        | AAT<br>Asn        | 1104 |
| ACC<br>Thr        | CCT<br>Pro        | GTT<br>Val        | TAT<br>Tyr        | AAC<br>Asn        | TGG<br>Trp        | GAT<br>Asp<br>375 | AGG<br>Arg        | GGC<br>Gly        | TAC<br>Tyr        | TTT<br>Phe        | GGC<br>Gly<br>380 | AAC<br>Asn        | ATC<br>Ile        | ACC<br>Thr        | TTG<br>Leu        | 1152 |
| CAA<br>Gln<br>385 | TAC<br>Tyr        | GCC<br>Ala        | CTG<br>Leu        | CAA<br>Gln<br>390 | CAA<br>Gln        | TCG<br>Ser        | CGA<br>Arg        | AAC<br>Asn        | GTC<br>Val<br>395 | CCA<br>Pro        | GCC<br>Ala        | GTG<br>Val        | GAA<br>Glu        | ACT<br>Thr        | CTA<br>Leu<br>400 | 1200 |
| AAC<br>Asn        | AAG<br>Lys        | GTC<br>Val        | GGA<br>Gly        | CTC<br>Leu<br>405 | AAC<br>Asn        | CGC<br>Arg        | GCC<br>Ala        | AAG<br>Lys        | ACT<br>Thr<br>410 | TTC<br>Phe        | CTA<br>Leu        | AAT<br>Asn        | GGT<br>Gly        | CTC<br>Leu<br>415 | GGA<br>Gly        | 1248 |
| ATC<br>Ile        | GAC<br>Asp        | TAC<br>Tyr        | CCA<br>Pro<br>420 | AGT<br>Ser        | ATT<br>Ile        | CAC<br>His        | TAC<br>Tyr        | TCA<br>Ser<br>425 | AAT<br>Asn        | GCC<br>Ala        | ATT<br>Ile        | TCA<br>Ser        | AGT<br>Ser        | AAC<br>Asn        | ACA<br>Thr        | 1296 |
| ACC<br>Thr        | GAA<br>Glu        | TCA<br>Ser        | GAC<br>Asp<br>435 | AAA<br>Lys        | AAA<br>Lys        | TAT<br>Tyr        | GGA<br>Gly<br>440 | GCA<br>Ala        | AGT<br>Ser        | AGT<br>Ser        | GAA<br>Glu        | AAG<br>Lys        | ATG<br>Met        | GCT<br>Ala        | GCT<br>Ala        | 1344 |
| GCT<br>Ala        | TAC<br>Tyr<br>450 | GCT<br>Ala        | GCC<br>Ala        | TTT<br>Phe        | GCA<br>Ala        | AAT<br>Asn<br>455 | GGT<br>Gly        | GGA<br>Gly        | ACT<br>Thr        | TAC<br>Tyr        | TAT<br>Tyr<br>460 | AAA<br>Lys        | CCA<br>Pro        | ATG<br>Met        | TAT<br>Tyr        | 1392 |
| ATC<br>Ile<br>465 | CAT<br>His        | AAA<br>Lys        | GTC<br>Val        | GTC<br>Val        | TTT<br>Phe<br>470 | AGT<br>Ser        | GAT<br>Asp        | GGG<br>Gly        | AGT<br>Ser        | GAA<br>Glu<br>475 | AAA<br>Lys        | GAG<br>Glu        | TTC<br>Phe        | TCT<br>Ser        | AAT<br>Asn<br>480 | 1440 |
| GTC<br>Val        | GGA<br>Gly        | ACT<br>Thr        | CGT<br>Arg        | GCC<br>Ala<br>485 | ATG<br>Met        | AAA<br>Lys        | GAA<br>Glu        | ACG<br>Thr        | ACA<br>Thr<br>490 | GCC<br>Ala        | TAT<br>Tyr        | ATG<br>Met        | ATG<br>Met        | ACC<br>Thr        | GAC<br>Asp<br>495 | 1488 |
| ATG<br>Met        | ATG<br>Met        | AAA<br>Lys        | ACA<br>Thr<br>500 | GTC<br>Val        | TTG<br>Leu        | AGT<br>Ser        | TAT<br>Tyr        | GGA<br>Gly<br>505 | ACT<br>Thr        | GGA<br>Gly        | CGA<br>Arg        | AAT<br>Asn        | GCC<br>Ala<br>510 | TAT<br>Tyr        | CTT<br>Leu        | 1536 |
| GCT<br>Ala        | TGG<br>Trp        | CTC<br>Leu<br>515 | CCT<br>Pro        | CAG<br>Gln        | GCT<br>Ala        | GGT<br>Gly        | AAA<br>Lys<br>520 | ACA<br>Thr        | GGA<br>Gly        | ACC<br>Thr        | TCT<br>Ser        | AAC<br>Asn        | TAT<br>Tyr        | ACA<br>Thr        | GAC<br>Asp        | 1584 |
| GAG<br>Glu<br>530 | GAA<br>Glu        | ATT<br>Ile        | GAA<br>Glu        | AAC<br>Asn        | CAC<br>His        | ATC<br>Ile<br>535 | AAG<br>Lys        | ACC<br>Thr        | TCT<br>Ser        | CAA<br>Gln        | TTT<br>Phe        | GTA<br>Val        | GCA<br>Ala        | CCT<br>Pro        | GAT<br>Asp        | 1632 |
| GAA<br>Glu<br>545 | CTA<br>Leu        | TTT<br>Phe        | GCT<br>Ala        | GGC<br>Gly        | TAT<br>Tyr        | ACG<br>Thr<br>550 | CGT<br>Arg        | AAA<br>Lys        | TAT<br>Tyr        | TCA<br>Ser        | ATG<br>Met        | GCT<br>Ala        | GTA<br>Val        | TGG<br>Trp        | ACA<br>Thr<br>560 | 1680 |
| GGC<br>Gly        | TAT<br>Tyr        | TCT<br>Ser        | AAC<br>Asn        | CGT<br>Arg<br>565 | CTG<br>Leu        | ACA<br>Thr        | CCA<br>Pro        | CTT<br>Leu        | GTA<br>Val<br>570 | GGC<br>Gly        | AAT<br>Asn        | GGC<br>Gly        | CTT<br>Leu        | ACG<br>Thr        | GTC<br>Val        | 1728 |
| GCT<br>Ala        | GCC<br>Ala        | AAA<br>Lys        | GTT<br>Val<br>580 | TAC<br>Tyr        | CGC<br>Arg        | TCT<br>Ser        | ATG<br>Met        | ATG<br>Met<br>585 | ACC<br>Thr        | TAC<br>Tyr        | CTG<br>Leu        | TCT<br>Ser        | GAA<br>Glu        | GGA<br>Gly        | AGC<br>Ser        | 1776 |

|                                                                 |      |
|-----------------------------------------------------------------|------|
| AAT CCA GAG GAT TGG AAT ATA CCA GAG GGG CTC TAC AGA AAT GGA GAA | 1824 |
| Asn Pro Glu Asp Trp Asn Ile Pro Glu Gly Leu Tyr Arg Asn Gly Glu |      |
| 595 600 605                                                     |      |
| TTC GTA TTT AAA AAT GGT GCT CGT TCT ACG TGG AGC TCA CCT GCT CCA | 1872 |
| Phe Val Phe Lys Asn Gly Ala Arg Ser Thr Trp Ser Ser Pro Ala Pro |      |
| 610 615 620                                                     |      |
| CAA CAA CCC CCA TCA ACT GAA AGT TCA AGC TCA TCA TCA GAT AGT TCA | 1920 |
| Gln Gln Pro Pro Ser Thr Glu Ser Ser Ser Ser Ser Ser Asp Ser Ser |      |
| 625 630 635 640                                                 |      |
| ACT TCA CAG TCT AGC TCA ACC ACT CCA AGC ACA AAT AAT AGT ACG ACT | 1968 |
| Thr Ser Gln Ser Ser Ser Thr Thr Pro Ser Thr Asn Asn Ser Thr Thr |      |
| 645 650 655                                                     |      |
| ACC AAT CCT AAC AAT AAT ACG CAA CAA TCA AAT ACA ACC CCT GAT CAA | 2016 |
| Thr Asn Pro Asn Asn Asn Thr Gln Gln Ser Asn Thr Thr Pro Asp Gln |      |
| 660 665 670                                                     |      |
| CAA AAT CAG AAT CCT CAA CCA GCA CAA CCA TA                      | 2049 |
| Gln Asn Gln Asn Pro Gln Pro Ala Gln Pro                         |      |
| 675 680                                                         |      |

## (2) INFORMATION FOR SEQ ID NO: 6:

- (i) SEQUENCE CHARACTERISTICS:  
 (A) LENGTH: 682 amino acids  
 (B) TYPE: amino acid  
 (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: protein

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 6:

|     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| Met | Ala | Pro | Ser | Leu | Ser | Glu | Ser | Lys | Leu | Val | Ala | Thr | Thr | Ser | Ser |
| 1   |     |     |     | 5   |     |     |     |     | 10  |     |     |     |     | 15  |     |
| Lys | Ile | Tyr | Asp | Asn | Lys | Asn | Gln | Leu | Ile | Ala | Asp | Leu | Gly | Ser | Glu |
|     |     |     | 20  |     |     |     |     | 25  |     |     |     |     | 30  |     |     |
| Arg | Arg | Val | Asn | Ala | Gln | Ala | Asn | Asp | Ile | Pro | Thr | Asp | Leu | Val | Lys |
|     |     | 35  |     |     |     |     | 40  |     |     |     |     | 45  |     |     |     |
| Ala | Ile | Val | Ser | Ile | Glu | Asp | His | Arg | Phe | Phe | Asp | His | Arg | Gly | Ile |
|     | 50  |     |     |     |     | 55  |     |     |     |     | 60  |     |     |     |     |
| Asp | Thr | Ile | Arg | Ile | Leu | Gly | Ala | Phe | Leu | Arg | Asn | Leu | Gln | Ser | Asn |
| 65  |     |     |     |     | 70  |     |     |     |     | 75  |     |     |     |     | 80  |
| Ser | Leu | Gln | Gly | Gly | Ser | Ala | Leu | Thr | Gln | Gln | Leu | Ile | Lys | Leu | Thr |
|     |     |     | 85  |     |     |     |     |     | 90  |     |     |     | 95  |     |     |
| Tyr | Phe | Ser | Thr | Ser | Thr | Ser | Asp | Gln | Thr | Ile | Ser | Arg | Lys | Ala | Gln |
|     |     |     | 100 |     |     |     |     | 105 |     |     |     |     | 110 |     |     |
| Glu | Ala | Trp | Leu | Ala | Ile | Gln | Leu | Glu | Gln | Lys | Ala | Thr | Lys | Gln | Glu |
|     |     | 115 |     |     |     |     | 120 |     |     |     |     | 125 |     |     |     |
| Ile | Leu | Thr | Tyr | Tyr | Ile | Asn | Lys | Val | Tyr | Met | Ser | Asn | Gly | Asn | Tyr |
|     | 130 |     |     |     |     | 135 |     |     |     |     | 140 |     |     |     |     |
| Gly | Met | Gln | Thr | Ala | Ala | Gln | Asn | Tyr | Tyr | Gly | Lys | Asp | Leu | Asn | Asn |
| 145 |     |     |     |     | 150 |     |     |     |     | 155 |     |     |     |     | 160 |
| Leu | Ser | Leu | Pro | Gln | Leu | Ala | Leu | Leu | Ala | Gly | Met | Pro | Gln | Ala | Pro |
|     |     |     |     | 165 |     |     |     |     | 170 |     |     |     |     | 175 |     |



Asn Gln Tyr Asp Pro Tyr S r His Pro Glu Ala Ala Gln Asp Arg Arg  
 180 185 190  
 Asn Leu Val Leu Ser Glu Met Lys Asn Gln Gly Tyr Ile Ser Ala Glu  
 195 200 205  
 Gln Tyr Glu Lys Ala Val Asn Thr Pro Ile Thr Asp Gly Leu Gln Ser  
 210 215 220  
 Leu Lys Ser Ala Ser Asn Tyr Pro Ala Tyr Met Asp Asn Tyr Leu Lys  
 225 230 235 240  
 Glu Val Ile Asn Gln Val Glu Glu Glu Thr Gly Tyr Asn Leu Leu Thr  
 245 250 255  
 Thr Gly Met Asp Val Tyr Thr Asn Val Asp Gln Glu Ala Gln Lys His  
 260 265 270  
 Leu Trp Asp Ile Tyr Asn Thr Asp Glu Tyr Val Ala Tyr Pro Asp Asp  
 275 280 285  
 Glu Leu Gln Val Ala Ser Thr Ile Val Asp Val Ser Asn Gly Lys Val  
 290 295 300  
 Ile Ala Gln Leu Gly Ala Arg His Gln Ser Ser Asn Val Ser Phe Gly  
 305 310 315 320  
 Ile Asn Gln Ala Val Glu Thr Asn Arg Asp Trp Gly Ser Thr Met Lys  
 325 330 335  
 Pro Ile Thr Asp Tyr Ala Pro Ala Leu Glu Tyr Gly Val Tyr Glu Ser  
 340 345 350  
 Thr Ala Thr Ile Val His Asp Glu Pro Tyr Asn Tyr Pro Gly Thr Asn  
 355 360 365  
 Thr Pro Val Tyr Asn Trp Asp Arg Gly Tyr Phe Gly Asn Ile Thr Leu  
 370 375 380  
 Gln Tyr Ala Leu Gln Gln Ser Arg Asn Val Pro Ala Val Glu Thr Leu  
 385 390 395 400  
 Asn Lys Val Gly Leu Asn Arg Ala Lys Thr Phe Leu Asn Gly Leu Gly  
 405 410 415  
 Ile Asp Tyr Pro Ser Ile His Tyr Ser Asn Ala Ile Ser Ser Asn Thr  
 420 425 430  
 Thr Glu Ser Asp Lys Lys Tyr Gly Ala Ser Ser Glu Lys Met Ala Ala  
 435 440 445  
 Ala Tyr Ala Ala Phe Ala Asn Gly Gly Thr Tyr Tyr Lys Pro Met Tyr  
 450 455 460  
 Ile His Lys Val Val Phe Ser Asp Gly Ser Glu Lys Glu Phe Ser Asn  
 465 470 475 480  
 Val Gly Thr Arg Ala Met Lys Glu Thr Thr Ala Tyr Met Met Thr Asp  
 485 490 495  
 Met Met Lys Thr Val Leu Ser Tyr Gly Thr Gly Arg Asn Ala Tyr Leu  
 500 505 510  
 Ala Trp Leu Pro Gln Ala Gly Lys Thr Gly Thr Ser Asn Tyr Thr Asp  
 515 520 525  
 Glu Glu Ile Glu Asn His Ile Lys Thr Ser Gln Phe Val Ala Pro Asp  
 530 535 540

Glu Leu Ph Ala Gly Tyr Thr Arg Lys Tyr Ser M t Ala Val Trp Thr  
545 550 555 560

Gly Tyr Ser Asn Arg Leu Thr Pro Leu Val Gly Asn Gly Leu Thr Val  
565 570 575

Ala Ala Lys Val Tyr Arg Ser Met Met Thr Tyr Leu Ser Glu Gly Ser  
580 585 590

Asn Pro Glu Asp Trp Asn Ile Pro Glu Gly Leu Tyr Arg Asn Gly Glu  
595 600 605

Phe Val Phe Lys Asn Gly Ala Arg Ser Thr Trp Ser Ser Pro Ala Pro  
610 615 620

Gln Gln Pro Pro Ser Thr Glu Ser Ser Ser Ser Ser Ser Asp Ser Ser  
625 630 635 640

Thr Ser Gln Ser Ser Ser Thr Thr Pro Ser Thr Asn Asn Ser Thr Thr  
645 650 655

Thr Asn Pro Asn Asn Asn Thr Gln Gln Ser Asn Thr Thr Pro Asp Gln  
660 665 670

Gln Asn Gln Asn Pro Gln Pro Ala Gln Pro  
675 680

(2) INFORMATION FOR SEQ ID NO: 7:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 844 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(vi) ORIGINAL SOURCE:

- (A) ORGANISM: Escherichia coli

(vii) IMMEDIATE SOURCE:

- (B) CLONE: pARC0438 PEF 1B QQAA

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 7:

Met Ala Gly Asn Asp Arg Glu Pro Ile Gly Arg Lys Gly Lys Pro Thr  
1 5 10 15

Arg Pro Val Lys Gln Lys Val Ser Arg Arg Arg Tyr Glu Asp Asp Asp  
20 25 30

Asp Tyr Asp Asp Tyr Asp Asp Tyr Glu Asp Glu Glu Pro Met Pro Arg  
35 40 45

Lys Gly Lys Gly Lys Gly Lys Gly Arg Lys Pro Arg Gly Lys Arg Gly  
50 55 60

Trp Leu Trp Leu Leu Leu Lys Leu Ala Ile Val Phe Ala Val Leu Ile  
65 70 75 80

Ala Ile Tyr Gly Val Tyr Leu Asp Gln Lys Ile Arg Ser Arg Ile Asp  
85 90 95

Gly Lys Val Trp Gln Leu Ala Ala Ala Val Tyr Gly Arg Met Val Asn  
100 105 110

Leu Glu Pro Asp Met Thr Ile S r Lys Asn Glu Met Val Lys Leu Leu  
 115 120 125  
 Glu Ala Thr Gln Tyr Arg Gln Val Ser Lys Met Thr Arg Pro Gly Glu  
 130 135 140  
 Phe Thr Val Gln Ala Asn Ser Ile Glu Met Ile Arg Arg Pro Phe Asp  
 145 150 155 160  
 Phe Pro Asp Ser Lys Glu Gly Gln Val Arg Ala Arg Leu Thr Phe Asp  
 165 170 175  
 Gly Asp His Leu Ala Thr Ile Val Asn Met Glu Asn Asn Arg Gln Phe  
 180 185 190  
 Gly Phe Phe Arg Leu Asp Pro Arg Leu Ile Thr Met Ile Ser Ser Pro  
 195 200 205  
 Asn Gly Glu Gln Arg Leu Phe Val Pro Arg Ser Gly Phe Pro Asp Leu  
 210 215 220  
 Leu Val Asp Thr Leu Leu Ala Thr Glu Asp Arg His Phe Tyr Glu His  
 225 230 235 240  
 Asp Gly Ile Ser Leu Tyr Ser Ile Gly Arg Ala Val Leu Ala Asn Leu  
 245 250 255  
 Thr Ala Gly Arg Thr Val Gln Gly Ala Ser Thr Leu Thr Ala Ala Leu  
 260 265 270  
 Val Lys Asn Leu Phe Leu Ser Ser Glu Arg Ser Tyr Trp Arg Lys Ala  
 275 280 285  
 Asn Glu Ala Tyr Met Ala Leu Ile Met Asp Ala Arg Tyr Ser Lys Asp  
 290 295 300  
 Arg Ile Leu Glu Leu Tyr Met Asn Glu Val Tyr Leu Gly Gln Ser Gly  
 305 310 315 320  
 Asp Asn Glu Ile Arg Gly Phe Pro Leu Ala Ser Leu Tyr Tyr Phe Gly  
 325 330 335  
 Arg Pro Val Glu Glu Leu Ser Leu Asp Gln Gln Ala Leu Leu Val Gly  
 340 345 350  
 Met Val Lys Gly Ala Ser Ile Tyr Asn Pro Trp Arg Asn Pro Lys Leu  
 355 360 365  
 Ala Leu Glu Arg Arg Asn Leu Val Leu Arg Leu Leu Gln Gln Gln Gln  
 370 375 380  
 Ile Ile Asp Gln Glu Leu Tyr Asp Met Leu Ser Ala Arg Pro Leu Gly  
 385 390 395 400  
 Val Gln Pro Arg Gly Gly Val Ile Ser Pro Gln Pro Ala Phe Met Gln  
 405 410 415  
 Leu Val Arg Gln Glu Leu Gln Ala Lys Leu Gly Asp Lys Val Lys Asp  
 420 425 430  
 Leu Ser Gly Val Lys Ile Phe Thr Thr Phe Asp Ser Val Ala Gln Asp  
 435 440 445  
 Ala Ala Glu Lys Ala Ala Val Glu Gly Il Pro Ala Leu Lys Lys Gln  
 450 455 460  
 Arg Lys Leu Ser Asp Leu Glu Thr Ala Ile Val Val Val Asp Arg Phe  
 465 470 475 480

S r Gly Glu Val Arg Ala Met Val Gly Gly Ser Glu Pro Gln Phe Ala  
 485 490 495  
 Gly Tyr Asn Arg Ala M t Gln Ala Arg Arg Ser Ile Gly Ser Leu Ala  
 500 505 510  
 Lys Pro Ala Thr Tyr Leu Thr Ala Leu Ser Gln Pro Lys Ile Tyr Arg  
 515 520 525  
 Leu Asn Thr Trp Ile Ala Asp Ala Pro Ile Ala Leu Arg Gln Pro Asn  
 530 535 540  
 Gly Gln Val Trp Ser Pro Gln Asn Asp Asp Arg Arg Tyr Ser Glu Ser  
 545 550 555 560  
 Gly Arg Val Met Leu Val Asp Ala Leu Thr Arg Ser Met Asn Val Pro  
 565 570 575  
 Thr Val Asn Leu Gly Met Ala Leu Gly Leu Pro Ala Val Thr Glu Thr  
 580 585 590  
 Trp Ile Lys Leu Gly Val Pro Lys Asp Gln Leu His Pro Val Pro Ala  
 595 600 605  
 Met Leu Leu Gly Ala Leu Asn Leu Thr Pro Ile Glu Val Ala Gln Ala  
 610 615 620  
 Phe Gln Thr Ile Ala Ser Gly Gly Asn Arg Ala Pro Leu Ser Ala Leu  
 625 630 635 640  
 Arg Ser Val Ile Ala Glu Asp Gly Lys Val Leu Tyr Gln Ser Phe Pro  
 645 650 655  
 Gln Ala Glu Arg Ala Val Pro Ala Gln Ala Ala Tyr Leu Thr Leu Trp  
 660 665 670  
 Thr Met Gln Gln Val Val Gln Arg Gly Thr Gly Arg Gln Leu Gly Ala  
 675 680 685  
 Lys Tyr Pro Asn Leu His Leu Ala Gly Lys Thr Gly Thr Thr Asn Asn  
 690 695 700  
 Asn Val Asp Thr Trp Phe Ala Gly Ile Asp Gly Ser Thr Val Thr Ile  
 705 710 715 720  
 Thr Trp Val Gly Arg Asp Asn Asn Gln Pro Thr Lys Leu Tyr Gly Ala  
 725 730 735  
 Ser Gly Ala Met Ser Ile Tyr Gln Arg Tyr Leu Ala Asn Gln Thr Pro  
 740 745 750  
 Thr Pro Leu Asn Leu Val Pro Pro Glu Asp Ile Ala Asp Met Gly Val  
 755 760 765  
 Asp Tyr Asp Gly Asn Phe Val Cys Ser Gly Gly Met Arg Ile Leu Pro  
 770 775 780  
 Val Trp Thr Ser Asp Pro Gln Ser Leu Cys Gln Gln Ser Glu Met Gln  
 785 790 795 800  
 Gln Gln Pro Ser Gly Asn Pro Phe Asp Gln Ser Ser Gln Pro Gln Gln  
 805 810 815  
 Gln Pro Gln Gln Gln Pro Ala Gln Gln Glu Gln Lys Asp Ser Asp Gly  
 820 825 830  
 Val Ala Gly Trp Ile Lys Asp Met Phe Gly Ser Asn  
 835 840

## (2) INFORMATION FOR SEQ ID NO: 8:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 844 amino acids  
 (B) TYPE: amino acid  
 (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

## (vi) ORIGINAL SOURCE:

(A) ORGANISM: Escherichia coli

## (vii) IMMEDIATE SOURCE:

(B) CLONE: pARC0468 PBP 1B QQLL

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 8:

Met Ala Gly Asn Asp Arg Glu Pro Ile Gly Arg Lys Gly Lys Pro Thr  
 1 5 10 15  
 Arg Pro Val Lys Gln Lys Val Ser Arg Arg Arg Tyr Glu Asp Asp Asp  
 20 25 30  
 Asp Tyr Asp Asp Tyr Asp Asp Tyr Glu Asp Glu Glu Pro Met Pro Arg  
 35 40 45  
 Lys Gly Lys Gly Lys Gly Lys Gly Arg Lys Pro Arg Gly Lys Arg Gly  
 50 55 60  
 Trp Leu Trp Leu Leu Leu Lys Leu Ala Ile Val Phe Ala Val Leu Ile  
 65 70 75 80  
 Ala Ile Tyr Gly Val Tyr Leu Asp Gln Lys Ile Arg Ser Arg Ile Asp  
 85 90 95  
 Gly Lys Val Trp Gln Leu Ala Ala Ala Val Tyr Gly Arg Met Val Asn  
 100 105 110  
 Leu Glu Pro Asp Met Thr Ile Ser Lys Asn Glu Met Val Lys Leu Leu  
 115 120 125  
 Glu Ala Thr Gln Tyr Arg Gln Val Ser Lys Met Thr Arg Pro Gly Glu  
 130 135 140  
 Phe Thr Val Gln Ala Asn Ser Ile Glu Met Ile Arg Arg Pro Phe Asp  
 145 150 155 160  
 Phe Pro Asp Ser Lys Glu Gly Gln Val Arg Ala Arg Leu Thr Phe Asp  
 165 170 175  
 Gly Asp His Leu Ala Thr Ile Val Asn Met Glu Asn Asn Arg Gln Phe  
 180 185 190  
 Gly Phe Phe Arg Leu Asp Pro Arg Leu Ile Thr Met Ile Ser Ser Pro  
 195 200 205  
 Asn Gly Glu Gln Arg Leu Phe Val Pro Arg Ser Gly Phe Pro Asp Leu  
 210 215 220  
 Leu Val Asp Thr Leu Leu Ala Thr Glu Asp Arg His Phe Tyr Glu His  
 225 230 235 240  
 Asp Gly Ile Ser Leu Tyr Ser Ile Gly Arg Ala Val Leu Ala Asn Leu  
 245 250 255  
 Thr Ala Gly Arg Thr Val Gln Gly Ala Ser Thr Leu Thr Leu Leu  
 260 265 270

Val Lys Asn L u Phe Leu Ser S r Glu Arg Ser Tyr Trp Arg Lys Ala  
 275 280 285  
 Asn Glu Ala Tyr Met Ala Leu Ile Met Asp Ala Arg Tyr Ser Lys Asp  
 290 295 300  
 Arg Ile Leu Glu Leu Tyr Met Asn Glu Val Tyr Leu Gly Gln Ser Gly  
 305 310 315 320  
 Asp Asn Glu Ile Arg Gly Phe Pro Leu Ala Ser Leu Tyr Tyr Phe Gly  
 325 330 335  
 Arg Pro Val Glu Glu Leu Ser Leu Asp Gln Gln Ala Leu Leu Val Gly  
 340 345 350  
 Met Val Lys Gly Ala Ser Ile Tyr Asn Pro Trp Arg Asn Pro Lys Leu  
 355 360 365  
 Ala Leu Glu Arg Arg Asn Leu Val Leu Arg Leu Leu Gln Gln Gln Gln  
 370 375 380  
 Ile Ile Asp Gln Glu Leu Tyr Asp Met Leu Ser Ala Arg Pro Leu Gly  
 385 390 395 400  
 Val Gln Pro Arg Gly Gly Val Ile Ser Pro Gln Pro Ala Phe Met Gln  
 405 410 415  
 Leu Val Arg Gln Glu Leu Gln Ala Lys Leu Gly Asp Lys Val Lys Asp  
 420 425 430  
 Leu Ser Gly Val Lys Ile Phe Thr Thr Phe Asp Ser Val Ala Gln Asp  
 435 440 445  
 Ala Ala Glu Lys Ala Ala Val Glu Gly Ile Pro Ala Leu Lys Lys Gln  
 450 455 460  
 Arg Lys Leu Ser Asp Leu Glu Thr Ala Ile Val Val Val Asp Arg Phe  
 465 470 475 480  
 Ser Gly Glu Val Arg Ala Met Val Gly Gly Ser Glu Pro Gln Phe Ala  
 485 490 495  
 Gly Tyr Asn Arg Ala Met Gln Ala Arg Arg Ser Ile Gly Ser Leu Ala  
 500 505 510  
 Lys Pro Ala Thr Tyr Leu Thr Ala Leu Ser Gln Pro Lys Ile Tyr Arg  
 515 520 525  
 Leu Asn Thr Trp Ile Ala Asp Ala Pro Ile Ala Leu Arg Gln Pro Asn  
 530 535 540  
 Gly Gln Val Trp Ser Pro Gln Asn Asp Asp Arg Arg Tyr Ser Glu Ser  
 545 550 555 560  
 Gly Arg Val Met Leu Val Asp Ala Leu Thr Arg Ser Met Asn Val Pro  
 565 570 575  
 Thr Val Asn Leu Gly Met Ala Leu Gly Leu Pro Ala Val Thr Glu Thr  
 580 585 590  
 Trp Ile Lys Leu Gly Val Pro Lys Asp Gln Leu His Pro Val Pro Ala  
 595 600 605  
 Met Leu Leu Gly Ala L u Asn Leu Thr Pro Ile Glu Val Ala Gln Ala  
 610 615 620  
 Phe Gln Thr Ile Ala Ser Gly Gly Asn Arg Ala Pro Leu Ser Ala Leu  
 625 630 635 640

Arg Ser Val Ile Ala Glu Asp Gly Lys Val Leu Tyr Gln Ser Phe Pro  
 645 650 655  
 Gln Ala Glu Arg Ala Val Pro Ala Gln Ala Ala Tyr L u Thr Leu Trp  
 660 665 670  
 Thr Met Gln Gln Val Val Gln Arg Gly Thr Gly Arg Gln Leu Gly Ala  
 675 680 685  
 Lys Tyr Pro Asn Leu His Leu Ala Gly Lys Thr Gly Thr Thr Asn Asn  
 690 695 700  
 Asn Val Asp Thr Trp Phe Ala Gly Ile Asp Gly Ser Thr Val Thr Ile  
 705 710 715 720  
 Thr Trp Val Gly Arg Asp Asn Asn Gln Pro Thr Lys Leu Tyr Gly Ala  
 725 730 735  
 Ser Gly Ala Met Ser Ile Tyr Gln Arg Tyr Leu Ala Asn Gln Thr Pro  
 740 745 750  
 Thr Pro Leu Asn Leu Val Pro Pro Glu Asp Ile Ala Asp Met Gly Val  
 755 760 765  
 Asp Tyr Asp Gly Asn Phe Val Cys Ser Gly Gly Met Arg Ile Leu Pro  
 770 775 780  
 Val Trp Thr Ser Asp Pro Gln Ser Leu Cys Gln Gln Ser Glu Met Gln  
 785 790 795 800  
 Gln Gln Pro Ser Gly Asn Pro Phe Asp Gln Ser Ser Gln Pro Gln Gln  
 805 810 815  
 Gln Pro Gln Gln Gln Pro Ala Gln Gln Glu Gln Lys Asp Ser Asp Gly  
 820 825 830  
 Val Ala Gly Trp Ile Lys Asp Met Phe Gly Ser Asn  
 835 840

## (2) INFORMATION FOR SEQ ID NO: 9:

- (i) SEQUENCE CHARACTERISTICS:
  - (A) LENGTH: 836 amino acids
  - (B) TYPE: amino acid
  - (D) TOPOLOGY: linear
- (ii) MOLECULE TYPE: peptide
- (vi) ORIGINAL SOURCE:
  - (A) ORGANISM: Escherichia coli
- (vii) IMMEDIATE SOURCE:
  - (B) CLONE: pARC0469 PBP 1B del 8

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 9:

Met Ala Gly Asn Asp Arg Glu Pro Ile Gly Arg Lys Gly Lys Pro Thr  
 1 5 10 15  
 Arg Pro Val Lys Gln Lys Val Ser Arg Arg Arg Tyr Glu Asp Asp  
 20 25 30  
 Asp Tyr Asp Asp Tyr Asp Asp Tyr Glu Asp Glu Glu Pro Met Pro Arg  
 35 40 45  
 Lys Gly Lys Gly Lys Gly Lys Gly Arg Lys Pro Arg Gly Lys Arg Gly  
 50 55 60

Trp Leu Trp Leu L u Leu Lys Leu Ala Il Val Phe Ala Val Leu Ile  
 65 70 75 80  
 Ala Ile Tyr Gly Val Tyr Leu Asp Gln Lys Ile Arg Ser Arg Ile Asp  
 85 90 95  
 Gly Lys Val Trp Gln Leu Ala Ala Ala Val Tyr Gly Arg Met Val Asn  
 100 105 110  
 Leu Glu Pro Asp Met Thr Ile Ser Lys Asn Glu Met Val Lys Leu Leu  
 115 120 125  
 Glu Ala Thr Gln Tyr Arg Gln Val Ser Lys Met Thr Arg Pro Gly Glu  
 130 135 140  
 Phe Thr Val Gln Ala Asn Ser Ile Glu Met Ile Arg Arg Pro Phe Asp  
 145 150 155 160  
 Phe Pro Asp Ser Lys Glu Gly Gln Val Arg Ala Arg Leu Thr Phe Asp  
 165 170 175  
 Gly Asp His Leu Ala Thr Ile Val Asn Met Glu Asn Asn Arg Gln Phe  
 180 185 190  
 Gly Phe Phe Arg Leu Asp Pro Arg Leu Ile Thr Met Ile Ser Ser Pro  
 195 200 205  
 Asn Gly Glu Gln Arg Leu Phe Val Pro Arg Ser Gly Phe Pro Asp Leu  
 210 215 220  
 Leu Val Asp Thr Leu Leu Ala Thr Glu Asp Arg His Phe Tyr Glu His  
 225 230 235 240  
 Asp Gly Ile Ser Leu Tyr Ser Ile Gly Arg Ala Val Leu Ala Asn Leu  
 245 250 255  
 Thr Ala Gly Arg Thr Val Gln Leu Val Lys Asn Leu Phe Leu Ser Ser  
 260 265 270  
 Glu Arg Ser Tyr Trp Arg Lys Ala Asn Glu Ala Tyr Met Ala Leu Ile  
 275 280 285  
 Met Asp Ala Arg Tyr Ser Lys Asp Arg Ile Leu Glu Leu Tyr Met Asn  
 290 295 300  
 Glu Val Tyr Leu Gly Gln Ser Gly Asp Asn Glu Ile Arg Gly Phe Pro  
 305 310 315 320  
 Leu Ala Ser Leu Tyr Tyr Phe Gly Arg Pro Val Glu Glu Leu Ser Leu  
 325 330 335  
 Asp Gln Gln Ala Leu Leu Val Gly Met Val Lys Gly Ala Ser Ile Tyr  
 340 345 350  
 Asn Pro Trp Arg Asn Pro Lys Leu Ala Leu Glu Arg Arg Asn Leu Val  
 355 360 365  
 Leu Arg Leu Leu Gln Gln Gln Gln Ile Ile Asp Gln Glu Leu Tyr Asp  
 370 375 380  
 Met Leu Ser Ala Arg Pro Leu Gly Val Gln Pro Arg Gly Gly Val Ile  
 385 390 395 400  
 Ser Pro Gln Pro Ala Phe Met Gln Leu Val Arg Gln Glu Leu Gln Ala  
 405 410 415  
 Lys L u Gly Asp Lys Val Lys Asp Leu Ser Gly Val Lys Ile Phe Thr  
 420 425 430



Thr Phe Asp Ser Val Ala Gln Asp Ala Ala Glu Lys Ala Ala Val Glu  
 435 440 445  
 Gly Il Pro Ala Leu Lys Lys Gln Arg Lys Leu Ser Asp Leu Glu Thr  
 450 455 460  
 Ala Ile Val Val Val Asp Arg Phe Ser Gly Glu Val Arg Ala Met Val  
 465 470 475 480  
 Gly Gly Ser Glu Pro Gln Phe Ala Gly Tyr Asn Arg Ala Met Gln Ala  
 485 490 495  
 Arg Arg Ser Ile Gly Ser Leu Ala Lys Pro Ala Thr Tyr Leu Thr Ala  
 500 505 510  
 Leu Ser Gln Pro Lys Ile Tyr Arg Leu Asn Thr Trp Ile Ala Asp Ala  
 515 520 525  
 Pro Ile Ala Leu Arg Gln Pro Asn Gly Gln Val Trp Ser Pro Gln Asn  
 530 535 540  
 Asp Asp Arg Arg Tyr Ser Glu Ser Gly Arg Val Met Leu Val Asp Ala  
 545 550 555 560  
 Leu Thr Arg Ser Met Asn Val Pro Thr Val Asn Leu Gly Met Ala Leu  
 565 570 575  
 Gly Leu Pro Ala Val Thr Glu Thr Trp Ile Lys Leu Gly Val Pro Lys  
 580 585 590  
 Asp Gln Leu His Pro Val Pro Ala Met Leu Leu Gly Ala Leu Asn Leu  
 595 600 605  
 Thr Pro Ile Glu Val Ala Gln Ala Phe Gln Thr Ile Ala Ser Gly Gly  
 610 615 620  
 Asn Arg Ala Pro Leu Ser Ala Leu Arg Ser Val Ile Ala Glu Asp Gly  
 625 630 635 640  
 Lys Val Leu Tyr Gln Ser Phe Pro Gln Ala Glu Arg Ala Val Pro Ala  
 645 650 655  
 Gln Ala Ala Tyr Leu Thr Leu Trp Thr Met Gln Gln Val Val Gln Arg  
 660 665 670  
 Gly Thr Gly Arg Gln Leu Gly Ala Lys Tyr Pro Asn Leu His Leu Ala  
 675 680 685  
 Gly Lys Thr Gly Thr Thr Asn Asn Asn Val Asp Thr Trp Phe Ala Gly  
 690 695 700  
 Ile Asp Gly Ser Thr Val Thr Ile Thr Trp Val Gly Arg Asp Asn Asn  
 705 710 715 720  
 Gln Pro Thr Lys Leu Tyr Gly Ala Ser Gly Ala Met Ser Ile Tyr Gln  
 725 730 735  
 Arg Tyr Leu Ala Asn Gln Thr Pro Thr Pro Leu Asn Leu Val Pro Pro  
 740 745 750  
 Glu Asp Ile Ala Asp Met Gly Val Asp Tyr Asp Gly Asn Phe Val Cys  
 755 760 765  
 Ser Gly Gly Met Arg Ile Leu Pro Val Trp Thr Ser Asp Pro Gln Ser  
 770 775 780  
 Leu Cys Gln G Ser Glu Met Gln Gln Gln Pro Ser Gly Asn Pro Phe  
 785 790 795 800

|            |           |            |            |            |            |            |            |            |            |            |            |            |            |                   |
|------------|-----------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|-------------------|
| Met<br>1   | Lys       | Phe        | Val        | Lys<br>5   | Tyr        | Phe        | Leu        | Ile<br>10  | Ala        | Val        | Cys        | Cys        | Ile<br>15  | Leu               |
| Leu        | Gly       | Ala        | Gly<br>20  | Ser        | Ile        | Tyr        | Gly        | Leu<br>25  | Tyr        | Arg        | Tyr        | Ile        | Glu<br>30  | Pro<br>Gln        |
| Leu        | Pro       | Asp<br>35  | Val        | Ala        | Thr        | Leu        | Lys<br>40  | Asp        | Val        | Arg        | Leu        | Gln<br>45  | Ile        | Pro<br>Met        |
| Gln        | Ile<br>50 | Tyr        | Ser        | Ala        | Asp        | Gly<br>55  | Glu        | Leu        | Ile        | Ala        | Gln<br>60  | Tyr        | Gly        | Glu<br>Lys        |
| Arg<br>65  | Arg       | Ile        | Pro        | Val        | Thr<br>70  | Leu        | Asp        | Gln        | Ile        | Pro<br>75  | Pro        | Glu        | Met        | Val<br>Lys<br>80  |
| Ala        | Phe       | Ile        | Ala        | Thr<br>85  | Glu        | Asp        | Ser        | Arg        | Phe<br>90  | Tyr        | Glu        | His        | His        | Gly<br>Val<br>95  |
| Asp        | Pro       | Val        | Gly<br>100 | Ile        | Phe        | Arg        | Ala        | Ala<br>105 | Ser        | Val        | Ala        | Leu        | Phe<br>110 | Ser<br>Gly        |
| His        | Ala       | Ser<br>115 | Gln        | Gly        | Ala        | Ser        | Thr<br>120 | Ile        | Thr        | Ala        | Ala        | Leu<br>125 | Ala        | Arg<br>Asn        |
| Phe<br>130 | Phe       | Leu        | Ser        | Pro        | Glu        | Arg<br>135 | Thr        | Leu        | Met        | Arg        | Lys<br>140 | Ile        | Lys        | Glu<br>Val        |
| Phe<br>145 | Leu       | Ala        | Ile        | Arg        | Ile<br>150 | Glu        | Gln        | Leu        | Leu        | Thr<br>155 | Lys        | Asp        | Glu        | Ile<br>Leu<br>160 |
| Glu        | Leu       | Tyr        | Leu        | Asn<br>165 | Lys        | Ile        | Tyr        | Leu        | Gly<br>170 | Tyr        | Arg        | Ala        | Tyr        | Gly<br>Val<br>175 |
| Gly        | Ala       | Ala        | Ala<br>180 | Gln        | Val        | Tyr        | Phe        | Gly<br>185 | Lys        | Thr        | Val        | Asp        | Gln<br>190 | Leu<br>Thr        |
| Leu        | Asn       | Glu<br>195 | Met        | Ala        | Val        | Ile        | Ala<br>200 | Gly        | Leu        | Pro        | Lys        | Ala<br>205 | Pro        | Ser<br>Thr        |
| Phe<br>210 | Asn       | Pro        | L u        | Tyr        | Ser        | Met<br>215 | Asp        | Arg        | Ala        | Val        | Ala<br>220 | Arg        | Arg        | Asn<br>Val        |

Val L u Ser Arg Met Leu Asp Glu Gly Tyr Il Thr Gln Gln Gln Phe  
 225 230 235 240  
 Asp Gln Thr Arg Thr Glu Ala Ile Asn Ala Asn Tyr His Ala Pro Glu  
 245 250 255  
 Ile Ala Phe Ser Ala Pro Tyr Leu Ser Glu Met Val Arg Gln Glu Met  
 260 265 270  
 Tyr Asn Arg Tyr Gly Glu Ser Ala Tyr Glu Asp Gly Tyr Arg Ile Tyr  
 275 280 285  
 Thr Thr Ile Thr Arg Lys Val Gln Gln Ala Ala Gln Gln Ala Val Arg  
 290 295 300  
 Asn Asn Val Leu Asp Tyr Asp Met Arg His Gly Tyr Arg Gly Pro Ala  
 305 310 315 320  
 Asn Val Leu Trp Lys Val Gly Glu Ser Ala Trp Asp Asn Asn Lys Ile  
 325 330 335  
 Thr Asp Thr Leu Lys Ala Leu Pro Thr Tyr Gly Pro Leu Leu Pro Ala  
 340 345 350  
 Ala Val Thr Ser Ala Asn Pro Gln Gln Ala Thr Ala Met Leu Ala Asp  
 355 360 365  
 Gly Ser Thr Val Ala Leu Ser Met Glu Gly Val Arg Trp Ala Arg Pro  
 370 375 380  
 Tyr Arg Ser Asp Thr Gln Gln Gly Pro Thr Pro Arg Lys Val Thr Asp  
 385 390 395 400  
 Val Leu Gln Thr Gly Gln Gln Ile Trp Val Arg Gln Val Gly Asp Ala  
 405 410 415  
 Trp Trp Leu Ala Gln Val Pro Glu Val Asn Ser Ala Leu Val Ser Ile  
 420 425 430  
 Asn Pro Gln Asn Gly Ala Val Met Ala Leu Val Gly Gly Phe Asp Phe  
 435 440 445  
 Asn Gln Ser Lys Phe Asn Arg Ala Thr Gln Ala Leu Arg Gln Val Gly  
 450 455 460  
 Ser Asn Ile Lys Pro Phe Leu Tyr Thr Ala Ala Met Asp Lys Gly Leu  
 465 470 475 480  
 Thr Leu Ala Ser Met Leu Asn Asp Val Pro Ile Ser Arg Trp Asp Ala  
 485 490 495  
 Ser Ala Gly Ser Asp Trp Gln Pro Lys Asn Ser Pro Pro Gln Tyr Ala  
 500 505 510  
 Gly Pro Ile Arg Leu Arg Gln Gly Leu Gly Gln Ser Lys Asn Val Val  
 515 520 525  
 Met Val Arg Ala Met Arg Ala Met Gly Val Asp Tyr Ala Ala Glu Tyr  
 530 535 540  
 Leu Gln Arg Phe Gly Phe Pro Ala Gln Asn Ile Val His Thr Glu Ser  
 545 550 555 560  
 Leu Ala Leu Gly Ser Ala Ser Phe Thr Pro Met Gln Val Ala Arg Gly  
 565 570 575  
 Tyr Ala Val Met Ala Asn Gly Gly Phe Leu Val Asp Pro Trp Phe Il  
 580 585 590

S r Lys Il Glu Asn Asp Gln Gly Gly Val Ile Phe Glu Ala Lys Pro  
 595 600 605  
 Lys Val Ala Cys Pro Glu Cys Asp Ile Pro Val Ile Tyr Gly Asp Thr  
 610 615 620  
 Gln Lys Ser Asn Val Leu Glu Asn Asn Asp Val Glu Asp Val Ala Ile  
 625 630 635 640  
 Ser Arg Glu Gln Gln Asn Val Ser Val Pro Met Pro Gln Leu Glu Gln  
 645 650 655  
 Ala Asn Gln Ala Leu Val Ala Lys Thr Gly Ala Gln Glu Tyr Ala Pro  
 660 665 670  
 His Val Ile Asn Thr Pro Leu Ala Phe Leu Ile Lys Ser Ala Leu Asn  
 675 680 685  
 Thr Asn Ile Phe Gly Glu Pro Gly Trp Gln Gly Thr Gly Trp Arg Ala  
 690 695 700  
 Gly Arg Asp Leu Gln Arg Arg Asp Ile Gly Gly Lys Thr Gly Thr Thr  
 705 710 715 720  
 Asn Ser Ser Lys Asp Ala Trp Phe Ser Gly Tyr Gly Pro Gly Val Val  
 725 730 735  
 Thr Ser Val Trp Ile Gly Phe Asp Asp His Arg Arg Asn Leu Gly His  
 740 745 750  
 Thr Thr Ala Ser Gly Ala Ile Lys Asp Gln Ile Ser Gly Tyr Glu Gly  
 755 760 765  
 Gly Ala Lys Ser Ala Gln Pro Ala Trp Asp Ala Tyr Met Lys Ala Val  
 770 775 780  
 Leu Glu Gly Val Pro Glu Gln Pro Leu Thr Pro Pro Pro Gly Ile Val  
 785 790 795 800  
 Thr Val Asn Ile Asp Arg Ser Thr Gly Gln Leu Ala Asn Gly Gly Asn  
 805 810 815  
 Ser Arg Glu Glu Tyr Phe Ile Glu Gly Thr Gln Pro Thr Gln Gln Ala  
 820 825 830  
 Val His Glu Val Gly Thr Thr Ile Ile Asp Asn Gly Glu Ala Gln Glu  
 835 840 845  
 Leu Phe  
 850

## (2) INFORMATION FOR SEQ ID NO: 11:

- (i) SEQUENCE CHARACTERISTICS:
  - (A) LENGTH: 553 amino acids
  - (B) TYPE: amino acid
  - (D) TOPOLOGY: linear
- (ii) MOLECULE TYPE: peptide
- (vi) ORIGINAL SOURCE:
  - (A) ORGANISM: Escherichia coli
- (vii) IMMEDIATE SOURCE:
  - (B) CLONE: pARC 0592 truncated PBP 1B

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 11:

|     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| Met | Ala | Gly | Asn | Asp | Arg | Glu | Pro | Ile | Gly | Arg | Lys | Gly | Lys | Pro | Thr | 1   | 5   | 10  | 15  |
| Arg | Pro | Val | Lys | Gln | Lys | Val | Ser | Arg | Arg | Arg | Tyr | Glu | Asp | Asp | Asp | 20  | 25  | 30  |     |
| Asp | Tyr | Asp | Asp | Tyr | Asp | Asp | Tyr | Glu | Asp | Glu | Glu | Pro | Met | Pro | Arg | 35  | 40  | 45  |     |
| Lys | Gly | Lys | Gly | Lys | Gly | Lys | Gly | Arg | Lys | Pro | Arg | Gly | Lys | Arg | Gly | 50  | 55  | 60  |     |
| Trp | Leu | Trp | Leu | Leu | Leu | Lys | Leu | Ala | Ile | Val | Phe | Ala | Val | Leu | Ile | 65  | 70  | 75  | 80  |
| Ala | Ile | Tyr | Gly | Val | Tyr | Leu | Asp | Gln | Lys | Ile | Arg | Ser | Arg | Ile | Asp | 85  | 90  | 95  |     |
| Gly | Lys | Val | Trp | Gln | Leu | Ala | Ala | Ala | Val | Tyr | Gly | Arg | Met | Val | Asn | 100 | 105 | 110 |     |
| Leu | Glu | Pro | Asp | Met | Thr | Ile | Ser | Lys | Asn | Glu | Met | Val | Lys | Leu | Leu | 115 | 120 | 125 |     |
| Glu | Ala | Thr | Gln | Tyr | Arg | Gln | Val | Ser | Lys | Met | Thr | Arg | Pro | Gly | Glu | 130 | 135 | 140 |     |
| Phe | Thr | Val | Gln | Ala | Asn | Ser | Ile | Glu | Met | Ile | Arg | Arg | Pro | Phe | Asp | 145 | 150 | 155 | 160 |
| Phe | Pro | Asp | Ser | Lys | Glu | Gly | Gln | Val | Arg | Ala | Arg | Leu | Thr | Phe | Asp | 165 | 170 | 175 |     |
| Gly | Asp | His | Leu | Ala | Thr | Ile | Val | Asn | Met | Glu | Asn | Asn | Arg | Gln | Phe | 180 | 185 | 190 |     |
| Gly | Phe | Phe | Arg | Leu | Asp | Pro | Arg | Leu | Ile | Thr | Met | Ile | Ser | Ser | Pro | 195 | 200 | 205 |     |
| Asn | Gly | Glu | Gln | Arg | Leu | Phe | Val | Pro | Arg | Ser | Gly | Phe | Pro | Asp | Leu | 210 | 215 | 220 |     |
| Leu | Val | Asp | Thr | Leu | Leu | Ala | Thr | Glu | Asp | Arg | His | Phe | Tyr | Glu | His | 225 | 230 | 235 | 240 |
| Asp | Gly | Ile | Ser | Leu | Tyr | Ser | Ile | Gly | Arg | Ala | Val | Leu | Ala | Asn | Leu | 245 | 250 | 255 |     |
| Thr | Ala | Gly | Arg | Thr | Val | Gln | Gly | Ala | Ser | Thr | Leu | Thr | Gln | Gln | Leu | 260 | 265 | 270 |     |
| Val | Lys | Asn | Leu | Phe | Leu | Ser | Ser | Glu | Arg | Ser | Tyr | Trp | Arg | Lys | Ala | 275 | 280 | 285 |     |
| Asn | Glu | Ala | Tyr | Met | Ala | Leu | Ile | Met | Asp | Ala | Arg | Tyr | Ser | Lys | Asp | 290 | 295 | 300 |     |
| Arg | Ile | Leu | Glu | Leu | Tyr | Met | Asn | Glu | Val | Tyr | Leu | Gly | Gln | Ser | Gly | 305 | 310 | 315 | 320 |
| Asp | Asn | Glu | Ile | Arg | Gly | Phe | Pro | Leu | Ala | Ser | Leu | Tyr | Tyr | Phe | Gly | 325 | 330 | 335 |     |
| Arg | Pro | Val | Glu | Glu | Leu | Ser | Leu | Asp | Gln | Gln | Ala | Leu | Leu | Val | Gly | 340 | 345 | 350 |     |

[illegible]

## (2) INFORMATION FOR SEQ ID NO: 12:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 532 amino acids

(B) TYPE: amino acid

(D) TOPOLOGY: linear

## (ii) MOLECULE TYPE: peptide

(vi) ORIGINAL SOURCE:

(A) ORGANISM: *Escherichia coli*

(vii) IMMEDIATE SOURCE:

(B) CLONE: pARC 0593      truncated soluble PBP 1B

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 12:

|     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |     |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| Met | Ala | Gly | Asn | Asp | Arg | Glu | Pro | Ile | Gly | Arg | Lys | Gly | Lys | Pro | Thr |
| 1   |     |     |     | 5   |     |     |     |     | 10  |     |     |     |     | 15  |     |
| Arg | Pro | Val | Lys | Gln | Lys | Val | Ser | Arg | Arg | Tyr | Glu | Asp | Asp | Asp |     |
|     |     |     | 20  |     |     |     | 25  |     |     |     |     | 30  |     |     |     |
| Asp | Tyr | Asp | Asp | Tyr | Asp | Asp | Tyr | Glu | Asp | Glu | Glu | Pro | Met | Pro | Arg |
|     |     | 35  |     |     |     |     | 40  |     |     |     |     | 45  |     |     |     |
| Lys | Gly | Lys | Gly | Lys | Gly | Lys | Gly | Arg | Lys | Pro | Arg | Gly | Lys | Arg | Gly |
|     | 50  |     |     |     |     | 55  |     |     |     |     | 60  |     |     |     |     |

Ser Ile Asp Gln Lys Ile Arg Ser Arg Ile Asp Gly Lys Val Trp Gln  
 65 70 75 80  
 Leu Ala Ala Ala Val Tyr Gly Arg Met Val Asn Leu Glu Pro Asp Met  
 85 90 95  
 Thr Ile Ser Lys Asn Glu Met Val Lys Leu Leu Glu Ala Thr Gln Tyr  
 100 105 110  
 Arg Gln Val Ser Lys Met Thr Arg Pro Gly Glu Phe Thr Val Gln Ala  
 115 120 125  
 Asn Ser Ile Glu Met Ile Arg Arg Pro Phe Asp Phe Pro Asp Ser Lys  
 130 135 140  
 Glu Gly Gln Val Arg Ala Arg Leu Thr Phe Asp Gly Asp His Leu Ala  
 145 150 155 160  
 Thr Ile Val Asn Met Glu Asn Asn Arg Gln Phe Gly Phe Phe Arg Leu  
 165 170 175  
 Asp Pro Arg Leu Ile Thr Met Ile Ser Ser Pro Asn Gly Glu Gln Arg  
 180 185 190  
 Leu Phe Val Pro Arg Ser Gly Phe Pro Asp Leu Leu Val Asp Thr Leu  
 195 200 205  
 Leu Ala Thr Glu Asp Arg His Phe Tyr Glu His Asp Gly Ile Ser Leu  
 210 215 220  
 Tyr Ser Ile Gly Arg Ala Val Leu Ala Asn Leu Thr Ala Gly Arg Thr  
 225 230 235 240  
 Val Gln Gly Ala Ser Thr Leu Thr Gln Gln Leu Val Lys Asn Leu Phe  
 245 250 255  
 Leu Ser Ser Glu Arg Ser Tyr Trp Arg Lys Ala Asp Glu Ala Tyr Met  
 260 265 270  
 Ala Leu Ile Met Asp Ala Arg Tyr Ser Lys Asp Arg Ile Leu Glu Leu  
 275 280 285  
 Tyr Met Asn Glu Val Tyr Leu Gly Gln Ser Gly Asp Asn Glu Ile Arg  
 290 295 300  
 Gly Phe Pro Leu Ala Ser Leu Tyr Tyr Phe Gly Arg Pro Val Glu Glu  
 305 310 315 320  
 Leu Ser Leu Asp Gln Gln Ala Leu Leu Val Gly Met Val Lys Gly Ala  
 325 330 335  
 Ser Ile Tyr Asn Pro Trp Arg Asn Pro Lys Leu Ala Leu Glu Arg Arg  
 340 345 350  
 Asn Leu Val Leu Arg Leu Leu Gln Gln Gln Ile Ile Asp Gln Glu  
 355 360 365  
 Leu Tyr Asp Met Leu Ser Ala Arg Pro Leu Gly Val Gln Pro Arg Gly  
 370 375 380  
 Gly Val Ile Ser Pro Gln Pro Ala Phe Met Gln Leu Val Arg Gln Glu  
 385 390 395 400  
 Leu Gln Ala Lys Leu Gly Asp Lys Val Lys Asp Leu Ser Gly Val Lys  
 405 410 415  
 Ile Phe Thr Thr Phe Asp Ser Val Ala Gln Asp Ala Ala Glu Lys Ala  
 420 425 430

Ala Val Glu Gly Ile Pro Ala Leu Lys Lys Gln Arg Lys Leu Ser Asp  
 435 440 445

Leu Glu Thr Ala Ile Val Val Val Asp Arg Phe Ser Gly Glu Val Arg  
 450 455 460

Ala Met Val Gly Gly Ser Glu Pro Gln Phe Ala Gly Tyr Asn Arg Ala  
 465 470 475 480

Met Gln Ala Arg Arg Ser Ile Gly Ser Leu Ala Lys Pro Ala Thr Tyr  
 485 490 495

Leu Thr Ala Leu Ser Gln Pro Lys Ile Tyr Arg Leu Asn Thr Trp Ile  
 500 505 510

Ala Asp Ala Pro Ile Ala Leu Arg Gln Pro Asn Gly Gln Val Trp Ser  
 515 520 525

Pro Gln Asn Asp  
 530

## (2) INFORMATION FOR SEQ ID NO: 13:

- (i) SEQUENCE CHARACTERISTICS:
  - (A) LENGTH: 159 amino acids
  - (B) TYPE: amino acid
  - (C) STRANDEDNESS: double
  - (D) TOPOLOGY: linear
- (ii) MOLECULE TYPE: DNA (genomic)
- (vi) ORIGINAL SOURCE:
  - (A) ORGANISM: Escherichia coli
- (vii) IMMEDIATE SOURCE:
  - (B) CLONE: pARC 0392
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 13:

Gly Glu Gln Arg Leu Phe Val Pro Arg Ser Gly Phe Pro Asp Leu Leu  
 1 5 10 15

Val Asp Thr Leu Leu Ala Thr Glu Asp Arg His Phe Tyr Glu His Asp  
 20 25 30

Gly Ile Ser Leu Tyr Ser Ile Gly Arg Ala Val Leu Ala Asn Leu Thr  
 35 40 45

Ala Gly Arg Thr Val Gln Gly Ala Ser Thr Leu Thr Gln Gln Leu Val  
 50 55 60

Lys Asn Leu Phe Leu Ser Ser Glu Arg Ser Tyr Trp Arg Lys Ala Asn  
 65 70 75 80

Glu Ala Tyr Met Ala Leu Ile Met Asp Ala Arg Tyr Ser Lys Asp Arg  
 85 90 95

Ile Leu Glu Leu Tyr Met Asn Glu Val Tyr Leu Gly Gln Ser Gly Asp  
 100 105 110

Asn Glu Ile Arg Gly Phe Pro Leu Ala Ser Leu Tyr Tyr Phe Gly Arg  
 115 120 125

Pro Val Glu Glu L u Ser Leu Asp Gln Gln Ala Leu Leu Val Gly Met  
 130 135 140



Val Lys Gly Ala Ser Ile Tyr Asn Pro Trp Arg Asn Pro Lys Leu  
145 150 155

## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                                                                                                                                                                                                     |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>A.</b> The indications made below relate to the microorganism referred to in the description<br>on page <u>19</u> , line <u>18</u>                                                                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                                                     |
| <b>B. IDENTIFICATION OF DEPOSIT</b> <span style="float: right;">Further deposits are identified on an additional sheet <input type="checkbox"/></span>                                                                                                                                                                                                                                                                                                             |                                                                                                                                                                                                                                                                                                                     |
| Name of depositary institution<br><b>The National Collections of Industrial and Marine Bacteria Limited (NCIMB)</b>                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                     |
| Address of depositary institution (including postal code and country)<br><br><b>23 St Machar Drive<br/>Aberdeen AB2 1RY<br/>Scotland, UK</b>                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                                     |
| Date of deposit<br><b>28 June 1994</b>                                                                                                                                                                                                                                                                                                                                                                                                                             | Accession Number<br><b>NCIMB 40666</b>                                                                                                                                                                                                                                                                              |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) <span style="float: right;">This information is continued on an additional sheet <input type="checkbox"/></span>                                                                                                                                                                                                                                                                                  |                                                                                                                                                                                                                                                                                                                     |
| In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state. |                                                                                                                                                                                                                                                                                                                     |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                  |                                                                                                                                                                                                                                                                                                                     |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                                                                                                                                                                                                     |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                                     |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                               |                                                                                                                                                                                                                                                                                                                     |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                                                                                                                                                                                                     |
| <div style="border: 1px solid black; padding: 2px;">For receiving Office use only</div> <div style="border: 1px solid black; padding: 2px;"><input checked="" type="checkbox"/> This sheet was received with the international application</div> <div style="border: 1px solid black; padding: 2px;">Authorized officer </div>                                                                                                                                     | <div style="border: 1px solid black; padding: 2px;">For International Bureau use only</div> <div style="border: 1px solid black; padding: 2px;"><input type="checkbox"/> This sheet was received by the International Bureau on:</div> <div style="border: 1px solid black; padding: 2px;">Authorized officer</div> |

## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                             |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| <b>A. The indications made below relate to the microorganism referred to in the description</b><br>on page <u>24</u> , line <u>10</u>                                                                                                                                                                                                                                                                                                                              |                                                                                                                                             |
| <b>B. IDENTIFICATION OF DEPOSIT</b>                                                                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                             |
| Further deposits are identified on an additional sheet <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                             |
| Name of depositary institution<br>The National Collections of Industrial and Marine Bacteria Limited (NCIMB)                                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                             |
| Address of depositary institution (including postal code and country)<br>23 St Machar Drive<br>Aberdeen AB2 1RY<br>Scotland, UK                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                             |
| Date of deposit<br>28 June 1994                                                                                                                                                                                                                                                                                                                                                                                                                                    | Accession Number<br>NCIMB 40667                                                                                                             |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) This information is continued on an additional sheet <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                     |                                                                                                                                             |
| In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state. |                                                                                                                                             |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                  |                                                                                                                                             |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                             |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                             |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                               |                                                                                                                                             |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                             |
| For receiving Office use only<br><input checked="" type="checkbox"/> This sheet was received with the international application<br>Authorized officer <u>C. H. W.</u>                                                                                                                                                                                                                                                                                              | For International Bureau use only<br><input type="checkbox"/> This sheet was received by the International Bureau on:<br>Authorized officer |

## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                        |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------|
| <b>A.</b> The indications made below relate to the microorganism referred to in the description<br>on page <u>27</u> , line <u>2</u>                                                                                                                                                                                                                                                                                                                               |                                        |
| <b>B. IDENTIFICATION OF DEPOSIT</b> Further deposits are identified on an additional sheet <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                                                |                                        |
| Name of depositary institution<br><b>The National Collections of Industrial and Marine Bacteria Limited (NCIMB)</b>                                                                                                                                                                                                                                                                                                                                                |                                        |
| Address of depositary institution (including postal code and country)<br><br><b>23 St Machar Drive<br/>Aberdeen AB2 1RY<br/>Scotland, UK</b>                                                                                                                                                                                                                                                                                                                       |                                        |
| Date of deposit<br><b>28 June 1994</b>                                                                                                                                                                                                                                                                                                                                                                                                                             | Accession Number<br><b>NCIMB 40665</b> |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) This information is continued on an additional sheet <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                     |                                        |
| In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state. |                                        |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                  |                                        |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                        |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                       |                                        |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                               |                                        |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                        |

|                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                  |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <div style="text-align: right; font-size: small;">For receiving Office use only</div> <div style="padding: 5px;"><input checked="" type="checkbox"/> This sheet was received with the international application</div> <div style="border-top: 1px solid black; padding-top: 5px;">Authorized officer <u>CSM</u></div> | <div style="text-align: right; font-size: small;">For International Bureau use only</div> <div style="padding: 5px;"><input type="checkbox"/> This sheet was received by the International Bureau on:</div> <div style="border-top: 1px solid black; padding-top: 5px;">Authorized officer</div> |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|
| <b>A.</b> The indications made below relate to the microorganism referred to in the description<br>on page <u>29</u> , line <u>19</u>                                                                                                                                                                                                                                                                                                                                     |                                 |
| <b>B. IDENTIFICATION OF DEPOSIT</b> <span style="float: right;">Further deposits are identified on an additional sheet <input type="checkbox"/></span>                                                                                                                                                                                                                                                                                                                    |                                 |
| Name of depositary institution<br><p style="text-align: center;">The National Collections of Industrial and Marine Bacteria Limited (NCIMB)</p>                                                                                                                                                                                                                                                                                                                           |                                 |
| Address of depositary institution (including postal code and country)<br><p>23 St Machar Drive<br/>         Aberdeen AB2 1RY<br/>         Scotland, UK</p>                                                                                                                                                                                                                                                                                                                |                                 |
| Date of deposit<br>28 June 1994                                                                                                                                                                                                                                                                                                                                                                                                                                           | Accession Number<br>NCIMB 40661 |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) <span style="float: right;">This information is continued on an additional sheet <input type="checkbox"/></span>                                                                                                                                                                                                                                                                                         |                                 |
| <p>In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state.</p> |                                 |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                         |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                              |                                 |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                                      |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |

|                                                                                                                                                                                                               |                                                                                                                                                                                        |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p style="text-align: center;">For receiving Office use only</p> <p><input checked="" type="checkbox"/> This sheet was received with the international application</p> <p>Authorized officer <u>CSind</u></p> | <p style="text-align: center;">For International Bureau use only</p> <p><input type="checkbox"/> This sheet was received by the International Bureau on:</p> <p>Authorized officer</p> |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                    |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------|
| <b>A.</b> The indications made below relate to the microorganism referred to in the description<br>on page <u>29</u> , line <u>23</u>                                                                                                                                                                                                                                                                                                                                     |                                                                    |
| <b>B. IDENTIFICATION OF DEPOSIT</b> <span style="float: right;">Further deposits are identified on an additional sheet <input type="checkbox"/></span>                                                                                                                                                                                                                                                                                                                    |                                                                    |
| Name of depositary institution<br><p style="text-align: center;">The National Collections of Industrial and Marine Bacteria Limited (NCIMB)</p>                                                                                                                                                                                                                                                                                                                           |                                                                    |
| Address of depositary institution (including postal code and country)<br><p style="text-align: center;">23 St Machar Drive<br/>Aberdeen AB2 1RY<br/>Scotland, UK</p>                                                                                                                                                                                                                                                                                                      |                                                                    |
| Date of deposit<br><p style="text-align: center;">28 June 1994</p>                                                                                                                                                                                                                                                                                                                                                                                                        | Accession Number<br><p style="text-align: center;">NCIMB 40662</p> |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) <span style="float: right;">This information is continued on an additional sheet <input type="checkbox"/></span>                                                                                                                                                                                                                                                                                         |                                                                    |
| <p>In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state.</p> |                                                                    |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                         |                                                                    |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                    |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                              |                                                                    |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                                      |                                                                    |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                    |

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <div style="border: 1px solid black; padding: 5px;"> <div style="display: flex; justify-content: space-between; align-items: center;"> <span>For receiving Office use only</span> <span>_____</span> </div> <div style="display: flex; align-items: center;"> <input checked="" type="checkbox"/> This sheet was received with the international application         </div> <div style="border-top: 1px solid black; padding-top: 5px;">           Authorized officer <span style="float: right; margin-right: 50px;"><i>CS7nd</i></span> </div> </div> | <div style="border: 1px solid black; padding: 5px;"> <div style="display: flex; justify-content: space-between; align-items: center;"> <span>For International Bureau use only</span> <span>_____</span> </div> <div style="display: flex; align-items: center;"> <input type="checkbox"/> This sheet was received by the International Bureau on:         </div> <div style="border-top: 1px solid black; padding-top: 5px;">           Authorized officer         </div> </div> |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|
| <b>A.</b> The indications made below relate to the microorganism referred to in the description<br>on page <u>29</u> , line <u>27</u>                                                                                                                                                                                                                                                                                                                                     |                                 |
| <b>B. IDENTIFICATION OF DEPOSIT</b> <div style="text-align: right;">Further deposits are identified on an additional sheet <input type="checkbox"/></div>                                                                                                                                                                                                                                                                                                                 |                                 |
| Name of depositary institution<br>The National Collections of Industrial and Marine Bacteria Limited (NCIMB)                                                                                                                                                                                                                                                                                                                                                              |                                 |
| Address of depositary institution (including postal code and country)<br><br>23 St Machar Drive<br>Aberdeen AB2 1RY<br>Scotland, UK                                                                                                                                                                                                                                                                                                                                       |                                 |
| Date of deposit<br>28 June 1994                                                                                                                                                                                                                                                                                                                                                                                                                                           | Accession Number<br>NCIMB 40663 |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) <div style="text-align: right;">This information is continued on an additional sheet <input type="checkbox"/></div>                                                                                                                                                                                                                                                                                      |                                 |
| <p>In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state.</p> |                                 |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                         |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                              |                                 |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                                      |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |

|                                                                                                                                                                                                                                                                                                                                                                                                                          |                                                                                                                                                                                                                                                                                                                                                               |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <div style="text-align: right;">For receiving Office use only</div> <div> <input checked="" type="checkbox"/> This sheet was received with the international application         </div> <div style="border-top: 1px solid black; height: 40px; display: flex; align-items: center;"> <div style="flex: 1;">Authorized officer</div> <div style="flex: 1; text-align: center; font-family: cursive;">C. Lind</div> </div> | <div style="text-align: right;">For International Bureau use only</div> <div> <input type="checkbox"/> This sheet was received by the International Bureau on:         </div> <div style="border-top: 1px solid black; height: 40px; display: flex; align-items: center;"> <div style="flex: 1;">Authorized officer</div> <div style="flex: 1;"></div> </div> |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                                                                             |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| <b>A.</b> The indications made below relate to the microorganism referred to in the description<br>on page <u>34</u> , line <u>3</u>                                                                                                                                                                                                                                                                                                                                      |                                                                                                                                             |
| <b>B. IDENTIFICATION OF DEPOSIT</b> <span style="float: right;">Further deposits are identified on an additional sheet <input type="checkbox"/></span>                                                                                                                                                                                                                                                                                                                    |                                                                                                                                             |
| Name of depositary institution<br>The National Collections of Industrial and Marine Bacteria Limited (NCIMB)                                                                                                                                                                                                                                                                                                                                                              |                                                                                                                                             |
| Address of depositary institution (including postal code and country)<br><br>23 St Machar Drive<br>Aberdeen AB2 1RY<br>Scotland, UK                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                             |
| Date of deposit<br>28 June 1994                                                                                                                                                                                                                                                                                                                                                                                                                                           | Accession Number<br>NCIMB 40668                                                                                                             |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) <span style="float: right;">This information is continued on an additional sheet <input type="checkbox"/></span>                                                                                                                                                                                                                                                                                         |                                                                                                                                             |
| <p>In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state.</p> |                                                                                                                                             |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                             |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                                                                             |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                              |                                                                                                                                             |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                                      |                                                                                                                                             |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                                                                             |
| For receiving Office use only<br><input checked="" type="checkbox"/> This sheet was received with the international application<br>Authorized officer <i>CSH</i>                                                                                                                                                                                                                                                                                                          | For International Bureau use only<br><input type="checkbox"/> This sheet was received by the International Bureau on:<br>Authorized officer |



## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|
| <b>A.</b> The indications made below relate to the microorganism referred to in the description<br>on page <u>36</u> , line <u>7</u>                                                                                                                                                                                                                                                                                                                                      |                                 |
| <b>B. IDENTIFICATION OF DEPOSIT</b> <span style="float: right;">Further deposits are identified on an additional sheet <input type="checkbox"/></span>                                                                                                                                                                                                                                                                                                                    |                                 |
| Name of depositary institution<br>The National Collections of Industrial and Marine Bacteria Limited (NCIMB)                                                                                                                                                                                                                                                                                                                                                              |                                 |
| Address of depositary institution (including postal code and country)<br><br>23 St Machar Drive<br>Aberdeen AB2 1RY<br>Scotland, UK                                                                                                                                                                                                                                                                                                                                       |                                 |
| Date of deposit<br>28 June 1994                                                                                                                                                                                                                                                                                                                                                                                                                                           | Accession Number<br>NCIMB 40669 |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) <span style="float: right;">This information is continued on an additional sheet <input type="checkbox"/></span>                                                                                                                                                                                                                                                                                         |                                 |
| <p>In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state.</p> |                                 |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                         |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                              |                                 |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                                      |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |

|                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                                                                                                                                                                                                                                                                        |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p style="text-align: center;">For receiving Office use only</p> <div style="border: 1px solid black; padding: 5px;"> <input checked="" type="checkbox"/> This sheet was received with the international application       </div> <div style="border: 1px solid black; padding: 5px; margin-top: 10px;">         Authorized officer <span style="float: right;"><i>ASnd</i></span> </div> | <p style="text-align: center;">For International Bureau use only</p> <div style="border: 1px solid black; padding: 5px;"> <input type="checkbox"/> This sheet was received by the International Bureau on:       </div> <div style="border: 1px solid black; padding: 5px; margin-top: 10px;">         Authorized officer       </div> |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                                                                                                                                                                                                     |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>A.</b> The indications made below relate to the microorganism referred to in the description<br>on page <u>36</u> , line <u>15</u>                                                                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                                                     |
| <b>B. IDENTIFICATION OF DEPOSIT</b> Further deposits are identified on an additional sheet <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                     |
| Name of depositary institution<br><b>The National Collections of Industrial and Marine Bacteria Limited (NCIMB)</b>                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                     |
| Address of depositary institution (including postal code and country)<br><br><b>23 St Machar Drive<br/>Aberdeen AB2 1RY<br/>Scotland, UK</b>                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                                     |
| Date of deposit<br><b>28 June 1994</b>                                                                                                                                                                                                                                                                                                                                                                                                                             | Accession Number<br><b>NCIMB 40670</b>                                                                                                                                                                                                                                                                              |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) This information is continued on an additional sheet <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                     |                                                                                                                                                                                                                                                                                                                     |
| In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state. |                                                                                                                                                                                                                                                                                                                     |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                  |                                                                                                                                                                                                                                                                                                                     |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                                                                                                                                                                                                     |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                                     |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                               |                                                                                                                                                                                                                                                                                                                     |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                                                                                                                                                                                                     |
| <div style="border: 1px solid black; padding: 2px;">For receiving Office use only</div> <div style="border: 1px solid black; padding: 2px;"><input checked="" type="checkbox"/> This sheet was received with the international application</div> <div style="border: 1px solid black; padding: 2px;">Authorized officer <span style="float: right; font-family: cursive;">C. J. N. D.</span></div>                                                                 | <div style="border: 1px solid black; padding: 2px;">For International Bureau use only</div> <div style="border: 1px solid black; padding: 2px;"><input type="checkbox"/> This sheet was received by the International Bureau on:</div> <div style="border: 1px solid black; padding: 2px;">Authorized officer</div> |

## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|
| <b>A.</b> The indications made below relate to the microorganism referred to in the description<br>on page <u>37</u> , line <u>11</u>                                                                                                                                                                                                                                                                                                                                     |                                 |
| <b>B. IDENTIFICATION OF DEPOSIT</b> <div style="text-align: right;">Further deposits are identified on an additional sheet <input type="checkbox"/></div>                                                                                                                                                                                                                                                                                                                 |                                 |
| Name of depositary institution<br>The National Collections of Industrial and Marine Bacteria Limited (NCIMB)                                                                                                                                                                                                                                                                                                                                                              |                                 |
| Address of depositary institution (including postal code and country)<br><br>23 St Machar Drive<br>Aberdeen AB2 1RY<br>Scotland, UK                                                                                                                                                                                                                                                                                                                                       |                                 |
| Date of deposit<br>28 June 1994                                                                                                                                                                                                                                                                                                                                                                                                                                           | Accession Number<br>NCIMB 40659 |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) <div style="text-align: right;">This information is continued on an additional sheet <input type="checkbox"/></div>                                                                                                                                                                                                                                                                                      |                                 |
| <p>In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state.</p> |                                 |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                         |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                              |                                 |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                                      |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                 |

|                                                                                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                                                        |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <div style="text-align: center;">For receiving Office use only</div> <div style="border: 1px solid black; padding: 5px;"> <input checked="" type="checkbox"/> This sheet was received with the international application </div> <div style="border: 1px solid black; padding: 5px; margin-top: 10px;"> Authorized officer <u>CSMD</u> </div> | <div style="text-align: center;">For International Bureau use only</div> <div style="border: 1px solid black; padding: 5px;"> <input type="checkbox"/> This sheet was received by the International Bureau on: </div> <div style="border: 1px solid black; padding: 5px; margin-top: 10px;"> Authorized officer </div> |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                 |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|
| <b>A.</b> The indications made below relate to the microorganism referred to in the description<br>on page <u>38</u> , line <u>23</u>                                                                                                                                                                                                                                                                                                                              |                                 |
| <b>B. IDENTIFICATION OF DEPOSIT</b> Further deposits are identified on an additional sheet <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                                                |                                 |
| Name of depositary institution<br>The National Collections of Industrial and Marine Bacteria Limited (NCIMB)                                                                                                                                                                                                                                                                                                                                                       |                                 |
| Address of depositary institution (including postal code and country)<br><br>23 St Machar Drive<br>Aberdeen AB2 1RY<br>Scotland, UK                                                                                                                                                                                                                                                                                                                                |                                 |
| Date of deposit<br>28 June 1994                                                                                                                                                                                                                                                                                                                                                                                                                                    | Accession Number<br>NCIMB 40664 |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) This information is continued on an additional sheet <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                     |                                 |
| In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state. |                                 |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                  |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                 |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                       |                                 |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                               |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                 |

|                                                                                                                                                                                                                                                                                                             |                                                                                                                                                                                                                                                                                                  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <div style="text-align: right; font-size: small;">For receiving Office use only</div> <div style="padding: 5px;"><input checked="" type="checkbox"/> This sheet was received with the international application</div> <div style="border-top: 1px solid black; padding-top: 5px;">Authorized officer </div> | <div style="text-align: right; font-size: small;">For International Bureau use only</div> <div style="padding: 5px;"><input type="checkbox"/> This sheet was received by the International Bureau on:</div> <div style="border-top: 1px solid black; padding-top: 5px;">Authorized officer</div> |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## INDICATIONS RELATING TO A DEPOSITED MICROORGANISM

(PCT Rule 13bis)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                 |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|
| <b>A.</b> The indications made below relate to the microorganism referred to in the description<br>on page <u>40</u> , line <u>13</u>                                                                                                                                                                                                                                                                                                                              |                                 |
| <b>B. IDENTIFICATION OF DEPOSIT</b>                                                                                                                                                                                                                                                                                                                                                                                                                                |                                 |
| Further deposits are identified on an additional sheet <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                                                                                    |                                 |
| Name of depositary institution<br>The National Collections of Industrial and Marine Bacteria Limited (NCIMB)                                                                                                                                                                                                                                                                                                                                                       |                                 |
| Address of depositary institution (including postal code and country)<br>23 St Machar Drive<br>Aberdeen AB2 1RY<br>Scotland, UK                                                                                                                                                                                                                                                                                                                                    |                                 |
| Date of deposit<br>28 June 1994                                                                                                                                                                                                                                                                                                                                                                                                                                    | Accession Number<br>NCIMB 40660 |
| <b>C. ADDITIONAL INDICATIONS</b> (leave blank if not applicable) This information is continued on an additional sheet <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                     |                                 |
| In respect of all designated states in which such action is possible and to the extent that it is legally permissible under the law of the designated state, it is requested that a sample of the deposited micro-organism be made available only by the issue thereof to an independent expert, in accordance with the relevant patent legislation, e.g. Rule 28(4) EPC, and generally similar provisions <i>mutatis mutandis</i> for any other designated state. |                                 |
| <b>D. DESIGNATED STATES FOR WHICH INDICATIONS ARE MADE</b> (if the indications are not for all designated States)                                                                                                                                                                                                                                                                                                                                                  |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                 |
| <b>E. SEPARATE FURNISHING OF INDICATIONS</b> (leave blank if not applicable)                                                                                                                                                                                                                                                                                                                                                                                       |                                 |
| The indications listed below will be submitted to the International Bureau later (specify the general nature of the indications e.g., "Accession Number of Deposit")                                                                                                                                                                                                                                                                                               |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                 |

|                                                                                                                                                                         |                                                                                                                                                          |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>For receiving Office use only</b><br><input checked="" type="checkbox"/> This sheet was received with the international application<br>Authorized officer <u>CSM</u> | <b>For International Bureau use only</b><br><input type="checkbox"/> This sheet was received by the International Bureau on:<br>Authorized officer _____ |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|

## CLAIMS

1. A polypeptide which is a water-soluble active derivative of a bacterial bifunctional penicillin binding protein, said penicillin binding protein being bound to the cell membrane when expressed in a bacterial cell and being capable of exhibiting both transglycosylase and transpeptidase activities and said derivative lacking a membrane anchoring sequence but retaining the capability to exhibit one or both of said enzymic activities.
2. A polypeptide according to claim 1 which amino acid sequence is identical to, or substantially similar to, SEQ ID NO: 2, 4, 6, 12 or 13 in the Sequence Listing.
3. A polypeptide which is a transglycosylase deficient derivative of a bacterial bifunctional penicillin binding protein, said penicillin binding protein being bound to the cell membrane when expressed in a bacterial cell and being capable of exhibiting both transglycosylase and transpeptidase activities and said derivative lacking the capability to exhibit transglycosylase activity but retaining the capability to exhibit transpeptidase activity.
4. A polypeptide according to claim 3 wherein the said derivative is lacking transglycosylase activity because of a mutation or deletion in the second conserved region of the gene coding for said polypeptide.
5. A polypeptide according to claim 3 which amino acid sequence is identical to, or substantially similar to, SEQ ID NO: 7, 8, 9, or 10 in the Sequence Listing.

6. A polypeptide according to claim 1 or 3 wherein the bacterial cell is an *Escherichia coli* cell or a *Streptococcus pneumoniae* cell.
- 5 7. A polypeptide comprising (a) a first polypeptide according to claim 1 or 3; and (b) an additional polypeptide which allows binding to an affinity matrix; there being a cleavage site between said polypeptides.
- 10 8. A polypeptide according to claim 7 wherein the additional polypeptide is glutathione-S-transferase or a polypeptide substantially similar to glutathione-S-transferase.
- 15 9. A polypeptide according to claim 7 wherein the additional polypeptide is a polypeptide rich in histidine residues.
- 20 10. An isolated and purified DNA molecule which has a nucleotide sequence coding for a polypeptide according to claim 1, 3 or 7.
- 25 11. A DNA molecule according to claim 10, which nucleotide sequence is identical to, or substantially similar to, SEQ ID NO: 1, 3 or 5 in the Sequence Listing.
- 30 12. A replicable expression vector which carries and is capable of mediating the expression of a DNA molecule according to claim 10.
13. A vector according to claim 12 which is the vector  
pARC0558 (NCIMB No. 40666),  
pARC0559 (NCIMB No. 40667),  
pARC0512 (NCIMB No. 40665),  
pARC0438 (NCIMB No. 40661),  
pARC0468 (NCIMB No. 40662),  
pARC0469 (NCIMB No. 40663),

5 pARC0571 (NCIMB No. 40668),  
pARC0593 (NCIMB No. 40670),  
pARC0392 (NCIMB No. 40659),  
pARC0499 (NCIMB No. 40664), or  
pARC0400 (NCIMB No. 40660).

14. A cell harbouring a vector according to claim 12.
- 10 15. A process for production of a polypeptide which is a derivative of penicillin binding protein, comprising growing a cell according to claim 14 in or on a culture medium for expression of the polypeptide and optionally recovering the polypeptide.
- 15 16. A process for the production of a water soluble polypeptide according to claim 1 which comprises culturing *Escherichia coli* cells harbouring an expression vector wherein a DNA coding sequence for said polypeptide is under the control of an isopropyl thiogalactoside (IPTG) - inducible promoter, said culturing being carried out in the presence of a sub-optimal concentration of IPTG  
20 for induction of the said promoter and at a temperature in the range of 20 to 24°C.
- 25 17. A method of identifying an antibody capable of binding a bacterial bifunctional penicillin binding protein which includes the step of employing a polypeptide according to claim 1 or 3 in an antibody binding assay and selecting antibodies that bind to the polypeptide.
- 30 18. A method of assaying for compounds which bind to a penicillin binding protein, said method comprising (a) contacting a polypeptide according to claim 1, 3 or 7 with a compound to be investigated; and (b) detecting whether said compound binds to the penicillin binding protein.



19. A method of assaying for compounds which bind to a penicillin binding protein, said method comprising (a) culturing cells according to claim 14; (b) lysing the said cells and isolating the crude cell extract; (c) exposing the said cell extract to potential inhibitors of a penicillin binding protein; (d) introducing an agent, known to bind a penicillin binding protein, to the said cell extract; (e) removing the unbound fraction of said agent; and (f) assaying the presence of said agent remaining in the cell extract.
20. A method of assaying for compounds which bind to a penicillin binding protein, said method comprising (a) exposing a polypeptide according to claim 1, 3 or 7, immobilised on a solid support, to a potential inhibitor of a penicillin binding protein; (b) exposing an agent, known to bind a penicillin binding protein, to the immobilised polypeptide; (c) removing the unbound fraction of said agent; and (d) assaying the presence of said agent bound to the immobilised polypeptide.
21. A method of assaying for compounds which bind to a penicillin binding protein, said method comprising (a) exposing a polypeptide according to claim 1, 3 or 7 to a potential inhibitor of a penicillin binding protein; (b) exposing the said polypeptide to an agent, known to bind a penicillin binding protein, which agent is immobilised on a solid support; and (c) assaying the presence of polypeptide bound to the immobilised agent.
22. A method of assaying for compounds which bind to the transglycosylase domain of a penicillin binding protein, said method comprising (a) exposing the transglycosylase domain of a polypeptide according to claim 1 or 7, said polypeptide being immobilised on a solid support, to a potential inhibitor of the transglycosylase activity of a penicillin binding protein; (b) exposing

an agent, known to bind the transglycosylase domain of a penicillin binding protein, to the immobilised polypeptide; (c) removing the unbound fraction of said agent; and (d) assaying the presence of said agent bound to the immobilised polypeptide.

5

23. A method of assaying for compounds which bind to the transglycosylase domain of a penicillin binding protein, said method comprising (a) exposing the transglycosylase domain of a polypeptide according to claim 1 or 7 to a potential inhibitor of a penicillin binding protein; (b) exposing the said polypeptide to an agent, known to bind to the transglycosylase domain of a penicillin binding protein, which agent is immobilised on a solid support; and (c) assaying the presence of polypeptide bound to the immobilised agent.
- 10
- 15
24. A method according to any one of claims 19 to 23 wherein the agent known to bind a penicillin binding protein is a monoclonal antibody.
- 20
25. A method according to any one of claims 19 to 23 wherein the agent known to bind a penicillin binding protein is a labelled antibiotic compound.
- 25
26. A method of determining the protein structure of a penicillin binding protein, characterized in that a polypeptide according to claim 1 or 3 is utilized in X-ray crystallography.

1 / 14

Fig. 1





Fig. 2



2 / 14

Fig. 3




Fig. 4




Panel A

Panel B

3 / 14


Fig. 5



Panel A

Panel B

Fig. 6



4 / 14






Fig. 8

NcoI - NruI fragment from pARC0547 cloned into  
NcoI -EcoRV sites of pARC0534.

5 / 14

**Fig. 7**Hydrophobicity plot *E. coli* PBP 1B (1-150 a.a)**Fig. 9**

6 / 14

Fig. 10





Fig. 11

Hydrophobicity Plot *S. pneumoniae* PBP1A (1-100 a.a)



7 / 14

Fig. 12






Fig. 13



1 2 3 4 5

Panel A



1 2 3 4 5

Panel B

8 / 14

E.1B 190 DPRLITMISSPNGEQR.LFVPRSGFPDLLVDTLLATEDRHFYEHGDISLYSIGRAVLANLTAGRTVQGASTLTQQLVKNLFLS...  
 S.1A 55 IYDNKNQLIADLGSERRVNAQANDIPTDLVKAIVSIEDHRFPDHRGIDTIRILGAFRLNLSN.SLQGGHALTQQLIKLITYPSTST  
 E.1A 50 IYSADGELIAQYGEKRRIPVTLQDIPPEMVKAFIATEDSRFYEHHGVDPVGIFRAASVALPSGHASQGASTITQQLARNFFLS...  
 H.inf 50 IYTADGKLGGEVGEQRRIPVKLDVPQRLIDAFLATEDSRFPYDHHGLDPIGIARALFVYFYSNGGASQGASTITQQLARNFFLI...

REGION 2

REGION 1

E.1B SERSYWRKANEAYMALINDARYSKDRILELYMNEVYLGQSGDNEIRGFPLASLYYFGRPVEELSLDQOALLVGMVKGASIYNPWR  
 S.1A SDQTSRKAQEAWLAIQLEQKATQKQEIILTYYNKVMS....NGNYGMQTAQNYYGKDLNNLSLPQLALLAGMPQAPNQYDPYS  
 E.1A PERTLMRKIKEVELAIRIEQLLTDEILELYLNKIYLG....YRAYGVGAAQVYFGKTVQDLTNEMAVIAGLPKAPSTFNPLY  
 H.inf SEKTIKKFKEAVLAVEIENTLNKQEIILELYLNKIFLG....YRSYGVAAAQTYFGKSLNLLTSEMALTAGLPFAFSTMNALY

REGION 3

E.1B NPKLALERRNLVLRLLQOQQII 386  
 S.1A HPEAAQDRRNVLVLSMKNQGYI 242  
 E.1A SMDRAVARRNVLVLSRMLDEGYI 235  
 H.inf SLKPSERRNVVLSRMLDEKYI 236

Fig. 14

9/14

Fig. 15





Fig. 17

10 / 14





Fig. 18



11/14

Fig. 19





Fig. 20



12 / 14

Fig. 21





Fig. 22



13 / 14

Fig. 23






Fig. 24



14 / 14

Fig. 25





## INTERNATIONAL SEARCH REPORT

International application N .

PCT/SE 95/00761

## A. CLASSIFICATION OF SUBJECT MATTER

IPC6: C07K 14/245, C12N 9/24, C12N 9/52, G01N 33/53  
According to International Patent Classification (IPC) or to both national classification and IPC

## B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC6: C07K, C12N

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

SE,DK,FI,NO classes as above

Electronic data base consulted during the international search (name of data base and, where practicable, search terms used)

WPI, EPOC, BIOSIS, MEDLINE, SCISEARCH, PATENT CITATION INDEX, EMBL/GENSEQ/D  
DBJ

## C. DOCUMENTS CONSIDERED TO BE RELEVANT

| Category* | Citation of document with indication, where appropriate, of the relevant passages                                                                                                                                                                                                                                                 | Relevant to claim No. |
|-----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------|
| X         | JOURNAL OF BACTERIOLOGY, Volume 172, No 1, January 1990, TANNEKE den BLAAUWEN et al, "Interaction of Monoclonal antibodies with the Enzymatic Domains of Penicillin-Binding Protein 1b of Escherichia coli", page 63 - page 70, the whole document<br>--                                                                          | 17-25                 |
| X         | THE JOURNAL OF BIOLOGICAL CHEMISTRY, Volume 8, 1993, Robert A. Nicholas et al, "Penicillin-binding Protein 1B from Escherichia coli Contains a Membrane Association Site in Addition to Its Transmembrane Anchor", page 5632 - page 5641, page 5633, column 2, line 48-51; page 634, column 2, line 64 - page 5635, line 13<br>-- | 1-2,6-10,<br>12-16    |

☒ Further documents are listed in the continuation of Box C.☒ See patent family annex.

## \* Special categories of cited documents:

- \*A\* document defining the general state of the art which is not considered to be of particular relevance
- \*E\* earlier document but published on or after the international filing date
- \*L\* document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)
- \*O\* document referring to an oral disclosure, use, exhibition or other means
- \*P\* document published prior to the international filing date but later than the priority date claimed

- \*T\* later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
- \*X\* document of particular relevance: the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
- \*Y\* document of particular relevance: the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art
- \*Z\* document member of the same patent family

Date of actual completion of the international search

28 February 1996

Date of mailing of the international search report

28 -02- 1996

Name and mailing address of the ISA/  
Swedish Patent Office  
Box 5055, S-102 42 STOCKHOLM  
Facsimil No. +46 8 666 02 86

Authorized officer

Patrick Andersson  
Telephone N. +46 8 782 25 00

## INTERNATIONAL SEARCH REPORT

International application No.

PCT/SE 95/00761

## C (Continuation). DOCUMENTS CONSIDERED TO BE RELEVANT

| Category* | Citation of document, with indication, where appropriate, of the relevant passages                                                                                                                                                                                                                                                                                                                                       | Relevant to claim No. |
|-----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------|
| X         | EP 0505151 A2 (ELI LILLY AND COMPANY),<br>23 Sept 1992 (23.09.92), page 3, line 1 - line 3                                                                                                                                                                                                                                                                                                                               | 1,26                  |
| A         | --                                                                                                                                                                                                                                                                                                                                                                                                                       | 2-25                  |
| A         | EUROPEAN JOURNAL OF BIOCHEMISTRY, Volume 147, 1985,<br>Jenny K. BROOME-SMITH et al, "The nucleotide<br>sequences of the ponA and ponB genes encoding<br>penicillin-binding proteins 1A and 1B of<br>Escherichia coli K12" page 437 - page 446                                                                                                                                                                            | 1-26                  |
| A         | Dialog Information Services, file 154, MEDLINE,<br>Dialog accession no, 08300977, Medline accession<br>no. 93010977, Martin C et al: "Relatedness of<br>penicillin-resistant Streptococcus pneumoniae<br>isolated in south Africa and Spain"; & EMBO J<br>(ENGLAND) Nov 1992, 11 (11) p3831-6                                                                                                                            | 1-26                  |
| A         | Dialog Information Services, file 154, MEDLINE,<br>Dialog accession no. 08187042, Medline accession<br>no. 92325042, Martin C et al: "Nucleotide sequences<br>of genes encoding penicillin-binding proteins from<br>Streptococcus pneumoniae and Streptococcus oralis<br>with high homology to Escherichia coli penicillin-<br>binding proteins 1a and 1b"; & J Bacteriol (UNITED<br>STATES) Jul 1992, 174 (13) p4517-23 | 1-26                  |
|           | --<br>-----                                                                                                                                                                                                                                                                                                                                                                                                              |                       |

### Information on patent family members

**PCT/SE 95/00761**

**05/02/96**

Form PCT/ISA/210 (patent family annex) (July 1992)

**THIS PAGE BLANK (USPTO)**